

# Buffer Capacity

# Buffers:

- Buffer solutions are solutions that can resist changes in pH upon addition of small amounts of acid/base.
- Common buffer mixtures contain two substances: conjugate acid and a conjugate base .
- Together the two species (conjugate acid and conjugate base) resist large changes in pH by absorbing the H<sup>+</sup> ions or OH<sup>-</sup> ions added to the system.

# How buffers resist the change in pH:

1. When **H<sup>+</sup> ions** are added to the buffer system they will react with the **conjugate base** in the buffer as following:



2. When **OH<sup>-</sup> ions** are added they will react with the **conjugate acid** in the buffer as following:



**→ Thus, the buffer is effective as long as it does not run out of one of its components.**

(There are enough conjugated base and conjugated acid to absorb the H<sup>+</sup> ions or OH<sup>-</sup> ions added to the system respectively).

# Buffer Capacity (Theoretically):

- Quantitative measure of buffer resistance to pH changes is called **buffer capacity**.
- Buffer capacity can be defined in many ways, **it can be defined as:**
  - The number of moles of H<sup>+</sup>/OH<sup>-</sup> ions that must be added to one liter of the buffer in order to decrease /increase the pH by one unit respectively.
- The instantaneous buffer capacity is **expressed as  $\beta$**  and can be derived from Henderson Hasselbalch equation:

$$\beta = \frac{2.3 K_a [H^+][C]}{(K_a + [H^+])^2}$$

This equation will give you an overview of buffer capacity in both directions (when adding H and OH).

- **Where:**  $\beta$  = the buffer capacity, [H<sup>+</sup>] = the hydrogen ion concentration of the buffer, [C] = concentration of the buffer and  $K_a$  = acid dissociation constant.
- **Note:** The buffer capacity is **directly proportional** to the buffer concentration

# Buffer Capacity (Practical):

## □ Buffer capacity of acid and alkaline direction:

→ Buffer capacity<sub>a</sub> ( $BC_a$ ) = the number of moles of  $H^+$  that must be added to one liter of the buffer in order to decrease the pH by one unite.

This called buffer capacity in the ACID direction.

$$BC_a = \frac{9[HA] [A^-]}{10 [HA] + [A^-]}$$

→ Buffer capacity<sub>b</sub> ( $BC_b$ ) = the number of moles of  $OH^-$  that must be added to one liter of the buffer in order to increase the pH by one unite.

This called buffer capacity in the ALKAILNE direction.

$$BC_b = \frac{9[HA] [A^-]}{10 [A^-] + [HA]}$$


# Buffer capacity in acid and base direction:




$\text{H}_3\text{O}^+$  added, equilibrium position shifts to the left




$\text{OH}^-$  added, equilibrium position shifts to the right






Buffer solution  
after addition  
of strong acid




Buffer solution  
equimolar in  
acid and base



Buffer solution  
after addition  
of strong base




# Buffer capacity in acid and base direction:



**Figure 6** (a) Ethanoic acid buffer with a strong base added (b) Ethanoic acid buffer with a strong acid added. The pH changes quickly once all of the available buffer is depleted.

# Buffer capacity in acid and base direction:

- ❑ The figure represents the buffer capacity when the buffer is titrated in **both directions**.
- ❑ The buffer will react differently upon titration with acid and base **depending on the specific components of the buffer**.



The buffer capacity curve of 0.05M Tricine buffer, pH 7.5 ( $pK_a=8.15$ )



# Buffer capacity in acid and base direction:

The buffer will react differently upon titration with acid and base depending on the specific components of the buffer. This relation is shown in the table below:

Buffer	[HA]	[A <sup>-</sup> ]	Buffer capacity (resistance)
pH = Pka	Equal	Equal	Equal
pH lower than Pka	Higher	Lower	High buffer capacity in alkaline direction. [OH] ions
pH higher than Pka	Lower	Higher	High buffer capacity in acid. [H <sup>+</sup> ] ions

**Example:** Calculate the practical buffer capacity in the acid directions of a 0.1M and 0.2M acetate buffer, pH 5, pKa = 4.76.

First, calculate the concentration of the weak base and its conjugated acid that make up the buffer with 0.1M:

$$pH = pKa + \log \frac{[A^-]}{[HA]} \rightarrow 5 = 4.76 + \log \frac{[y]}{[0.1-y]} \rightarrow 0.24 = \log \frac{[y]}{[0.1-y]} \rightarrow \text{Anti log for both sides}$$

$$1.74 = \frac{y}{0.1-y} \rightarrow y=0.063\text{M. SO: } [A^-]=0.063\text{M} , [HA]=0.037\text{M}$$

Second: Calculate the practical buffer capacity in both directions

$$BC_a = \frac{9 [HA][A^-]}{10 [HA]+[A^-]} \rightarrow \frac{9 \times 0.037 \times 0.063}{(10 \times 0.037) + 0.063} \rightarrow \frac{0.021}{0.433} = 0.048\text{M } [H^+]$$

- **Note:** do the same calculation for the same buffer when its concentration equals 0.2M.

# Practical Part

# Objective:

---

- To understand the concept of buffer capacity.
- To determine the capacity of two acetate buffers in the acid directions.

# Method:

- You are provided **0.1M** and **0.2M** acetate buffer (pH=5).
- In one beaker add **10 ml** of the 0.1M acetate buffer, and in the other add 0.2M acetate buffer
- Titrate the two beakers by adding **0.5 ml of 0.1 M HCl** from the burette and determine the pH of the solution after each addition.
- Continue adding the acid/base until the pH drops by two units from your initial reading
- Record the values in the titration table.

# Results:

- Plot the capacity curve (pH against the volume (ml) of 0.1M HCl).
- For both buffers, determine the practical buffer capacity in the acid direction from the graph and the formula then summarize your value in the table:

Acetate buffer	Practical capacity (from the formula)	Practical capacity (from the curve)
0.1M	0.048M $[H^+]$	
0.2M		

# Results:

- **To determine the capacity from the graph:**
  - a) Find the ml of 0.1M HCl needed to drop the pH one unit from the initial reading value.
  - b) Then find the final concentration of the HCl.

## Example from the curve:

3.8 ml of 0.1M HCl is needed to drop the pH from 3.8 to 2.8 of 10 ml of acetate buffer. Thus:

$$C_1 \times V_1 = C_2 \times V_2$$

$$0.1\text{M} \times 3.8 \text{ ml} = ? \text{ M} \times 13.8 \text{ ml}$$

$$= 0.027\text{M} [\text{H}^+]$$




Figure: 0.1M acetate buffer capacity in the acid directions