



Characterization of DNA by Spectrophotometric Assay and Melting Temperature (T_m)



After DNA isolation:



Characterization of extracted DNA by spectrophotometric assay:

1. DNA concentration.
2. DNA purity.

By measuring the absorption of ultraviolet light (UV).

- DNA has maximum absorption at **260nm**. **WHY?**

1st. UV for quantification of nucleic acid concentration :

- Is determined by measuring absorbance at **260 nm**. **WHY?**
- At 260 nm double-stranded DNA has specific absorption coefficient of $0.02 (\mu\text{g/ml})^{-1}\text{cm}^{-1}$.
- So: **→ Concentration of DNA** = $(A_{260} / \epsilon L) \times \text{Dilution Factor (DF)}$.

Beer-Lambert Law:

$$A = \epsilon cl$$

2nd. DNA purity:

1. To detect nucleic acid purity from proteins contamination:

→ Calculate A_{260}/A_{280} WHY?

• Highly purified DNA samples have a A_{260}/A_{280} nm ratio of (1.8-1.9).

→ What if the ratio is below 1.8? What that means?


→ What if the ratio is higher than 1.9? What that means?

2. To detect nucleic acid purity from carbohydrates, peptides, ethanol or any organic compounds:

→ Calculate A_{260}/A_{230} WHY?

• Purified DNA samples have a A_{260}/A_{230} nm ratio of (2-2.2).

DNA and protein absorption spectrum:



***What is the effect of the contaminants on DNA concentration?**

*** What if the samples contaminated by proteins or organic compound?**

3^{ed}. GC content:

- The two strands of a DNA molecule can be dissociated ("**melted**") into single strands by **heat or altered pH**, which breaks the hydrogen bonds between complementary bases (A=T and G=C).
- What that process called ?
- Hyperchromic and hypochromic effect ?
- The **melting temperature (T_m)** is the temperature at which **50% of the DNA is unpaired** (denatured).
- Melting temperature profile ?
- Is the T_m same for all DNA molecules ? WHY?
- What is the important of knowing T_m of DNA?

DNA melting curve:




FIGURE 4.4 DNA melting curve. A melting curve of DNA showing T_m (the melting temperature) and possible molecular conformations for various degrees of melting.

3^{ed}. GC content:

- GC content can be calculated by generating T_m profile (DNA melting curve).

$$\%(G+C) = 2.44 (T_m - 69.3)$$

Relationship between T_m and GC%:





FIGURE 4.5 Effect of G-C content on DNA melting temperature. T_m increases with increasing percent of G + C.



What do you notice about the GC content in relation to T_m ?

Practical Part

Aim:

- Determination the concentration and purity of extracted DNA using UV spectrophotometer. [SEP]
- Determination of DNA melting temperature and GC content percentage.

Principle:

- DNA and proteins have a maximum absorbance at 260 and 280 respectively.
- dsDNA will be separated to ssDNA by heat (denaturation).
- O.D at 260 nm will increase during denaturation... Why?
- Temperature for midpoint of denaturation gives T_m .
- The DNA of each species has a specific melting curve.. Why?

Results:

1. Characterization of DNA by Spectrophotometric Assay:

Wavelength (nm)	Absorbance of DNA	
	Blood	Plant
230		
260		
280		

→ Find out the concentration of DNA using the following equation:

$$\text{Concentration of DNA} = (A_{260} / \epsilon L) \times \text{Dilution Factor.}$$

→ Determine the purity of the DNA.

Results:

2. Melting Temperature of DNA: [L] [SEP]

Temperature (°C)	DNA Absorbance at 260 nm	
	Blood	Plant
25		
50		
60		
70		
Boiling		

→ Plot the value of absorbance vs. temperature and calculate the T_m for sample DNA.

→ Find out the GC content of your sample using the following formula:

$$\text{Percent of G + C} = (T_m - 69.3) \times 2.44$$

Home Work:

- Watch the following videos :
 - <https://www.youtube.com/watch?v=wXiiTW3pflM>
 - https://www.youtube.com/watch?v=U2-5ukpKg_Q
- ➔ What is agarose gel electrophoresis?