
BCH312 [Practical]

### **Dilution of Prepared Solutions**

# **Dilution of Solution :**


□ Dilution of solution: means to add more <u>solvent</u> without the addition of more solute → To make it less concentrated.

- 1. Volume to volume dilutions (ratio).
- 2. Preparing dilutions by using the  $V_1XC_1=V_2XC_2$  formula.
- 3. Serial Dilutions.




# (1) Volume to volume dilutions (ratio):

- □ This type of dilutions describes the ratio of the **solute** to the **final volume** of the dilute solution.
- □ For example: to make 1:10 dilution of 1M NaCl solution, <u>one part</u> of the 1M NaCl solution, should be mixed with <u>nine parts</u> of water, for a <u>total of ten parts</u>.
- □ Therefore 1:10 dilution means  $\rightarrow$ 1 part of 1M NaCl + 9 parts of water.
- **Thus:**




- $\rightarrow$  if 10 ml of the 1:10 dilution was needed, then 1ml of 1M NaCl should be mixed with 9 ml of water.
- → if 100 ml of 1:10 dilution was needed, then 10 ml of the 1M NaCl should be mixed with 90 ml of water. [The final concentration of NaCl in both cases will be 0.1 M (1/10) = 0.1]



### **Example:** How to Prepare 20 ml of 2:10 dilution from 7 M solution (A) ?

how many ml of 7 M solution A we need to make 20 ml of 2:10 A solution?



= (2 X 20) / 10 = 4 ml

So, 4 ml from solution (A) of 7 M is needed and complete volume up to 20 ml (adding 16 ml water). Note: [16 ml water= 20 ml -4 ml].

> How to Know the concentration of solution A after dilution? First we will find the DILUTION FACTOR by the following : Dilution factor (D.F) = final volume / aliquot volume =10/2 = 5Then we will <u>divide</u> the stock concentration (before dilution) by the D.F: 7/5 = 1.4M

Note: To find out the stock concentration you will **<u>multiply</u>** the diluted concentration by the D.F

### (2) Preparing dilutions by using the V1XC1=V2XC2 formula:

- Sometimes it is necessary to use one solution to make a <u>specific amount</u> of a more dilute solution.
- To do this the following formula can be used:

$$V_1 X C_1 = V_2 X C_2$$

- □ Where:
- > V1= Volume of starting solution needed to make the new solution (volume of stock solution).
- >  $C_1$ = Concentration of starting solution (stock solution).
- $\triangleright$  V2= Final volume of new solution.
- $\succ$  C<sub>2</sub>= Final concentration of new solution.

## **Example:** Make 5ml of 0.25M solution from a 1.0M solution?

how many ml of 1M solution we need to make 5 ml of 0.25M solution?


→ V1XC1=V2XC2 Where: V1 = ? , C1= 1M , V2= 5ml , C2= 0.25M

So:  $(V1) \ge (1M) = (5ml) \ge (0.25M)$  $\rightarrow V1 = (5 \ge 0.25)/1 = 1.25 ml$ 

So 1.25ml of the 1M solution is needed (starting solution) then complete the volume up to 5 ml by diluent (generally water).

# (3) Serial Dilutions :

- □ It is a stepwise dilution of a solution, where the **<u>dilution factor is constant at each step.</u>**
- □ The source of dilution material for each step comes from the diluted material of the previous step.




#### 1:10

Dilution factor (D.F) = final volume / aliquot volume = 10/1 = 10 (for each step)

## Find out the concentration of the diluted solutions:

**Dilution factor (D.F) = final volume / aliquot volume = 10/1 = 10 (for each step)** 



From the **more** concentrated solution to the **lower** one



Starting with a 2.0 M stock solution of hydrochloric acid, prepare four standard solutions by serial dilution of the following Molarity respectively 1 M, 0.5 M, 0.25 M, 0.125 M. [with 1:2 dilution] ?


→ Dilution factor (D.F) = final volume / aliquot volume =  $2/1 = 2 \rightarrow 1:2$ 

#### -To prepare standard solution 1:

1 ml of the stock 2.0M solution is needed and volume made up to 2 ml with distilled water (never forget to mix properly).

#### -To prepare standard solutions 2-4:

1 ml of the previously diluted solution is taken, and the volume will made up to 2 ml by the addition of distilled water.



how to calculate the concentration of the diluted solutions if they unknown ?

First: find the D.F: Dilution factor (D.F) = final volume / aliquot volume = 2/1 = 2

→ Second: divide the previous solution concentration by the D.F:

-concentration of solution 1 = 2.0 M stock solution /2 = 1 M -concentration of solution 2 = 1M/2 = 0.5 M -concentration of solution 3 = 0.5M/2 = 0.25 M -concentration of solution 4 = 0.25/2 = 0.125 M

# Preschieel Pert



□ To dilute prepared solutions by different methods.



### Solution dilutions:

(1).....

**Prepare 50 ml with 1:20 dilution using the 0.08 M NaOH solution you previously prepared.** 

□ <u>Calculation:</u>

.....

.....

→ To prepare the 1:20 dilution ......ml of the starting solution (0.08M NaOH) is needed and volume made up to a final volume of ......ml.



(2).....

### **Prepare 50 ml of 0.2M HCl from the previously 0.4M HCl solution you previously prepared.**

□ <u>Calculation:</u>

.....

.....

→ To prepare the 0.2M HCl .....ml of the starting solution (0.4M HCl) is needed and volume made up to a total volume of .....ml by adding water.



(3).....

Starting with a 2.0 M stock solution of hydrochloric acid, prepare 4 ml of four standard solutions (1 to 4) of the following Molarity respectively (dilution 1:4):
(1) ...... M (2) ..... M (3) ..... M (4) ..... M.

**Calculation:** 

.....

.....

→ To prepare standard solution 1: ...... ml of the stock 2.0M solution is needed and volume made up to ...... ml with distilled water.

→ To prepare standard solution 2-4: ...... ml of the previously diluted solution  $(8.00 \times 10-2 \text{ M})$  is taken and volume is made up to a final volume of ...... ml by the addition of distilled water.