



# **Estimation of total protein in milk and egg using turbidimetric method**

# Diet protein sources:


- Proteins in human diet are derived from two main sources namely:
  1. **Animal proteins** (e.g. egg, milk, meat and fish).
  2. **Plant proteins** (e.g. pulses, cereals, nuts, beans and soy products).
- Animal proteins are more “**biologically complete**” than vegetable proteins.



- **Food analysts** are interested in knowing the total concentration, type, molecular structure and functional properties of the proteins in foods.
- Proteins are also the major structural components of many natural foods, often determining their overall texture.
- Isolated proteins are often used in foods as **ingredients** because of their unique functional properties, i.e., **their ability to provide desirable appearance, texture.**

# Milk proteins:

- Normal bovine milk contains **30–35 grams of protein per liter**.
- Primary group of milk proteins are the **caseins 80%**.
- **All other proteins** found in milk are grouped together under the name of whey proteins.
- The major whey proteins in cow milk are **beta-lactoglobulin and alpha-lactalbumin**.



# Egg proteins:



- They supply all essential amino acids for humans (a source of 'complete protein').
- **Egg white** consists primarily of about 90% water into which is dissolved 10% proteins (including albumins, mucoproteins, and globulins).
- Unlike the **yolk**, which is high in lipids (**fats**), egg **white** contains almost **no fat**, and the carbohydrate content is less than 1%.

Fat & Protein	Yolk	White
Cholesterol	389 mg	0 mg
O3 FA's	71.8 mg	0.0 mg
Saturated Fat	3.0 g	0.0 g
Mono Fat	3.7 g	0.0 g
Protein	5.0 g	22.7 g

# Practical Part


---

# Objective:

- Determine the total protein content in milk and egg using turbidimetric method (by sulfosalicylic acid).

# Turbidimetric method:


- Determination of total protein by measurement of **protein turbidity** produce by mixed with an anionic organic acid such as sulfosalicylic acid , TCA , or benzethonium chloride.
- These methods are sensitive, but the reagent **does not** react **equally** with each protein fraction.
- Proteins are precipitated as fine particals, turbidity is measured spectrophotometry.  
→ The higher protein concentration , the higher turbidity.






# Principle:

- **Sulfosalicylic acid** is an **anion(-)** which neutralizes the protein **cations(+)** leading to its precipitation (pH in highly acidic media, the protein will be positively charged, which is attracted to the acid anions that cause them to precipitate).
- Then the radiation of a wavelength which is **not absorbed** by the solution is made to pass through the suspension and the apparent absorption will be solely because of the scattering by the particles. → (The higher protein concentration, the lower transmittance value).




Increased concentration



# Principle:

- So, the transmitted light will have lower intensity as compared to that of the **incident light**.
- As a result, if the intensity of the transmitted light is measured, it will give an idea of the **number of particles in the suspension -proteins-**. (inversely related)



# Method:

- Set up a series of test tube as follows, label from 1- 6:

<b>Tube</b>	<b>Protein solution</b>	<b>Water</b>
1	4.5	<b>1.5</b>
2	3	<b>3</b>
3	2.4	<b>3.6</b>
4	1.5	<b>4.5</b>
5	0.9	<b>5.1</b>
6	0.3	<b>5.7</b>


# Method cont':

- Label a fresh set of test tubes 1 to 6. blank, egg sample and protein sample :
- Add 8 ml of sulphosalisalic acid to each test tube.
- Into tube 1 pipette 2 ml of protein solution 1(that you prepared before).
- Into tube 2 pipette 2 ml of protein solution 2 ... etc
- In the blank add 2ml water.
- Add 0.5 ml of egg sample and 1.5 ml water.
- Add 2 ml of milk sample.
- Mix the content of each tube well and allow to stand for five minutes.
- Using solution 7 (Blank) to set transmittance 100% at 500nm.
- Then use solutions from 1-6, to recorded respective transmittance of each suspension.
- record your results.

# Results :

- Plot transmittance against protein concentration on semi-logarithm paper (standard curve).
- Read the protein concentration of the “unknown samples” from the standard curve.
- Calculate the concentration of protein in the original sample (g/100 ml)

Tube	Transmittance at 500nm	Protein concentration (mg/dl)
B		
1		
2		
3		
4		
5		
6		
Milk sample		
Egg sample		



# Calculations :

- The concentration from the standard curve (mg/dl) x dilution factor= ----- mg/dl
- Dilution factor:
- Egg=.....
- Milk =.....