


Blood Biochemistry BCH 471 [Practical]


**Lab (4) Hemolyzing Agents & Detection of Blood**



# Blood Hemolysis

---

- **Hemolysis** (from the Greek **Hemo**: meaning blood, **lysis**, meaning to break open).
- It is the **breaking open** of red blood cells and the release of hemoglobin and the red cell contents into the surrounding fluid (plasma).
- Hemolysis may occur *in vivo* or *in vitro*.



# Hemolysis *in-vivo*

---

- **Conditions that can cause hemolysis include:**
  1. Immune reactions (e.g. during blood transfusion)
  2. Infections (e.g. malaria infection)
  3. Medications (e.g. cisplatin, a chemotherapy drug)
  4. Toxins and poisons
  
- Since the potassium concentration inside red blood cells is much higher than in the plasma, elevated potassium in plasma is usually found in biochemistry **tests of hemolyzed blood**.

# Hemolysis *in-vitro*


---

- **Factors that can cause hemolysis include:**

1. Improper technique during collection (e.g. incorrect needle size, excessive suction)
2. pH imbalance (addition acid or base)
3. Placing RBCs in a hypotonic solution

- **When blood hemolysis should be done?**

1. Breaking down RBCs to release their content
2. Estimation of haemoglobin
3. To obtain erythrocyte-free preparation of leukocyte and platelet



**Note:** In this lab blood hemolysis will be done by using hypotonic solutions and pH imbalance.

# Osmosis and Osmotic Pressure

## Osmosis:

It is the diffusion of solvent molecules across a semi-permeable membrane into a region of higher solute concentration.

- Once an *equilibrium* is reached the flow of water stops.


## Osmotic pressure:

The pressure exerted by a solvent passes through a semi-permeable membrane in osmosis.

## Tonicity:

Is the ability of an extracellular solution to make water move into or out of a cell by osmosis.

- Accordingly, there are three types of solutions: **Isotonic**, **Hypotonic**, and **Hypertonic**.





# Isotonic

---

- A solution that has the same solutes concentration as the normal cells of the body and the blood, having **equal osmotic pressure**.
- Example of Isotonic solution is **sodium chloride 0.9% (normal saline)**, have the same osmotic pressure as serum and they do not affect the membranes of the RBCs.
- In hospitals, intravenous fluids are isotonic.

Solute inside the cell = Solute outside the cell




# Hypotonic

---

- In a hypotonic solution, there is a lower concentration of solute outside a cell, creating an environment with **lower osmotic pressure** than what is contained within the cell.
- The RBCs will burst or hemolyzed.
- Any concentration of NaCl that is **lower than 0.9%**, will be considered hypotonic for cells.

Solute outside the cell < Solute inside the cell




# Hypertonic

---

- In a hypertonic solution, there is a higher concentration of solute outside a cell, creating an environment with **higher osmotic pressure** than what is contained within the cell.
- The RBCs will be shrink.
- Any concentration of NaCl that is **higher than 0.9%**, will be considered hypertonic for cells.

Solute outside the cell > Solute inside the cell





# Practical Part

---

# Objectives

---

1. To evaluate the effects of different solutions in blood hemolysis.
2. To detect the presence of blood in a biological sample.

## Calculate the concentration of an isotonic solution of a specific substance

---

For example, what is the concentration of NaCl (w/v%) that will make an isotonic solution for RBC? knowing that the osmolarity of RBC = 0.308 Osmolar.

**First:** Calculate the molarity from the osmolarity equation: [1]

Osmolarity = 0.308 Osmolar

No. of dissociation particles = 2, since  $\text{NaCl} \rightarrow \text{Na}^+ + \text{Cl}^-$

$$\rightarrow M = \frac{\text{Osmolarity}}{n} = \frac{0.308}{2} = 0.154 \text{ M}$$

$$[1] \text{ Osmolarity} = M \times n$$

**Where:**

M = molarity

n = No. of dissociation particles

**Pause and think** Why do you think it is important to prepare isotonic solutions?

# Calculations

---

To Calculate in w/v% expression:

**Second:** Calculate the No. of moles: [2]

→ No. of moles = M x V (in L) =

→ 0.154 M x 0.1 L = 0.0154 moles

$$[2] \text{ Molarity} = \frac{\text{No. of moles of solute}}{\text{Volume (L)}}$$

$$[3] \text{ No. of moles} = \frac{\text{Wt(g)}}{\text{Mwt}}$$

**Third:** Calculate weight in grams knowing that Mwt of NaCl = 58.5 g/mol: [3]

→ Wt (g) = No. of moles x Mwt

→ 0.0154 x 58.5 = 0.9 g in 100 ml then **0.9% w/v**

= 0.9 % → the concentration of NaCl that will make an isotonic solution

# Experiment (1): Hemolysis Test

## Method

1. Label 6 tubes (A → F). Then add as following:


	Tube A	Tube B	Tube C	Tube D	Tube E	Tube F
NaCl 0.45%	5 ml					
NaCl 1.2%		5 ml				
Sucrose 6%			5 ml			
NaOH 0.1 M				3 drops		
HCl 0.1 M					3 drops	
Dis. Water						5 ml
NaCl 0.9%				5 ml	5 ml	

2. Add 1ml of Sample A into each tube
3. Wait 10 - 30 min
4. Observe whether hemolysis has taken place



Pause and Think **What type of solution is distilled water considered?**

# Results



**A** Normal, non-hemolyzed sample

**B** Sedimented after one hour

**C** Hemolyzed sample

**Note:** the hemolyzed sample is transparent because there are no cells to scatter the light.


## Experiment (2): Detection of Blood by Benzidine Test

---

- It is often necessary to detect the presence of small quantities of blood in urine, stomach contents etc.

### Principle

- This method depends on the fact that the **heme group of hemoglobin** possesses a **peroxidase-like activity** which catalyzes the breakdown of hydrogen peroxide ( $\text{H}_2\text{O}_2$ ).
- The oxidizing species formed in this reaction can then react with benzidine giving **blue greenish color**.



**Note:** the test is not specific for blood as peroxidases present in milk, potatoes and pus, as well as the ions of  $\text{Fe}^{+3}$ ,  $\text{Cu}^{+2}$  and  $\text{K}^{+1}$  will give false positive results

# Experiment (2): Detection of Blood by Benzidine Test

---

## Method

### In a test tube add:

- Pinch of benzidine powder.
- 2ml glacial acetic acid.
- 3ml hydrogen peroxide ( $\text{H}_2\text{O}_2$ ).
- Add 1 – 2 ml of the sample B.

## Results

- If the test is **negative** → blood is absent from sample.
- If the test is **positive** → blood is probably not definitely present in sample.
- For this reason these tests are often described as “presumptive tests” .



Positive results