


### Exp.03: flavonoids extraction

#### Objectives:

- Extraction flavonoids from orange

#### Introduction:

Flavonoids are a class of plant secondary metabolites Consisting of 2 benzene rings linked via a heterocyclic pyrone or pyran ring they are benzo- $\gamma$ -pyrone derivatives so not open by alkalis.



#### The basic flavonoid skeleton can have a large number of substitutions on it:

- Hydroxyl groups
- Sugars-e.g. glucose, galactose, rhamnose.
- Sugars and hydroxyl groups: increase the water solubility of flavonoids
- Methyl and isopentyl groups: make flavonoids lipophile
- If no sugar: -AGLYCONE
- With sugar: -GLYCOSIDE

#### Classification of flavonoids:

According to the position of the linkage of the aromatic ring to the benzopyrano (chromano) moiety, it divided into three classes:

|   |  |
|---|--|
| 1- The flavonoids (2-phenylbenzopyrans) |  |
| 2- isoflavonoids (3-benzopyrans)        |  |
| 3- the neoflavonoids (4-benzopyrans)    |  |

|                    |  |
|--------------------|--|
| 4-Minor Flavonoids |  |
| 5-Anthocyanins     |  |
| 6-Bioflavonoids    |  |

All share a common chalcone precursor, and therefore are biogenetically and structurally related.

**Uses:**

**- For humane:**

Seem to have major health benefits for humans such as Antioxidant and other medicinal use

**- For plants:**

Protection from animal Breeding photosynthesis

**Tests of flavonoids:**

The extracts were dissolved in ethanol, filtered and subjected to following test.

**(a) Shinoda test:**

To the test solution, few drops of concentrated hydrochloric acid (HCL) were added. Then the magnesium turnings were put into the solution and observed for appearance of pink red color.

**(b) Zinc Hydrochloride reduction test**

To the test solution add a mixture of Zinc dust and conc. Hydrochloric acid. It gives red color after few minutes

**Experiment Procedure:**

| Step | Procedure   |
|------|---|
| 1    | - 10 g of orange in round- bottom flask   |
|      | - The flavonoids extraction from orangewill performed using soxhlet extractor (Fig. 1)  |
|      | - Add alcohol (Ethanol) in spherical flask.   |
|      | - Heat for 1 h  |
| 2    | Concentrated the solution.  |
| 3    | TLC using two system<br>- First system: Chloroform-methanol (95:5 v/v).<br>- Second system: Ethyl acetate- methanol-water (60:5:4 v/v/v).   |
| 4    | Notes the color of spot under UV light(brown).  |
| 5    | Use $\text{AlCl}_3$ as a reagent and notes the change in spot color under UV light (yellow).  |
| 6    | Calculate $R_f$ for flavonoid spot in two systems.  |
| 7    | <b>Performed:</b><br><b>(a) Shinoda test:</b><br>To the test Solution, add few fragments of Magnesium and concentrated Hydrochloric acid drop wise crimson red.<br><b>(b) Zinc Hydrochloride reduction test:</b><br>To the test solution add a mixture of Zinc dust and conc. Hydrochloric acid. It gives red color after few minutes |




Fig. 1