

**BCH 471**

**Hemolysing Agents  
& Detection of Blood**


# Blood Hemolysis

- **Hemolysis** (from the Greek Hemo: meaning blood, - lysis, meaning to break open)
- It is the breaking open of red blood cells and the release of hemoglobin and the red cell contents into the surrounding fluid (plasma)
- Hemolysis may occur *in vivo* or *in vitro*



# Hemolysis *in vivo*

- Conditions that can cause hemolysis include: Immune reactions, Infections, Medications. Toxins and poisons.
- Because the concentration of potassium inside red blood cells is much higher than in the plasma and so elevated potassium is usually found in biochemistry tests of hemolysed blood.



# Hemolysis *in vitro*

- Placing RBCs in a hypotonic solution
- Improper technique during collection (eg. incorrect needle size, excessive suction)
- pH imbalance (addition acid or base)

In this lab blood hemolysis will be done by using hypotonic solutions and pH imbalance.


# When Blood Hemolysis Should Be Done?

- Breaking down RBCs to release their content is often necessary for biochemistry
- Estimation of hemoglobin
- To obtain erythrocyte free preparation of leukocyte and platelet

# Osmosis and Osmotic Pressure


**Osmosis:** It is the diffusion of water across a selectivity preamble membrane into a region of higher solute conc. Once an equilibrium is reached the flow of water stops

**Osmotic pressure:** If the pressure that generate from osmosis



# Isotonic Solution


- A solution that has the same salt concentration as the normal cells of the body and the blood, having equal **osmotic pressure**.
- Example of Isotonic solution is **sodium chloride 0.9%**, have the same osmotic pressure as serum and they do not affect the membranes of the RBCs.
- In hospitals, intravenous fluids are isotonic.



# Hypotonic Solution

- In a hypotonic solution, there is a **lower concentration of solute outside a cell**, creating an environment with lower osmotic pressure than what is contained within the cell.

Solute inside the cell > Solute outside the cell



- The RBCs will burst or hemolyzed

- Any concentration of NaCl that is lower than 0.9%, will be considered hypotonic for cells




# Hypertonic Solution

- In a hypertonic solution, there is a **higher concentration of solute outside a cell**, creating an environment with higher osmotic pressure than what is contained within the cell.

Solute inside the cell < Solute outside the cell

- The RBCs will be shrink



- Any concentration of NaCl that is higher than 0.9%, will be considered hypotonic for cells

# How to Calculate the Concentration of an Isotonic Solution of a Specific Substance

For example you want to know the concentration of NaCl that will make an isotonic solution

Osmolarity of RBC = 0.308 Osmolar  $\rightarrow O = M \times \text{no. of dissociation particles}$

$\text{NaCl} \rightarrow \text{Na}^+ + \text{Cl}^-$  (no. of dissociation particles=2)

$$0.308 = M \times 2 \rightarrow M = 0.154 \text{ M}$$

**-To calculate in w/v % expression  $\rightarrow M = \text{mole}/V$  (in L)**

No. of moles =  $0.154 \times 0.1$  (100 ml, because you want it as %)

No. of moles = 0.0154 moles  $\rightarrow \text{No. of moles} = \text{Wt}/\text{Mwt}$

$$\text{Wt} = 0.0154 \times 58.5 = 0.9 \text{ g}$$

= 0.9 %  $\rightarrow$  the concentration of NaCl that will make an isotonic solution


# Method

	Tube A	Tube B	Tube C	Tube D	Tube E	Tube F
NaCl 0.45%	5 ml					
NaCl 1.2%		5 ml				
Sucrose 6%			5 ml			
NaOH 0.1M				3 drops		
HCl 0.1 M					3 drops	
Dis. Water						5 ml
NaCl 0.9%				5 ml	5 ml	

- Wait 30 min
- Observe wither hemolysis has taken place

**What type of solution is distilled water considered?**

# Results




Note that the hemolyzed sample is transparent, because there are no cells to scatter light.

# Detection of Blood by Benzidine Test

It is often necessary to detect the presence of small quantities of blood in urine, stomach contents etc.

## Principle

This method depends on the fact that the haem group of haemoglobin possesses a peroxidase-like activity which catalyses the breakdown of hydrogen peroxide ( $H_2O_2$ ).

The oxidising species formed in this reaction can then react with benzidine giving  blue greenish color

**However**, the test is not specific for blood as peroxidases present in milk, potatoes and pus, as well as the ions of  $Fe^{+3}$ ,  $Cu^{+2}$  and  $K^{+1}$  will give false positive results

**Method** → 3 ml of Sample + 2 ml Benzidine + 1 ml H<sub>2</sub>O<sub>2</sub>

## **Results**

- if the test is negative, blood is absent.

But

- if the test is positive, blood is probably, not definitely present.

→ For this reason the tests are often described as  
*"presumptive tests"*.