

Hemolysing Agents & Detection of Blood


BCH 471



Blood Hemolysis


Hemolysis (from the Greek Hemo: meaning blood, - lysis, meaning to break open)

- It is the breaking open of red blood cells and the release of hemoglobin and the red cell contents into the surrounding fluid (plasma)
- Hemolysis may occur *in vivo* or *in vitro*



Hemolysis *In Vivo*

- Conditions that can cause hemolysis include: Immune reactions, Infections, Medications, Toxins and poisons.
- Because the concentration of **potassium inside red blood cells is much higher** than in the plasma and so elevated potassium in plasma is usually found in biochemistry **tests of hemolysed blood**.



Hemolysis *In Vitro*

- Placing RBCs in a hypotonic solution.
- ⊙ Improper technique during collection (eg. incorrect needle size, excessive suction).
- ⊙ pH imbalance (addition acid or base).

In this lab blood hemolysis will be done by using hypotonic solutions and pH imbalance.


When Blood Hemolysis Should Be Done?

- Breaking down RBCs to release their content is often necessary for biochemistry
 - ⊙ Estimation of hemoglobin
 - ⊙ To obtain erythrocyte free preparation of leukocyte and platelet

Osmosis And Osmotic Pressure


Osmosis: It is the diffusion of water across a selectivity preamble membrane into a region of higher solute conc. Once an equilibrium is reached the flow of water stops

Osmotic pressure: Is the pressure that generate from osmosis.



Isotonic Solution


- A solution that has the **same salt concentration as the normal cells** of the body and the blood, having equal **osmotic pressure**.
- ⊙ Example of Isotonic solution is **sodium chloride 0.9% (Saline)**, have the same osmotic pressure as serum and they do not affect the membranes of the RBCs.
- ⊙ In hospitals, intravenous fluids are isotonic.



Hypotonic Solution

- In a hypotonic solution, there is a **lower concentration of solute outside a cell**, creating an environment **with lower osmotic pressure** than what is contained within the cell.
- The RBCs will burst or hemolyzed
- Any concentration of NaCl that is **lower than 0.9%**, will be considered hypotonic for cells


Solute inside the cell > Solute outside the cell



Hypertonic Solution

- In a hypertonic solution, there is a **higher concentration of solute outside a cell**, creating an environment with **higher osmotic pressure** than what is contained within the cell.
- The RBCs will be shrink
- Any concentration of NaCl that is higher than 0.9%, will be considered hypotonic for cells

Solute inside the cell < Solute outside the cell



How To Calculate The Concentration Of An Isotonic Solution Of A Specific Substance

For example you want to know the concentration of NaCl that will make an isotonic solution

Osmolarity of RBC = 0.308 Osmolar $\rightarrow O=M \times$ no.of dissociation partials

NaCl \rightarrow Na⁺ + Cl⁻ (no.of dissociation particles=2)

$$0.308 = M \times 2 \rightarrow M=0.154 \text{ M}$$

-To calculate in w/v % expression $\rightarrow M=\text{mole}/V$ (in L)

No.of moles = 0.154×0.1 (100 ml, because you want it as %)

No.of moles = 0.0154 moles \rightarrow No.of moles =

Wt/Mwt Wt= $0.0154 \times 58.5=0.9$ g

=0.9 % \rightarrow the concentration of NaCl that will make an isotonic solution

Method


- Label 6 Tubes (A → F). Then, Add 3 Drops Of Rbcs Suspended In Saline Into Each Tube

	Tube A	Tube B	Tube C	Tube D	Tube E	Tube F
NaCl 0.45%	5 ml					
NaCl 1.2%		5 ml				
Sucrose 6%			5 ml			
NaOH 0.1M				3 drops		
HCl 0.1 M					3 drops	
Dis. Water						5 ml
NaCl 0.9%				5 ml	5 ml	

- Wait 30 min
- Observe wither hemolysis has taken place

What type of solution is distilled water considered?

Results




Note that the hemolyzed sample is transparent, because there are no cells to scatter light.

Detection Of Blood By Benzidine Test

It is often necessary to detect the presence of small quantities of blood in urine, stomach contents etc.

Principle


This method depend on the fact that the haem group of haemoglobin possesses a peroxidase-like activity which catalyses the breakdown of hydrogen peroxide (H_2O_2).

The oxidising species formed in this reaction can then react with benzidine giving  blue greenish color

However, the test is not specific for blood as peroxidases present in milk, potatoes and pus, as well as the ions of Fe^{+3} , Cu^{+2} and K^{+1} will give false positive results

Method

- 3 ml of Sample heated in a boiling water bath for 3min
- 2 ml Benzidine + 1 ml H_2O_2



Results

- if the test is negative, blood is absent.

But

- if the test is positive, blood is probably, not definitely present.

→ For this reason the tests are often described as *"presumptive tests"*.