


# Quantitative estimation of protein in urine

By sulphosalicylic acid Method

BCH 472




- **In a healthy renal and urinary tract system**, the urine contains no protein or only trace amounts.
- **The presence of increased amounts of protein** in the urine can be an important indicator of **renal disease**. It may be the first sign of a serious problem and may appear before any other clinical symptoms.
- However, there are other physiologic conditions (eg, exercise, fever) that can lead to increased protein excretion in urine. Also, there are some renal disorders in which proteinuria is absent.



# Proteinuria:


- Protein in **normal** urine should be less than **150 mg/L**
- **Proteinuria** is defined as urinary protein excretion of greater than **150 mg** per day(per one liter).
- **Note: Dipsticks** (is the most common initial screening test for proteinuria) can only detect around 150 mg/L of albumin(The dipstick will not detect nonalbumin proteins).
- **Note:** /L = /24-hour = /day.



# Types of Proteinuria:

Type	Cause
Glomerular proteinuria	<ul style="list-style-type: none"><li>• Results from a <b>disruption of the glomerular filtration barrier</b> which <u>increased filtration of normal plasma protein</u> and because albumin has the highest concentration in the plasma it is called albuminuria eg. <u>Malignant hypertension</u></li></ul>
Tubular proteinuria	<ul style="list-style-type: none"><li>• <b>Defect in the reabsorption</b> eg, <u>Fanconi Syndrom</u> - low molecular weight protein that is found in urine</li></ul>
Overflow proteinuria	<ul style="list-style-type: none"><li>• <b>Overflow of high plasma</b> - high concentrations of low molecular weight protein found in urine eg, <u>Multiple myeloma</u>. In multiple myeloma excessive amounts of immunoglobulin light chains are produced.</li></ul>
Secretory proteinuria	<ul style="list-style-type: none"><li>• <b>Over secretion</b> of certain proteins in the tubules, most notably the over secretion of Tamm-Horsfall proteins eg, in <u>interstitial nephritis</u></li></ul>

# Types of Proteinuria:



# Cause of Proteinuria as Related to Quantity

Protein amount per 24-hour	Type of proteinuria
0.15 to 2.0 g	<ul style="list-style-type: none"><li>• Tubular proteinuria</li><li>• Overflow proteinuria (an increased proportion of low molecular weight proteins)</li></ul>
2.0 to 4.0 g	<ul style="list-style-type: none"><li>• Usually glomerular</li></ul>
➤ > 4.0 g	<ul style="list-style-type: none"><li>• Always glomerular (mainly albumin)</li></ul>

- **The quantitative estimation** of the daily excretion of protein is of value to the clinician in order to give a general idea of the type of renal disease, its severity and to monitor the results of treatment given.
  - The protein content can be determined by numerous methods eg, Biuret, Lowry, Bradford.
  - In this lab **turbidimetric method** will be used.
  - Determination of total protein by measurement of protein turbidity produce by mixed with an anionic organic acid such as sulfosalicylic acid , TCA , or benzethonium chloride.
  - Sulphosalicylic acid is used in this experiment to precipitate the protein in *a 24 hour* sample of urine. The turbidity is proportional to the concentration of the protein, and may be measured with a spectrophotometer at **500 nm**.
-



- **Practical Part**

---



# Sulfosalicylic acid (SSA) test :

- The sulfosalicylic acid (SSA) turbidity test quantitatively screens for proteinuria. The advantage of this **easily** performed test is its greater **sensitivity** for proteins such as Bence Jones.
- The SSA reaction will detect globulin and Bence-Jones proteins, in addition to albumin (although it is more sensitive to albumin).



# Principle:

- Sulfosalicylic acid is an **anion(-)** which neutralizes the protein **cations(+)** leading to its precipitation (pH in highly **acidic** media, the protein will be positively charged, which is attracted to the acid anions that cause them to precipitate.)
- Then the radiation of a wavelength which is not absorbed by the solution is made to pass through the suspension and the apparent absorption will be solely because of the **scattering** by the particles.
- The transmission decrease with increasing protein concentration.



# Method:

1-Set up a series of test tube as follows, label from 1- 7

Tube	Protein Stock Solution( 140 mg/dl)	0.85% Saline	Protein concentration mg/dl
1	4.5	1.5	
2	3	3	
3	2.4	3.6	
4	1.5	4.5	
5	0.9	5.1	
6	0.3	5.7	
7( Blank)	0	6	
Urine Sample	-	-	

2-Set another 8 test tube labeled 1-7 and pipette in each one Add 8 ml of sulfosalicylic acid


Tube	sulfosalicylic acid
1	8 ml
2	8 ml
3	8 ml
4	8 ml
5	8 ml
6	8 ml
7( Blank)	8 ml
Urine Sample	8 ml

- 3-Into tube 1 pipette 2 ml of protein solution 1, into tube 2 pipette 2 ml of protein solution 2 etc. For the Urine Sample pipette 2ml of the Sample.
- 4-**Mix** contents of each tube well and allow standing for 5 minutes.
- 5-Using solution 7 (Blank) to set transmittance at 100 at 500nm.
- 6-Then use solutions from 1-6, to recorded respective **transmittance** of each suspension.

## Results:

Tube	Transmittance at 500 nm	Protein concentration mg/dl
7( Blank)	100 %	
1		
2		
3		
4		
5		
6		
Urine Sample		

- Plot Transmittance against Protein concentration mg/dl.
- Determine the Protein concentration of Urine Sample from the standard curve.
- Compare the result you got with the normal range of protein excretion in 24 h urine specimen (you must convert the unit to g/L) if you know that the protein excretion in **healthy sample (0- less than 0.150g/24 h)**.
- Comment on the clinical conditions of the patient if it is present.
- Assuming that the **24 hour urine sample for the patient = 1000 ml**.



# Questions :

Which protein can be used as a marker of Multiple myeloma ? And how?

Explain how total protein can be determined by measuring protein turbidity ?

Why the resulting graph is a descending curve ?

---