




Preparation Of Buffer Solutions By Different Laboratory Ways



Dissociation of Triprotic acid

- **Triprotic acid** is acid that contain three hydrogens ions.
- It dissociates in solution in three steps, with three Ka values.
- **phosphoric acid** is an example of triprotic acid .
- It dissociates in solution as following:



preparation of Buffer by several ways:

For example if you asked to prepare sodium phosphate buffer [$\text{NaH}_2\text{PO}_4 / \text{Na}_2\text{HPO}_4$]

You can prepare it by:

1. By mixing NaH_2PO_4 (conjugate acid) and Na_2HPO_4 (conjugate base) in the proper proportions ,
2. By starting with H_3PO_4 and converting it to NaH_2PO_4 plus Na_2HPO_4 by adding the proper amount of **NaOH**.
3. By starting with NaH_2PO_4 and converting a portion of it to Na_2HPO_4 by adding **NaOH**.
4. By starting with Na_2HPO_4 and converting a portion of it to NaH_2PO_4 by adding a strong acid such as **HCL**.
5. By starting with Na_3PO_4 and converting it to Na_2HPO_4 plus NaH_2PO_4 by adding **HCL**.
6. By mixing Na_3PO_4 and NaH_2PO_4 in the proper proportions.

HCl
'donate H⁺'



NaOH
'accept H⁺'


Example: Prepare 0.1 liters of 0.045 M sodium phosphate buffer, pH=7.5, [pka₁= 2.12, pka₂ = 7.21 and pka₃ = 12.30]

- a) From concentrated (15M) H₃PO₄ and solution of 1.5 M NaOH .
- b) From solid NaH₂PO₄ and solid NaOH.

Calculations:


1st, write the equations of Dissociation of phosphoric acid and the pka of corresponding ones:

Because phosphoric acid [H₃PO₄] has (Triprotenation : it has 3 dissociation phases) so,



Regardless of which method is used, the first step involves calculating number of moles and amounts of the two ionic species in the buffer

2nd, choose the pka value which is near the pH value of the required buffer, to be able to know the ionic species involved in your buffer:



The pH of the required buffer [pH =7.5] is near the value of **pka₂**, consequently, the two major ionic species present are **H₂PO₄⁻ (conjugate acid)** and **HPO₄²⁻ (conjugate base)**. with the HPO₄²⁻ predominating { since the pH of the buffer is slightly basic }

3rd, calculate No. of moles for the two ionic species in the buffer:

$\text{pH} = \text{pKa}_2 + \log \left[\frac{\text{HPO}_4^{2-}}{\text{H}_2\text{PO}_4^-} \right]$ → Note that : $[\text{A}^-] = \text{HPO}_4^{2-}$, $[\text{HA}] = \text{H}_2\text{PO}_4^-$

• Since the buffer concentration is **0.045M**, so assume $[\text{A}^-] = y$, $[\text{HA}] = 0.045 - y$:

$$7.5 = 7.2 + \log (y / 0.045 - y)$$

$$7.5 - 7.2 = \log (y / 0.045 - y)$$

$$0.3 = \log(y / 0.045 - y) \rightarrow \text{antilog for both sides}$$

$$2 = (y / 0.045 - y) \rightarrow y = 0.09 - 2y \rightarrow 3y = 0.09 \rightarrow y = 0.09/3 = \mathbf{0.03M} \rightarrow \text{conc. of } [\text{HPO}_4^{2-}] = [\text{A}^-] = y$$



$$\text{So, conc. of } [\text{H}_2\text{PO}_4^-] = [\text{HA}] = 0.045 - y = 0.045 - 0.03 = \mathbf{0.015 M}$$

• **Now found the number of mole for the two ionic species in the buffer:**

$$\text{No. of moles of } \text{HPO}_4^{2-} (\text{A}^-) = M \times V = 0.03 \times 0.1 = \mathbf{0.003 \text{ moles.}}$$

$$\text{No. of moles of } \text{H}_2\text{PO}_4^- (\text{HA}) = M \times V = 0.015 \times 0.1 = \mathbf{0.0015 \text{ moles}}$$

Note :




Note: Total no. of moles of phosphate buffer = $M \times V = 0.045 \times 0.1 = \mathbf{0.0045 \text{ moles.}}$

Now, to prepare the required buffer:

a) From concentrated (15M) H_3PO_4 and solution of 1.5 M NaOH .

Remember that the two ionic species involved in the buffer are:



Calculations:


Start with **0.0045 mole** of H_3PO_4 and add **0.0045 moles** of NaOH to convert H_3PO_4 completely to H_2PO_4^- (HA) , then add **0.003 moles** of NaOH to convert H_2PO_4^- to give HPO_4^{2-} (A^-):

No. of moles needed of NaOH= $0.0045+0.003=$ **0.0075 moles**

Volume of NaOH needed(L)= no.of moles (of NaOH) / M (of NaOH) = $0.0075/ 1.5 = 0.005 \text{ L} =$ **5 ml**


Volume of H_3PO_4 needed(L) =no.of moles (of H_3PO_4) / M (of H_3PO_4) = $0.0045/ 15 =0.0003 \text{ L} =$ **0.3 ml**

→ Add **5ml** of NaOH to the **0.3 ml** of concentrate H_3PO_4 , mix ; then add sufficient water to bring the final volume to 0.1 liters (100 ml), and **check the pH**.




b) From solid NaH_2PO_4 and solid NaOH .

Remember that the two ionic species involved in the buffer are:



□ Calculations

- Start with **0.0045 mole** of NaH_2PO_4 (**HA**) and add **0.003 moles** of NaOH to convert NaH_2PO_4 to give Na_2HPO_4 (**A⁻**):



- Weight in grams of NaH_2PO_4 needed = no. of moles \times MW = $0.0045 \times 119.98 = 0.54 \text{ g}$
- Weight in grams of NaOH needed = no. of moles \times MW = $0.003 \times 40 = 0.12 \text{ g}$
- Dissolve the **0.54 g** of NaH_2PO_4 and **0.12 g** of NaOH in some water, mix ; then add sufficient water to bring the final volume to 0.1 liters (100 ml), and **check the pH**

Because we start with **solid**, We convert moles to **grams** not volume



Practical Part

Objective:

- To learn how to prepare a buffer by different laboratory ways.

Method:

Prepare 0.1 liters of 0.045 M sodium phosphate buffer, pH=7.5,
[pKa1= 2.12, pKa2 = 7.21 and pKa3 = 12.30]:

a) From concentrated (15M) H₃PO₄ and solution of 1.5 M NaOH:

Add **5ml** of NaOH to the **0.3 ml** of concentrate H₃PO₄, mix ; then add sufficient water to bring the final volume to 0.1 liters(100 ml), and check the pH.

b) From solid NaH₂PO₄ and solid NaOH:

Dissolve the **0.54g** of NaH₂PO₄ and **0.12g** of NaOH in some water, mix ; then add sufficient water to bring the final volume to 0.1 liters(100 ml), and check the pH.

H.W

1) Prepare 100 ml of 0.045 M sodium phosphate buffer, pH=7.5, [pka1= 2.12, pka2 = 7.21 and pka3 = 12.30]: (MW of Na₂HPO₄= 142, MW of Na₃PO₄= 164)

c) You are provided with solid Na₂HPO₄ and 2M solution of HCl.

d) You are provided with solid Na₃PO₄ and 2 M HCL .

2) Prepare 500 ml of 0.3 M sodium phosphate buffer, pH=2.5, [pka1= 2.12, pka2 = 7.21 and pka3 = 12.30]:

a) You are provided with solid Na₂HPO₄ and 5M (HCl or NaOH).

b) You are provided with 2M H₃PO₄ and 5M (HCL or NaOH) .