

Titration curve of amino acids




Titration curve

- Titration Curves are produced by monitoring the pH of a given volume of a sample solution after successive **addition of acid or alkali**.
- The curves are usually plots of pH against the volume of titrant added (acid or base).
- Each **dissociation group** represent **one stage** in the titration curve.

Amino acid general formula:

Amino acids consist of:

- ❑ A basic amino group ($-\text{NH}_2$)
- ❑ An acidic carboxyl group ($-\text{COOH}$)
- ❑ A hydrogen atom ($-\text{H}$)
- ❑ A distinctive side chain ($-\text{R}$).



Titration of amino acid:

- When an amino acid is dissolved in water it exists predominantly in the **isoelectric form**.
- Amino acid is an **amphoteric** compound → It act as either an acid or a base:
 - **Upon titration with acid** → it acts as a **BASE** (accept a proton).
 - **Upon titration with base** → it acts as an **ACID** (donate a proton)


□ Amino acids are example of **weak acid** which contain **more than one dissociate group**.

□ **Examples:**

(1) Alanine:

-Contain COOH ($pK_{a1}= 2.34$) and NH_3^+ ($pK_{a2}= 9.69$) groups (it has one pI value =6.010). **[Diprotic]**


-The COOH will dissociate first then NH_3^+ dissociate later . (Because $pK_{a1}<pK_{a2}$)




Full protonated alanine

(2) Arginine:

-Contain COOH ($pK_{a1}= 2.34$) , NH_3^+ ($pK_{a2}= 9.69$) groups and basic group ($pK_{a3}=12.5$) (it has one pI value=11). **[Triprotic]**



Titration curve of Alanine



pK_1 carboxylic acid = 2.34
 pK_2 amino group = 9.69
 $pI = (pK_1 + pK_2) / 2$

Titration curve of alanine or glycine [diprotic]:

[1] In starting point:


- Alanine is full protonated.
- $[\text{NH}_3^+-\text{CH}-\text{CH}_3-\text{COOH}]$.

[2] COOH will dissociate first:

- $[\text{NH}_3^+-\text{CH}-\text{CH}_3-\text{COOH}] > [\text{NH}_3^+-\text{CH}-\text{CH}_3-\text{COO}^-]$
- $\text{pH} < \text{pK}_{a1}$.


[3] In this point the component of alanine act as buffer:

- $[\text{NH}_3^+-\text{CH}-\text{CH}_3-\text{COOH}] = [\text{NH}_3^+-\text{CH}-\text{CH}_3-\text{COO}^-]$.
- $\text{pH} = \text{pK}_{a1}$



Cont.

[4] In this point:



$$\text{pH} > \text{pKa}_1.$$

[5] Isoelectric point:

The COOH is full dissociate to COO⁻.




Con. of -ve charge = Con. of +ve charge.

The amino acid present as **Zwitter ion** (neutral form) .


Remember that :PI (isoelectric point) is the pH value at which the net charge of amino acid equal to zero.

$$\text{pI} = (\text{pKa}_1 + \text{pKa}_2) / 2 = (2.32 + 9.96) / 2 = 6.01$$

[6] The NH₃⁺ start dissociate:




$$\text{pH} < \text{pKa}_2.$$




Cont.

[7] In this point the component of alanine act as **buffer**:



$$\text{pH} = \text{pK}_{a2}.$$

[8] In this point:



$$\text{pH} > \text{pK}_{a2}$$


[9] End point:

The alanine is full dissociated.



$$\text{pOH} = (\text{pK}_b + \text{P}[\text{A}^-]) / 2$$

$$\rightarrow \text{pK}_b = \text{pK}_w - \text{pK}_{a2}$$



Calculating the pH at different point of the titration curve :

The pH calculated by different way :

[1] at starting point :

$$\text{pH} = (\text{pK}_a + \text{P}[\text{HA}]) / 2$$

[2] At any point within the curve (before or in or after middle titration):

$$\text{pH} = \text{pK}_a + \log([\text{A}^-] / [\text{HA}])$$

[3] At end point:

$$\text{pOH} = (\text{pK}_b + \text{P}[\text{A}^-]) / 2$$

$$\text{pH} = \text{pK}_w - \text{pOH}$$

$$\text{pK}_b = \text{pK}_w - \text{pK}_a$$



Practical Part

Objectives

- To study titration curves of amino acid.
- To use this curve to estimate the pK_a values of the ionizable groups of the amino acid.
- To determine pI.
- To determine the buffering region.
- To understand the acid base behaviour of an amino acid.

Method:

- a) You are provided with 10 ml of a 0.1M alanine solution, titrate it with 0.1M NaOH adding the base drop wise mixing, and recording the pH after each 0.5 ml NaOH added until you reach a pH=11.


Measured pH value	Amount of 0.1M NaOH added [ml]

- a) Take another 10 ml of a 0.1M alanine solution, titrate it with 0.1 M HCL adding the acid drop wise mixing, and recording the pH after each 0.5 ml HCL added until you reach a pH=2.17.


Measured pH value	Amount of 0.1M HCl added [ml]

Results


- Record the titration table and plot a curve of pH versus ml of titrant added.
- Calculate the pH of the alanine solution after the addition of 0 ml, 5ml, of 0.1M NaOH, and calculate pH after addition of 0.5 ml , 2ml of HCl.
- Determine the pKa of ionizable groups of amino acids.
- Determine the PI value from your result
- Compare alanine pka and pI values with those obtained from Curve.
- Compare your calculated pH values with those obtained from Curve.
- Determine the buffering points.




Titration curve of alanine with 0.1M NaOH




Titration curve of alanine with 0.1M HCl





- Note: in calculating the pH:
- At any point within the curve $\text{pH} = \text{pK}_a + \log\left(\frac{[\text{A}^-]}{[\text{HA}]}\right)$



If a base is added:

The amino acid will be treated as an **acid**

The pK_a used is of the **amine group**.

The upper stage

If acid is added:

The amino acid will be treated as a **base**

The pK_a used is of the **carboxyl group**

The lower stage