

# Determination of food acidity

BCH 445



Lab 1



# Food acidity

Food acids are usually **organic acids**, with citric, malic, lactic, tartaric, and acetic acids being the most common.

However, **inorganic acids** such as phosphoric and carbonic acids (arising from carbon dioxide in solution) often play an important and even predominant role in food acidulation.



## The organic acids present in foods influence :

- **flavor** (i.e., tartness)
- **Color** (through their impact on anthocyanin and other pH-influenced pigments)
- prevent/retard the growth of **microorganisms** or inhibit the germination of spores
- Providing the proper environment for **metal ion chelation**, an important phenomenon in the minimization of lipid oxidation

## Organic acids may present :

- Naturally,
- By Fermentation
- Added as part of a specific food formulation



# The importance of determining food acidity

## 1. Determine the degree of maturity of fruits and vegetables

The titratable acidity of fruits is used, along with sugar content, as an indicator of **maturity**, generally

→ the higher the maturity, the **lower** the acid content.

e.g. in the ripening process, such as tomatoes from green to mature stage, there is an **increase** in sugar content.

## 2. To determine the freshness of foods

for example in milk, the **more** the lactic acid levels, means that milk is **rotten**.



### **3. Acidity indicators reflect the quality of food**


The amount of organic acids in food directly affects the food flavor, color, stability, and the level of quality.

### **4- Determination of acid on the microbial fermentation process**

Such as: fermentation products in soy sauce, vinegar and other acids is an important indicator of quality.



There are two ways to express food acidity:



**Titratable acidity**

- Simple estimate of the total acid content of food
- Better predictor of acid impact on flavor

**Hydrogen concentration pH**


- Depend on the strength of acid condition

# Titratable acidity

- **Titratable acidity** provide a simple estimate of acid in food, it is a **routine titration** that cannot differentiate between individual acids. Therefore, titratable acidity is usually stated in terms of **predominant acid**
- It is determined by **neutralizing** the acid present in a known quantity (weight or volume) of food sample using a standard base.
- **The endpoint for titration is** determined usually by the color change of a pH-sensitive dye, typically phenolphthalein.
- It is calculated by the following formula:

$$\text{TA}\% = (\text{wt of acid} / \text{wt of sample}) \times 100$$

▫ **Wt of acid is calculated** as no. of moles of NaOH that neutralizes the acid x mwt of acid



Note the color at end point:





# Objective

To determine total acidity of milk, juice, vinegar and oil acid value.



# 1-Determination of Milk Acidity:

Measuring milk acidity is an important test used to determine milk **quality**.

**The Natural acidity of fresh milk is due to :**

phosphates, casein and whey proteins, citrates and carbon dioxide dissolved during the process of milking.

**Developed acidity is due to:**

lactic acid produced by the action of bacteria on lactose in milk.





# Milk Acidity

**TA%**=0.12% - 0.16% the average 0.14%

If it **increased** more than 0.16% , is an indication of **lactic acid** by bacteria.

**Acidity is expressed as percentage of lactic acid** because lactic acid is the principal acid produced by fermentation.



## Method:

1. Mix the milk sample thoroughly by avoiding incorporation of air.
2. Transfer 10 ml milk to conical flask or beaker .
3. Add equal quantity of distilled water .
4. Add 2 drops of phenolphthalein indicator and stir.
5. Rapidly titrate the contents with 0.1 N NaOH solution, continue to add alkali drop by the drop and stirring the content till first definite **change to pink colour** .
6. **Note down the final burette reading.**

## Result and Calculation:

$$\text{Lactic acid \%} = \frac{(0.1\text{M NaOH} \times \text{vol. of NaOH (in liter)} \times \underline{90.08})}{\text{Weight of the sample}} \times 100$$

**90.08 g/ mol** is the molecular weight of Lactate.


**Normal range** = 0.12% - 0.16%

## 2-Determination of total acidity in juice :


- The acidity of natural fruit juices is the result mainly of their content of **organic acids**.

▫ For example,

▫ most fruits contain the **tricarboxylic acid** (citric acid)



▫ **whereas** grapes are rich in tartaric acid & peaches, apricots and plums in malic acids, both are **dicarboxylic acids**.



## Method:

- 1- Weight 10 gm juice in beaker.
- 2- Add 25 ml of distilled water.
- 3- Titrate with 0.1M NaOH , using 2 drops of **phenolphthalein** as an indicator.

# Calculations:

Calculate percent acidity of fruit juice ( citric acid):


$$1- \text{Wt. of citric acid} = \frac{0.1\text{M NaOH} \times \text{vol. of NaOH (in liter)} \times 192.43}{3}$$

\*192.43 g/mol is the molecular weight of citric acid

$$2- \% \text{ of total acidity} = (\text{wt. of acid} / \text{wt. of sample}) \times 100$$

Normal range for citric acid = 0.39 - 1.1 %


\*Why when calculation of the weight of citric acid it divided by 3?





### 3-Determination of total acidity in vinegars:

- The acidity of vinegars is derived by the fermentation of ethanol by acetic acid bacteria which produce **acetic acid**.
- It may be determined titrimetrically using **phenolphthalein** as an indicator .
- The **natural acidity** of vinegar is mainly due to the presence of **acetic acid (CH<sub>3</sub>COOH)** , which is volatile .



# Method:

## **Determination of total acidity**

1- Weight 1 gm vinegar.

2- Add 10 ml of distilled water.

3- Titrate with 0.1M NaOH , using 2 drops of **phenolphthalein** as an indicator.

## Calculations:

Calculate percent acidity as acetic acid (MW=60.05)

1- **Wt. of acetic acid**= (0.1M NaOH X volume of NaOH in liter X MW)

2- **% of total acidity**= (wt. of acid / wt. of sample) X 100


**Normal range**= 4-6 %

## 4-Acid value :

**The acid value** is defined as the number of milligrams of sodium hydroxide required to neutralize the **free fatty acids** present in one gram of fat.

It is a relative measure of rancidity as free fatty acids are normally formed during decomposition of oil glycerides.

The value is also expressed as percent of free fatty acids calculated as **oleic acid (main fatty acid in olive oil)**.



## Principle:

- The value is a measure of the **amount of fatty acids** which have been liberated by hydrolysis from the glycerides due to the action of **moisture, temperature and/or lipolytic enzyme lipase**.
- The acid value is determined by **directly titrating** the oil/fat in an alcoholic medium against standard sodium hydroxide solution.

## Method:

1. Mix the oil or melted fat thoroughly before weighting.
2. Weight accurately about 5 g of cooled oil sample in a 250 ml conical flask.
3. Add 50 ml of freshly neutralized hot ethanol.
4. Add one ml of phenolphthalein indicator solution.
5. Boil the mixture( in water bath) for about 5 minutes and titrate while hot against standard alkali solution shaking vigorously during the titration.

## Calculation:

**Acid value** =  $40 \times (V \times N) / \text{weight of sample}$

\*40 g/mol is the molecular weight of NaOH

Where **V** = Volume in **ml** of standard potassium hydroxide or sodium hydroxide used

**N** = Normality of the Sodium hydroxide solution = **0.1 N** .

**W** = Weight in **g** of the sample

→ The **maximum level** allowed for **acid value** of edible fats and oils is **0.6 mg NaOH/g**

# Discussion

Discuss the result you got for each sample and compare it to the normal range