



# **Estimation of Serum Creatinine, Urine Creatinine , and Creatinine Clearance**

BCH 472




# Kidney functions:

- **The kidneys serve three essential functions:**
  1. They function as filters, removing metabolic products and toxins from the blood and excreting them through the urine.
  2. They regulate the body's fluid status, electrolyte balance, and acid-base balance.
  3. The kidneys produce or activate hormones that are involved in erythropoiesis,  $\text{Ca}^{2+}$  metabolism, and the regulation of blood pressure and blood flow.





- **Renal function tests** are used to detect the presence of renal diseases and assess their progress.
- The most widely used test is to measure the **glomerular filtration rate (GFR)**, that is, the rate of filtrate formation by the kidneys.




# Glomerular Filtration Rate:

- Under **normal** conditions, approximately 625 mL of plasma flow through the kidneys each minute and the volume of plasma **filtered** is 125 mL/ min which is called the glomerular filtration rate.
- **Glomerular filtration rate(GFR)**, is the volume of plasma filtered by the kidneys in per unit of time.
- **GFR** is an important and the best overall measurement in the evaluation of kidney function.



# Measuring the GFR:


- Accurate measurement of the GFR by **clearance tests** requires determination of the concentration, in plasma and urine, of a substance is known to be completely filtered from the plasma at the glomerulus.
- This substance must not be reabsorbed nor secreted by renal tubules, broken down, or accumulated by the tubules and must remain at a constant concentration in the plasma throughout the period of urine collection.
- It's clearance is given by ,**clearance = U.V/ P**
- ( U= urine creatinine, p= plasma creatinine, v= volume of urine)



# Substances clearance used for Measuring GFR



	Inulin Clearance	Creatinine Clearance	Urea Clearance
Source	Non-toxic fructose polymer	End-product of skeletal muscle creatine metabolism	endproduct of protein Metabolism
Advantages	Not reabsorbed or secreted	An endogenous product of muscle metabolism; near constant production	An endogenous product of protein
Disadvantages	Not made by body; must be injected	Small amount is secreted	Partially reabsorbed synthesis varies with diet

\* Creatinine clearance is preferred because it is a normal constituent of blood and **no infusion** is needed unlike inulin. Moreover it is not reabsorbed by the tubules as in the case of urea.



# Creatinine:

- In the muscles “**creatine**” is converted to creatine phosphate which becomes the source of a high energy phosphate bond for the immediate reformation of ATP.
- “**Creatinine**” is the byproduct of muscle energy metabolism and is produced at a constant rate according to the muscle mass of the individual. It is a substance that, in health, is easily excreted by the kidney.
- Endogenous creatinine production is **constant** as long as the muscle mass remains constant



# Serum Creatinine:

## High plasma creatinine


Plasma creatinine tends to be higher in subjects with a large muscle mass.

**Other non-renal causes of increased plasma creatinine include the following :**

- A high meat intake can cause a temporary increase.
- Transient, small increases may occur after vigorous exercise .



\* If non-renal cause can be excluded, an **increased** plasma creatinine indicates a fall in GFR





# Urine Creatinine

## Decreased urine creatinine is found in:

- There are many cases leads to a decrease in urine creatinine eg. Advanced renal disease, renal stenosis.

## Increased urine creatinine is found in:

There are many cases eg. Diabetes mellitus



# Creatinine clearance :

A measure of the amount of creatinine **eliminated** from the blood by the kidneys.

- Creatinine is cleared from the body fluids almost **entirely by glomerular filtration** (small amount is secreted by kidney tubules).
  - Therefore, the clearance of creatinine can be used to assess **GFR**.
  - Because measurement of creatinine clearance does **not** require intravenous infusion into the patient, this method is much more **widely** used than inulin clearance for estimating GFR clinically.
  - Tubules, to variable degree, **secrete creatinine**, which, by itself, would lead to an ~20% overestimate of GFR in humans.
-

# Clinical Implications:

1. **Decreased creatinine clearance** is found in any condition that decreases renal blood flow:

- Impaired kidney function.
- Shock, dehydration.
- Hemorrhage.

2. **Increased creatinine clearance** is found in:

- Pregnancy
-

# Reference Values:

- Urine creatinine :1- 2 g/ 24h
- (serum) creatinine: 0.6–1.2 mg/dL
- Normal creatinine clearance= 100-130 ml/min/1.73m<sup>2</sup>

**Note:** What 1.73 m<sup>2</sup> means?

Kidney function is proportional to kidney size, which is proportional to body surface area. A body surface area of 1.73 m<sup>2</sup> is the normal mean value for young adults.

Adjustment for body surface area is necessary when comparing a patient's estimated GFR to normal values or to the levels defining the stages of CKD.

---

# Practical Part

## Experiments

```
graph TD; A[Experiments] --- B[1-Estimation of Serum Creatinine]; A --- C[2-Estimation of Urine Creatinine]; A --- D[3-Calculation of Creatinine Clearance];
```

1-Estimation of Serum Creatinine

2-Estimation of Urine Creatinine

3-Calculation of Creatinine  
Clearance

---



# Objective:

- 1- To estimate creatinine in serum and urine.
- 2- To calculate creatinine clearance value.

# Principle:

## (Jaffe's method):

Colorimetric estimation of creatinine using the alkaline picrate method.





- 2- Immerse the Tubes carefully in the boiling water bath for 40 seconds.
- 4- Pipette 0.6 ml of NaOH to all tube
- 5- Let the tubes stand for 20 min.
- 6- Read the absorbance at **520 nm**.

## Results:

Tube	Standard (serum)		Test (serum)		Test (urine)		Standard(Urine)	
	(A)	(B)	(C )	(D)	(E)	(F)	(G )	(H)
Absorbance at 520 nm								
Average (Mean of Absorbance)								



# Calculation:

## 1-Serum creatinine =

(Mean Absorbance of sample serum ÷ Mean Absorbance of Standard serum) X  
concentration of standard serum(3 mg/dl) = ..... mg / dl

## 2-Urine creatinine =

(Mean Absorbance of sample urine ÷ Mean Absorbance of Standard) X  
concentration of standard urine(0.75 mg/dl) X **DF** (100) = ..... mg / dl

To compare with normal range, convert from mg/dl to g/24 h

## 3- Creatinine Clearance : = U.V/ P

$$= \frac{\text{Urinary creatinine (mg/dl)}}{\text{plasma creatinine (mg/dl)}} \times \text{Urine volume(ml/min)} = A$$

**Note:** 24h urine volume = 100ml, body surface= 1.6m<sup>2</sup>

- A-----→ 1.6 m<sup>2</sup>
  - ?-----→ 1.73 m<sup>2</sup>
  - **-Corrected for surface area=                      ml/min/1.73 m<sup>2</sup>**
-

# Creatinine Clearance:

Find the Creatinine Clearance = if you know that the Urine creatinine  $U = 488$  mg/dl , Serum creatinine  $P = 2.32$  mg/dl, Volume of urine in 24 h  $V = 100$  ml and  $A$  (surface area) =  $1.6 \text{ m}^2$

**A- Creatinine Clearance:** =  $U.V / P$

For example:  $V = 100 \text{ ml} / 24 \text{ h}$  >> to convert hours to minutes /1440 “24x60”

$(488 \text{ mg/dl} \div 2.32 \text{ mg/dl}) \times (100 \div 1440) = 14.6 \text{ ml/min}$

14.6 ml/ min in  $1.6 \text{ m}^2$ , find Creatinine clearance in  $1.73 \text{ m}^2$  :

$= (14.6 \times 1.73) \div 1.6 = 15.8 \text{ ml/min} / 1.73 \text{ m}^2$

**B- Creatinine Clearance:** =  $(U \times V \times 1.73) / (P \times 1440 \times A)$

$(488 \text{ mg/dl} \times 100 \times 1.73) / (2.32 \times 1440 \times 1.6)$

$= 15.8 \text{ ml/min} / 1.73 \text{ m}^2$

---

## Discussion:

- Comment on the concentration of **creatinine in serum**.
- Comment on the concentration of **creatinine in urine**.
- Comment on the value of **Creatinine Clearance** .



## Question-Home work:

A man aged 35 years has a serum creatinine of 3 mg/dl. A 24 h urine of 2160 ml is collected and found to a creatinine concentration of 400 mg/dl

Calculate the Creatinine Clearance

The logo consists of the letters 'H' and 'W' in a bold, serif font. The 'H' is black and the 'W' is red, with the 'W' overlapping the 'H'.