Exp#6 Bilirubin

Total and Direct bilirubin

Daheeya Alenazi (Cromatest kit (CHS 264

Definition

Bilirubin: is the yellow breakdown product of .normal haemoglobin & excreted in bile and urine

(DIRECT (Conjugated	(INDIRECT(Unconjugated
bound to albumin and sent to the liver. In the liver it is conjugated with glucuronic acid soluble in water.	Erythrocytes damaged Hb The heme is then turned into unconjugated bilirubin in the reticulo-endothelial cells of .the spleen .not soluble in water

Clinical Significance

.An increase in bilirubin conc. Is called jaundice-

- High level of direct type means that the bile is not being .properly excreted
- .High levels of indirect mean more Hb is being damaged-

Hyberbilirubinemia, Elevation of total bilirubin* :occur due to

Over production of bilirubin-

Excessive hemolysis of RBC-

Liver diseases e.g. hepatitis &cirrhosis-

Obstruction of biliary tract e.g. gall stones-

Both conjugated (direct) and unconjugated-.(indirect) bilirubin are increased in hepatitis

The relative proportion of the conjugated fractionincreases with progression of the disease until eventually the liver loses its ability to carry out the conjugation reaction

Principle

Direct Bilirubin:

Bilirubin in the serum is coupled with diazotized sulfanilic acid to form azobilirubin (purple color).

The intensity of absorption of purple color that is formed is proportional to the bilirubin concentration in the serum.

This is the calculation of Direct Bilirubin

Total Bilirubin:

Add methanol to accelerates the reaction of unconjugated bilirubin in the serum, and a value of total bilirubin is obtained after letting the sample stand for 5 minutes.

***This is the calculation of Total Bilirubin**

The total bilirubin value represents the sum of the bilirubin glucoronide (direct) and the unconjugated (indirect) bilirubin.

Specimen

serum

While processing samples, protect from direct light as loss of bilirubin may occur.

Procedure

:Total Bilirubin-1

Tubes	Reagent blank	Test blank	Test	CAL		
DW	μl 100					
Test		μl 100	μl 100			
CAL				μl 100		
RT		ml 1				
WR	1ml		1ml	1ml		
.Mix and let stand for 2 min. at RT Wavelength= 540 nm Read A. of test blank against DW Read A. of samples against reagent blank						

Direct Bilirubin -2

Tubes	Reagent blank	Test blank	Test	CAL
DW	μ 100			
Test		μ 100	μ 100	
CAL				μ 100
RD		ml 1		
WR	ml 1		1ml	ml 1

.Mix and stand for 5 min. at 37 C Wavelength= 540 nm read A of test blank against DW-1 read A of samples against reagent blank-2

Calculation

:Same calculation for total and direct

Abs. test – Abs. test blankX Conc. Of CAL = mg/dlAbs. CAL

Normal Value

:Adult Total up to 1.0 mg/dl Direct up to 0.2 mg/dl