


BCH 445- Biochemistry of Nutrition [Practical]  
**Lab (9) Determination of Calcium in  
Milk**



# Calcium

---

- Calcium is an important component of a healthy diet and a **mineral** necessary for life.
- It is a mineral that people need to build and maintain strong bones and teeth.
- It is also very important for other physical functions, such as **muscle control and blood circulation**.
- If we do not have enough calcium in our diets to keep our bodies functioning, calcium is removed from where it is stored in our bones.
- Overtime, this causes our bones to grow weaker and may lead to **osteoporosis** (a disorder in which bones become very fragile).



# Milk and Calcium

---

- Milk is a heterogeneous mixture of proteins, sugar, fat, vitamins and minerals.
- Milk and milk products are some of the **natural sources of calcium**.
- The calcium concentration in bovine milk is about **1g/L**.
- Cow's milk has **good bioavailability** of calcium (about 30 to 35%).
- Milk is an excellent source of dietary calcium for those whose bodies tolerate it because it has a high concentration of calcium and the calcium in milk is **excellently absorbed**.
- It is estimated that without milk and milk products in the diet, less than half of the calcium requirements would be met.



# Practical Part


---

Objective:

- Determination of Calcium in milk sample.



# Principle

- In this experiment, The determination of calcium in milk is based on a **complexometric titration** of calcium with an aqueous solution of the disodium salt of EDTA at high pH value=12.
- **Complexometric titration** is a type of titration based on complex formation between the analyte and titrant.
- Such compounds (EDTA) are capable of forming **chelate complex** with many cations (**metal ion**) in which the cation is bound in a ring structure.
- The ring results from the formation of a salt-like bond between the cation and the carboxyl groups together with a coordinate bond through the lone pair of electrons of the nitrogen atom.



# Principle


- The common form of the agent is disodium salt  $\text{Na}_2\text{H}_2\text{EDTA}$ .
- It is colorless and can be weighed and dissolve in water to form a stable solution.
- At high pH ( $> 10$ ) the remaining protons leave EDTA forming  $\text{EDTA}^{4-}$  anion.




## Solochrome indicator

---

- The Solochrome dark blue indicator is a suitable indicator in this case.
- The dye **itself has a blue color**.
- This blue dye also forms a complex with the calcium ions changing color from **blue** to **pink/red** in the process, but the dye–metal ion complex is less stable than the EDTA–metal ion complex.
- As a result, when the calcium ion–dye complex is titrated with EDTA the  $\text{Ca}^{2+}$  ions react to form a stronger complex with the EDTA changing the dye color to **blue**.



At the beginning of reaction



Excess  $\text{Ca}^{2+}$   
ions present  
to complex  
with indicator

At the end of reaction



$\text{Ca}^{2+}$  ions  
almost all  
complexed by  
EDTA




All  $\text{Ca}^{2+}$  ions  
complexed by EDTA,  
indicator completely  
un-complexed



# How to determine calcium in the presence of Mg?

- This method for determining  $\text{Ca}^{2+}$  concentration in the presence of  $\text{Mg}^{2+}$  relies on the fact that the pH of the solution is sufficiently high ((The pH will be approximately 12.5 due to the addition of concentrated NaOH solution)) to ensure that all magnesium ions precipitate as magnesium hydroxide before the indicator is added.
- In this condition, magnesium ions are precipitated as hydroxide and do not interfere with the determination of calcium.



## Method

---

1. Combine 10mL of **sample**, 40mL distilled water, and 4mL of 8M sodium hydroxide solution into an Erlenmeyer flask and allow solution to stand for about 5 minutes with occasional swirling.
2. A small of magnesium hydroxide may precipitate during this time. Do not add the indicator until you have given this precipitate a chance to form.
3. Then add 6 drops of the **Solochrome** dark blue solution.
4. After that start to titrate with EDTA solution.
5. Repeat titration for three trials.

# Results

---

	<b>EDTA volume (ml)</b>
1	
2	
3	
Average	

$$\text{molarity (M)} = \frac{\text{number of moles}}{\text{liters of solution}}$$

## Calculation

1. Calculate the moles of EDTA required to complex the  $\text{Ca}^{2+}$  ions in the sample:

→ Number of moles (for EDTA) = Molarity of EDTA x volume of EDTA in L

Note : Ratio  $\text{Ca}^{2+}:\text{EDTA} = 1 : 1$  (i.e moles of EDTA = moles of  $\text{Ca}^{2+}$ )

Molarity of EDTA = 0.03408 M

1. Calculate weight of  $\text{Ca}^{2+}$  :

→ Weight of  $\text{Ca}^{2+}$  = Number of moles x molecular weight (40.78)

• % of  $\text{Ca}^{2+}$  = (weight of  $\text{Ca}^{2+}$  / weight of sample) x 100

OR

-  $\text{Ca}^{2+} \% = \frac{[\text{Molarity of EDTA} \times \text{vol. of EDTA (in liter)}] \times 40.78}{\text{Weight of the sample}} \times 100$

\* 40.78 g/mol is the molecular weight of  $\text{Ca}^{2+}$