




BCH 447

Triglyceride Determination in Serum

Introduction:

- Triglycerides are esters of fatty acids and are hydrolyzed by lipase to glycerol and free fatty acids.







- Triglyceride determinations when performed in conjunction with other lipid assays are useful in the diagnosis of **primary and secondary hyperlipidemia** (abnormally elevated of fat in blood)



Principle:

- Standard methods for the measurement of TG concentrations involved either enzymatic or alkaline hydrolysis to liberate glycerol.
- The glycerol concentration is then determined by enzymatic assay coupled with Trinder reaction that terminates in the formation of a quinoneimine dye.
- **The amount of the dye formed is** determined by its absorption at 505 nm, it is directly proportional to the concentration of triglycerides in the samples.



Trinder reaction: It is the reaction between hydrogen peroxide and the phenol and aminoantipyrine (AAP) to form quinoneimine (red-violet dye), catalyzed by the presence of a peroxidase

Specimen collection and storage:

1. Fresh, non-hemolyzed serum from fasting patients is recommended.
2. Triglycerides in serum appears stable for three days when stored at 2-8 °C.
3. Prolonged storage of the samples at room temperature is not recommended since **other glycerol containing compounds may hydrolyze**, releasing free glycerol with an apparent increase in total triglycerides content.

Method :

- By Triglyceride reagent kit.


-Follow the table:

| | Blank | Standard | Test |
|--|-------|-----------------|-----------------|
| Reconstituted Reagent | 1 ml | 1 ml | 1 ml |
| Pre-warm at 37°C for 2 min and add: | | | |
| Standard (200 mg/dl) | --- | 0.01 ml (10 µl) | --- |
| Sample | --- | --- | 0.01 ml (10 µl) |
| Mix and incubate at 37°C for 10 min | | | |
| ↓ | | | |
| Read the absorbance of standard and sample at 505 nm against blank | | | |

Calculation:

$$\text{Conc. of TG} = \frac{\text{Ab Test}}{\text{Ab Std.}} \times \text{conc. of Std. (200 mg/dl)}$$

Normal range of TG: 10 -190 mg/dl




TG → **Hyperthyroidism**

Malnutrition

Low-fat diet

Malabsorption



TG → **Hyperlipidemia**

Hypothyroidism

Obesity

Kidney Disease

Diabetes




HDL-Cholesterol determination




Introduction:

- Cholesterol is a fatty substance found in blood, bile and brain tissue.
- It serves as a precursor to **bile acids, steroids and vitamin D**.




•In the plasma, cholesterol is transported by three lipoproteins:

- high density lipoprotein (HDL-Cholesterol),
 - low density lipoprotein (LDL-Cholesterol),
 - and very low density lipoprotein (VLDL- Cholesterol).
-
- The concentration of **total cholesterol** in serum has been associated with metabolic, infectious and coronary heart diseases.

- The concentration of HDL-cholesterol in serum has important in diagnosis of the how the level of **risk to get coronary heart diseases**.
 - There is an **inverse relationship** between serum HDL-Cholesterol and the risk of coronary heart disease
 - **More HDL-Chol** → indicate **low risk** of coronary heart disease.
 - The measurement of **HDL Cholesterol and triglyceride** provides valuable information for the prediction of coronary heart disease and for lipoprotein phenotyping.
- 

Principle:

- Enzymatic methods, involving cholesterol esterase and oxidase and Trinders color system.
- The enzymatic reaction sequence employed in the assay of cholesterol is as follows:



- The amount of the dye formed is** determined by its absorption at $505 \pm 5 \text{ nm}$, it is directly proportional to the concentration of cholesterol in the samples.

Method :

- By HDL-Cholesterol reagent kit.

-Follow the table:

| | Blank | Standard | Test |
|---|--------|----------|--------|
| Cholesterol Enzymatic Reagent | 1ml | 1ml | 1ml |
| Pre-warm at 37°C for 2 min . | | | |
| Dis. water | 100 µl | _____ | _____ |
| Standard (50mg/dL) | _____ | 100 µl | _____ |
| Sample | _____ | _____ | 100 µl |
| Mix and incubate at 37°C for 10 min | | | |
| ↓ | | | |
| Read the absorbance of standard and sample at 505 nm against blank | | | |

Calculation :

$$\text{Conc. of HDL} = \frac{\text{Abs Test}}{\text{Abs Standard}} \times \text{conc. of Standard (50 mg/dl)} = \text{mg/dl}$$

- Normal value of HDL-Cholesterol :

- Low level (risk factor) < 40 mg/dl
- High HDL (protector factor) \geq 60 mg/dl

Questions :

- In a lipid profile test, what is being tested?**
- Is any test preparation needed to ensure the quality of the sample?**