



# The effect of substrate concentration on the rate of an enzyme catalyzed reaction





The rate of an enzyme catalyzed reaction depends directly on the concentration of an enzyme. With a fixed concentration of an enzyme and **with increasing substrate concentration , a rapid increase in the rate of the reaction is observed at first .**

But as substrate concentration continues to increase , the increase in rate of reaction *being slow down until , with a large substrate concentration, no further increase in the rate of the reaction obtained.*



*As the substrate level is increased* , the velocity increases in a hyperbolic fashion.

*At relatively low concentration* of substrate, the rate of reaction increase linearly with an increase in substrate concentration.

*At higher substrate concentration* the rate of reaction increase *smaller and smaller* amount in response to increase in substrate concentration.

*Finally* , a point is reached beyond which there is only *small increase* in the rate of the reaction with increasing substrate concentration. **This plateau is called maximum velocity,  $V_{max}$ .**

The hyperbolic shape of this curve can be expressed algebraically by the Michaelis -Menten equation:

$$V_o = V_{max} [S] / K_m + [S]$$

**$V_i$**  = initial velocity

**$V_{max}$**  = maximum velocity

**$[S]$**  = substrate concentration

**$K_m$**  = Michaelis -Menten constant

**$K_m$** : is a substrate concentration at half  $V_{max}$  ( $K_m$  indicate the affinity of an enzyme for its substrate)


Note that the larger the  $K_M$  (the weaker the binding), the larger the  $[S]$  needed to reach the half maximum rate.

The  $K_m$  can vary greatly from enzyme to enzyme , and even for different substrates of the same enzyme. **The Michaelis -Menten equation can be algebraically transformed** into forms that are useful in the practical determination of  $K_m$  an  $V_{max}$ . One common transformation is derived simply **by taking the reciprocal of both sides of the Michaelis -Menten equation to give Lineweaver – Burk equation:**

$$\frac{1}{V_i} = \frac{1}{V_{max}} + \frac{K_m}{V_{max}} \cdot \frac{1}{[S]}$$

By plotting  $1/v$  against  $1/[S]$   
a straight line plot,  
Lineweaver – Burk plot is  
obtained.

**Both  $V_{max}$  and  $K_m$  can be  
obtained accurately from  
intercepts of the straight  
line with the  $y$  – axis and  $x$ -  
axis**




# Objectives:

---

- 1) To establish the relationship between substrate concentration and the rate of an enzyme catalyzed reaction.
- 2) To determine the  $K_m$  and  $V_{max}$  of the enzyme for a particular substrate.





The object of this experiment is to demonstrate the effect of performing the standard 5 minutes assay in the presence of different substrate concentration ranging from 0 to 5.  $\mu\text{M}$  p-nitrophenyle phosphate. The results should provide classic Michaelis -Menten data from which approximations of  $V_{\text{max}}$  and  $K_m$  can be made.

**Double-reciprocal plots of the same data should be done to arrive at even more exact values for  $K_m$  and  $V_{\text{max}}$ .**

# Materials:



Test tubes

Pipettes

Cuvettes

Water bath

Stopwatch

Spectrophotometer

Chemicals:

Acid phosphatase enzyme

**0.05 M PNPP Different substrate concentration**

1 M sodium acetate buffer

0.1 M  $MgCl_2$

0.5 M KOH

# Method:

1. Prepare a series of substrate dilutions according to the protocol outlined in Tables I .

Reagent	A	B	C	D	E	F	G	H
0.05M PNP P ml	0	0.0 5	0.1	0.2 5	0.5	1	1.5	2
H <sub>2</sub> O(ml)	5	4.9 5	4.9	4.7 5	4.5	4	3.5	3
[S] M	...	...	...	...	...	...	.....	.....
	...	...	...	...	....	....	.	.

Concentration of stock Substrate 0.05M

1. Set up eight assay tubes. To each of these tubes **add 0.5 ml of 1.0 M sodium acetate buffer (pH 5.7), 0.5 ml of 0.1 M MgCl<sub>2</sub>, and**


- 3. add 0.5ml of the enzyme, run each reaction for 5 minutes**
- 4. Place the tubes in a test tube rack situated in a 37°C water bath and let stand for 5 minutes. stop it by adding 0.5ml of 0.5M KOH.**
- 5. Let the tube stand for 20 min. before determine the absorbance at 405 nm for each reaction mixture against the blank. (The tube containing no substrate should be used as the blank).**

# Result:

Tube	Absorbance at 405 nm	Substrate Conc.(x-values)	Velocity of reaction* (y-values)
B		.....	.....
C		.....	.....
D		.....	.....
E		.....	.....
F		.....	.....
G		.....	.....
H		.....	.....

- 1- Calculate the velocity ( M of p-nitrophenol/minute)**
- 2- Plot velocity against substrate concentration in the standard manner of Michaelis and Menten.**

# I) Michaelis –Menten Curve



**3-Determine  $V_{max}$  and  $K_m$  for acid phosphatase.**

**\*Velocity/Rate of reaction =  $A \times 10^3 / (\epsilon \times T)$**

Where:

A = Absorbance

$\epsilon$  = extinction coefficient for ACP enzyme =  $18.8 \times 10^3$

T = Time = 5 min.





**4- Calculate the reciprocals of velocity ( $1/v$ ) and substrate concentration ( $1/[S]$ ) and present these data as a table.**

Tube	Substrate Conc.	<u><math>1/(\text{Substrate Conc.})</math></u> (x-values)	Velocity of reaction	<u><math>1/(\text{Velocity of reaction})</math></u> (y-values)
B		.....		.....
C		.....		.....
D		.....		.....
E		.....		.....
F		.....		.....
G		.....		.....
H		.....		.....

## II) Lineweaver and Burk Curve

**5-Prepare the double –reciprocal plot of Lineweaver and Burk and determine the  $K_m$  and  $V_{max}$  from the x and y intercepts.**







**Thank You**