

# Buffer Capacity

# Buffers:

- Buffer solutions are solutions that can resist changes in pH upon addition of small amounts of acid/base.
- Common buffer mixtures contain two substances: conjugate acid and a conjugate base .
- Together the two species (conjugate acid and conjugate base) resist large changes in pH by absorbing the H<sup>+</sup> ions or OH<sup>-</sup> ions added to the system.

Buffer capacity = resistance

# How buffers resist the change in pH:

1. When **H<sup>+</sup> ions** are added to the buffer system they will react with the **conjugate base** in the buffer as following:



2. When **OH<sup>-</sup> ions** are added they will react with the **conjugate acid** in the buffer as following:



**→ Thus, the buffer is effective as long as it does not run out of one of its components.**

(There are enough conjugated base and conjugated acid to absorb the H<sup>+</sup> ions or OH<sup>-</sup> ions added to the system respectively).

# Buffer Capacity (Theoretically):

- Quantitative measure of buffer resistance to pH changes is called **buffer capacity**.
- Buffer capacity can be defined in many ways, **it can be defined as:**
  - The number of moles of H<sup>+</sup>/OH<sup>-</sup> ions that must be added to one liter of the buffer in order to decrease /increase the pH by one unit respectively.

- The instantaneous buffer capacity is **expressed as β** and can be derived from **Henderson Hasselbalch equation**:

$$\beta = \frac{2.3 K_a [H^+][C]}{(K_a + [H^+])^2}$$

- What the **relationship** between buffer capacity (β) and buffer concentration [C]?
- 0.1 M vs 0.2 M acetate buffer, which buffer has the **highest resistance**? Why?

- Where:** β = the buffer capacity, [H<sup>+</sup>] = the hydrogen ion concentration of the buffer, [C] = concentration of the buffer and K<sub>a</sub> = acid dissociation constant.
- Note:** The buffer capacity is **directly proportional** to the buffer concentration

# Buffer Capacity (Practical):

□ **Buffer capacity of acid and alkaline direction:**

→ Buffer capacity  $_a$  ( $BC_a$ ) = is the **concentration of  $H^+$**  that must be added to decrease the pH by one unit.

This called buffer capacity in the **ACID** direction.

$$BC_a = \frac{9[HA] [A^-]}{10 [HA] + [A^-]}$$

→ Buffer capacity  $_b$  ( $BC_b$ ) = is the **concentration of  $OH^-$**  that must be added to increase the pH by one unit.

This called buffer capacity in the **ALKALINE** direction.

$$BC_b = \frac{9[HA] [A^-]}{10 [A^-] + [HA]}$$


# Buffer capacity in acid and base direction:




$\text{H}_3\text{O}^+$  added, equilibrium position shifts to the left




$\text{OH}^-$  added, equilibrium position shifts to the right




Buffer solution  
after addition  
of strong acid



Buffer solution  
equimolar in  
acid and base  
**pH=pKa**



Buffer solution  
after addition  
of strong base




← Add  $\text{H}_3\text{O}^+$   
**Acid direction**

→ Add  $\text{OH}^-$   
**Base direction**

# Buffer capacity in acid and base direction:

## Base direction


Titration of Ethanoic Acid Buffer with Sodium Hydroxide Solution



(a) Volume of NaOH(aq)

## Acid direction

Titration of Ethanoic Acid Buffer with Hydrochloric Acid



(b) Volume of HCl(aq)


Figure 6 (a) Ethanoic acid buffer with a strong base added (b) Ethanoic acid buffer with a strong acid added. The pH changes quickly once all of the available buffer is depleted.

# Buffer capacity in acid and base direction:

The figure represents the buffer capacity when the buffer is titrated in **both directions**.

## Remember:

- Buffer capacity  $a$  ( $BC_a$ ) = is the **concentration of  $H^+$**  that must be added to decrease the pH by one unit.
- Buffer capacity  $b$  ( $BC_b$ ) = is the **concentration of  $OH^-$**  that must be added to increase the pH by one unit.



The buffer capacity curve of 0.05M Tricine buffer, pH 7.5 ( $pK_a=8.15$ )

💡 Can we calculate **practical buffer capacity** from the graph?

**Example:** Calculate the practical buffer capacity in the acid directions of a 0.1M and 0.2M acetate buffer, pH 5, pKa = 4.76.

First, calculate the concentration of the weak base and its conjugated acid that make up the buffer with 0.1M:

$$pH = pKa + \log \frac{[A^-]}{[HA]} \rightarrow 5 = 4.76 + \log \frac{[y]}{[0.1-y]} \rightarrow 0.24 = \log \frac{[y]}{[0.1-y]} \rightarrow \text{Anti log for both sides}$$

$$1.74 = \frac{y}{0.1-y} \rightarrow y=0.063\text{M. SO: } [A^-]=0.063\text{M} , [HA]=0.037\text{M}$$

Second: Calculate the practical buffer capacity in both directions

$$BC_a = \frac{9 [HA][A^-]}{10 [HA]+[A^-]} \rightarrow \frac{9 \times 0.037 \times 0.063}{(10 \times 0.037) + 0.063} \rightarrow \frac{0.021}{0.433} = 0.048\text{M } [H^+]$$

□ **Note:** do the same calculation for the same buffer when its concentration equals 0.2M.

# Practical Part



# Objective:

---

- To understand the concept of buffer capacity.
- To determine the capacity of two acetate buffers in the acid directions.

# Method:

- You are provided **0.1M** and **0.2M** acetate buffer (pH=5).
- In one beaker add **10 ml** of the 0.1M acetate buffer, and in the other add 0.2M acetate buffer
- Titrate the two beakers by adding **1 ml of 0.1 M HCl** from the burette and determine the pH of the solution after each addition.
- Continue adding the acid/base until the pH drops by two units from your initial reading
- Record the values in the titration table.

## Practical notes:

- Wash all glassware with distilled water followed by washing with the solution used.
- Check the flow of your burette and ensure accurate meniscus.
- No need to wash electrode after each addition, since the same solutions are used.

# Results:

- Plot the capacity curve (pH against the volume (ml) of 0.1M HCl).
- For both buffers, determine the practical buffer capacity in the acid direction from the graph and the formula then summarize your value in the table:

Acetate buffer	Practical capacity (from the formula)	Practical capacity (from the curve)
0.1M	0.048M [H <sup>+</sup> ]	
0.2M		

# Results:

- Buffer capacity  $a$  ( $BC_a$ ) = is the **concentration of  $H^+$**  that must be added to decrease the pH by one unit.

- To determine the capacity from the graph:
  - a) Find the ml of 0.1M HCl needed to drop the pH one unit from the initial reading value.
  - b) Then find the final concentration of the HCl.

## Example from the curve:

3.8 ml of 0.1M HCl is needed to drop the pH from 3.8 to 2.8 of 10 ml of acetate buffer. Thus:

$$C_1 \times V_1 = C_2 \times V_2$$

$$0.1M \times 3.8 \text{ ml} = ? M \times 13.8 \text{ ml}$$

$$= 0.027M [H^+]$$




Figure: 0.1M acetate buffer capacity in the acid directions