

BCH 202 General Biochemistry [Practical]

Lab (1) Preparation of Buffer Solutions



Emtenan Mohammed Alkhudair

Office: Building 5, 3rd floor, Office No. 269

Email: Ealkhudair@ksu.edu.sa

Website: <http://fac.ksu.edu.sa/ealkhudair>

Office hours: Sun, Tue 9-11

Hydrogen number (pH):

- The **acidity** of certain solutions can be described by using hydrogen number (pH).
- **pH defined as:** The negative logarithm of the hydrogen ion concentration.



$$\text{pH} = -\text{Log}[\text{H}^+]$$

- When the pH increase the concentration of hydrogen ion decrease and vice versa.
- The pH range goes from 0 to 14.


Acidic → pH < 7

Neutral → pH = 7



Basic → pH > 7






Acidic solution : has higher concentration of hydrogen ions $[H^+]$ than hydroxyl ions $[OH^-]$



Neutral solution: has equal concentration of hydrogen ions $[H^+]$ and hydroxyl ions $[OH^-]$



Basic solution : has lower concentration of hydrogen ions $[H^+]$ than hydroxyl ions $[OH^-]$



Measuring of hydrogen number:


➤ To measure the hydrogen number in certain solution in very accurate way, we use a special instrument called **pH meter**. It's consist of two electrodes:

1. **Reference electrode:** contains silver-silver chloride wire immersed in saturated KCl solution.
2. **Glass electrode:** which contains a very thin bulb, that is sensitive to pH.

This device measures the difference between the electrodes, and converts it into a **pH from 0 to 14**.

➤ **Supporting materials:**


1. How pH meter works: shorturl.at/fkpL3
2. How to use the pH meter: shorturl.at/IKXZ4



Measuring of hydrogen number cont.:


Test strip (inaccurate):

A pH test strip is a strip of litmus paper with which you can measure the pH value of a liquid which show a different color at different acidities.




pH and biological system:

- All **biochemical reactions** occur under strict conditions of the concentration of hydrogen ion.
- Biological life can not withstand large changes in hydrogen ion concentrations which we measure as the pH.
- So how to resist changes in pH ? **Buffers.**



It is vitally important to maintain the body acid alkaline balance at the correct pH level to enjoy good health and avoid degenerative disease.



Buffers:

- So, **buffers defines as:** the solutions that have the ability **to resist changes in pH** upon the addition of **limited** amounts of acid or base.
- A buffer is made up of:

a weak acid and its conjugate base.



Or

a weak base and its conjugate acid .

الحمض الضعيف + قاعدته المقترنة



القاعدة الضعيفة + حمضها المقترن



Two types of buffers

A buffer is made up of a **weak acid** and its conjugate base.

Or

A **weak base** and its conjugate acid.

Acidic Buffer

Are made from weak acid and its conjugated base [its salt].

Example:

1. $\text{CH}_3\text{COOH} / \text{CH}_3\text{COONa}$ (Pka)

→ CH_3COOH (Weak acid)

→ CH_3COONa (conjugated base –its salt-)

2. $\text{NaH}_2\text{PO}_4 / \text{Na}_2\text{HPO}_4$ (Pka)

Basic Buffer

Are made from weak base and its conjugated acid [its salt].

Example:

1. $\text{NH}_3 / \text{NH}_4\text{Cl}$ (Pkb)

→ NH_3 (Weak base)

→ NH_4Cl (conjugated acid –its salt-)

Mechanism of action:


➤ How buffers can resist the change in pH?

-Example using [HA/A⁻] buffer:

➔ Where: HA is Weak acid and A⁻ is conjugated base [its salt].



If H⁺ (acid) is added to this buffer system → H⁺ will react with conjugated base → to give conjugate acid.



If OH⁻ (base) is added to this buffer system → OH⁻ will react with conjugated acid → to give conjugate base and H₂O.




Mechanism of action cont.:

➤ Example:

Buffer system: $\text{CH}_3\text{COOH} / \text{CH}_3\text{COO}^-$, (CH_3COOH :acid - CH_3COO^- : conjugated base)

-When acid [H^+] added: $\text{CH}_3\text{COO}^- + \text{H}^+ \rightarrow \text{CH}_3\text{COOH} \uparrow$

-When base [OH^-] added: $\text{CH}_3\text{COOH} + \text{OH}^- \rightarrow \text{CH}_3\text{COO}^- \uparrow + \text{H}_2\text{O}$




- **NOTE:** It resists pH changes when its two components are present in specific proportions.
- Thus, a buffer can protect against pH changes from added H^+ or OH^- ion as long as there is **sufficient** basic and acidic forms respectively ➔ **As soon as you run out of one of the forms you no longer have a buffer.**

Henderson-Hasselbalch equation:

- **The Henderson-Hasselbalch equation** is an equation that is often used to:
 1. To prepare Buffer.
 2. To calculate the pH of the Buffer.

$$\text{pH} = \text{pka} + \log \frac{[\text{A}^-]}{[\text{HA}]}$$

- It relates the **Ka** [dissociation constant] of a weak acid, **[HA]** concentration of weak acid, **[A-]** concentration of conjugate base [salt of the weak acid] components and the **pH of the buffer**.



Choosing the proper buffer:

- A buffer is best used close to its pKa [to act as a good buffer the pH of the solution must be within one pH unit of the pKa].



➔ The **buffer capacity** is optimal when the ratio of the weak acid to its salt is 1:1; that is, when $\text{pH} = \text{pKa}$.

$$\text{pH} = \text{pka} + \log \frac{[\text{A}^-]}{[\text{HA}]}$$

$$\text{pH} = \text{pka} + \log 1$$

$$\text{pH} = \text{pka} + 0$$

$$\boxed{\text{pH} = \text{pka}}$$



Buffer capacity:

- Quantitative measure of buffer resistance to pH changes is called **buffer capacity**.
- **Buffer capacity** can be defined in many ways, it can be defined as:


The number of moles of H^+/OH^- ions that must be added to one litre of the buffer in order to decrease /increase the pH by one unit respectively.


- Buffer capacity is **directly proportional** to the buffer concentration.

A buffer has concentration of **0.5 M**

A buffer has concentration of **0.9 M**

Which buffer has the highest capacity?






Practical Part



Objectives:

- To learn how to prepare buffers.
- To understand the behaviour of buffers solutions.



A. Preparation of phosphate buffer:

Prepare 50 ml from phosphate buffer with concentration 0.25M and pH=7.4, if you know that (pKa=7.2).

You are provided with buffer solution content: Monosodium dihydrogen phosphate NaH_2PO_4 and Disodium hydrogen phosphate Na_2HPO_4

Solution:

■ Provided:


$$\text{pKa} = 7.2 \quad \text{pH} = 7.4$$

Final volume of buffer = 50 ml

Concentration of buffer = 0.25 M \rightarrow [HA] + [A⁻]

■ Required:

Weight (g) of NaH_2PO_4 (as HA) and Na_2HPO_4 (as A⁻).



Calculations:

-To prepare a buffer Henderson-Hasselbalch equation is used:

$$\text{pH} = \text{pKa} + \log \frac{[\text{A}^-]}{[\text{HA}]}$$

1. First calculate the concentration of the weak acid and its conjugated base that make up the buffer with 0.25 M:

$$\rightarrow \text{Assume } [\text{A}^-] = y \quad \text{and} \quad [\text{HA}] = 0.25 - y$$

So:

$$7.4 = 7.2 + \log \frac{y}{0.25-y} \rightarrow 0.2 = \log \frac{y}{0.25-y}$$


By taking the “Anti log for both sides” :

$$1.6 = \frac{y}{0.25-y} \rightarrow y = 1.6 \times (0.25 - y) \rightarrow y = 0.4 - 1.6y \rightarrow y + 1.6y = 0.4 \rightarrow 2.6y = 0.4$$

y = 0.15 M [which is the concentration of **[A⁻]** in the buffer]

So, **[HA] = 0.25 - 0.15 = 0.1 M** [which is the concentration of **[HA]** in the buffer]

$$\left. \begin{array}{l} \\ \\ \end{array} \right\} 0.15 + 0.1 = 0.25\text{M}$$



Calculations cont.:

Molecular weight
(A) $\text{Na}_2\text{HPO}_4 = 142 \text{ g/mol}$
(HA) $\text{NaH}_2\text{PO}_4 = 120 \text{ g/mol}$

2. Calculate the **weight** in (g) needed from [A-] to prepare the buffer, so number of mole of [A-] should be calculated first :

→ **Calculate moles of A- in buffer: [1]**

$$\begin{aligned}\text{No. of mole (of A-)} &= \text{molarity (of A- calculated in the buffer)} \times \text{volume L (volume of the buffer)} \\ &= 0.15 \times 0.05 = 0.0075 \text{ mole}\end{aligned}$$

→ **Calculate weight of A- needed: [2]**

$$\begin{aligned}\text{Weight in (g) of [A-]} &= \text{No. of moles} \times \text{MW} \\ &= 0.0075 \times 142 = \boxed{1.065 \text{ g}}\end{aligned}$$

3. Calculate the **weight** in (g) needed from [HA] to prepare the buffer, so number of mole of [HA] should be calculated first :

→ **Calculate moles of HA in buffer: [1]**

$$\begin{aligned}\text{No. of mole (of HA)} &= \text{molarity (of HA calculated in the buffer)} \times \text{volume L (volume of the buffer)} \\ &= 0.1 \times 0.05 = 0.005 \text{ mole}\end{aligned}$$

→ **Calculate weight of HA needed: [2]**

$$\begin{aligned}\text{Weight in (g) of [HA]} &= \text{No. of moles} \times \text{MW} \\ &= 0.005 \times 120 = \boxed{0.6 \text{ g}}\end{aligned}$$

$$\text{[1] Molarity} = \frac{\text{No. of moles of solute}}{\text{Volume (L)}}$$

$$\text{[2] No. of moles} = \frac{\text{Wt(g)}}{\text{Mwt}}$$

Homework:

➤ You are provided with 0.15 M acetic acid and sodium acetate.

Prepare 100 ml of a 0.2M acetate buffer, pH =5.2 if you know that pK_a =4.76.

Hint: you will calculate ml of acetic acid and g of sodium acetate.

