



Qualitative Tests of Carbohydrate-I

Introduction:

- Carbohydrates are defined as the **polyhydroxy aldehydes or polyhydroxy ketones**.
- Most, but not all carbohydrate have a formula $(CH_2O)_n$ (hence the name hydrate of carbon).
- In human body, **the D-glucose** is used.
- Simple sugars ends with **-ose**.



- **Biological role:**
 - Carbohydrates are the **key source of energy** used by living things.
 - Also serve as extracellular structural elements as in cell wall of bacteria and plant.

Classification of carbohydrates:

- Several classifications of carbohydrates have proven useful, and are outlined in the following table:

Complexity	Simple Carbohydrates monosaccharides	Complex Carbohydrates disaccharides, oligosaccharides & polysaccharides		
Size	Tetrose C ₄ sugars	Pentose C ₅ sugars	Hexose C ₆ sugars	Heptose C ₇ sugars etc.
C=O Function	Aldose sugars having an aldehyde function or an acetal equivalent. Ketose sugars having a ketone function or an acetal equivalent.			
Reactivity	Reducing sugars oxidized by Tollens' reagent (or Benedict's or Fehling's reagents). Non-reducing sugars not oxidized by Tollens' or other reagents.			

Classification based on complexity:

1. Simple sugar (one unit) :

Monosaccharides contain one monosaccharide unit.

2. Complex sugar (more than one) :


1. **Disaccharides** contain **two** monosaccharide units.

2. **Oligosaccharides** contain **3-9** monosaccharide units.

3. **Polysaccharides** can contain **more than 9** monosaccharide units.

➤ Complex carbohydrates can be broken down into smaller sugar units through a process known as **hydrolysis**.


Classification based on complexity:




Classification based on reactivity:

➤ Reducing and non Reducing sugar :


If the oxygen on the anomeric carbon of a sugar is **not attached** to any other structure, that sugar can act as a **reducing agent** and is termed a **reducing sugar**.



FREE anomeric carbon
Reducing sugar



attached anomeric carbon
Non-Reducing sugar



FREE anomeric carbon
Reducing sugar

Classification based on reactivity cont':

- **All monosaccharides** are **reducing sugars**; they all have a free reactive carbonyl group.
- **Some disaccharides** have exposed carbonyl groups and are also **reducing sugars** like **lactose**. While other disaccharides such as **sucrose** are **non-reducing sugars** and will not react with Benedict's solution.
- **Large polymers** of glucose, such as starch, are **not reducing sugars**, since the concentration of hemiacetal groups is very low.

Solubility of sugars [physical property]:

- **Monosaccharide** and **disaccharide** can be **dissolved freely** in water because water is a polar substance.
- **Polysaccharide** cannot be dissolved easily in water, because, it has **high molecular weight**, which give colloidal solutions in water.



Practical Part

Chemical properties of carbohydrates:

1

Molisch test: To identify the carbohydrate from other macromolecules.

2

Benedict test: For the presence of reducing sugars.

3

Barfoed's Test: for to distinguish between reducing monosaccharides, reducing disaccharides and non reducing di-polysaccharides (detect reducing monosaccharides).

4

Bial's Test: To distinguish between pentose monosaccharide and hexose monosaccharide (to detect pentoses).

5

Seliwanoff's Test: To distinguish between aldoses and ketoses (to detect ketoses).

Experiment I : Molisch test

➤ Objective:

To identify the carbohydrate from other macromolecules lipids and proteins (**this test is specific for all carbohydrates**).


➤ Principle:

Two solutions are added : H_2SO_4 and α -naphthol

1. The test reagent (H_2SO_4) dehydrates pentose to form **furfural** and dehydrates hexoses to form **5-hydroxymethyl furfural**.
2. The furfural and 5- hydroxymethyl furfural further react with α -naphthol present in the test reagent to produce a **purple product**.





(1)



Pentose sugar

(2)



Hexose sugar



Experiment I : Molisch test

➤ Method:


1. Two ml of a sample solution is placed in a test tube.
2. 0.5 ml of the Molisch reagent (which α -naphthol in 95% ethanol) is added.
3. The solution is then poured slowly into a tube containing two ml of concentrated sulfuric acid so that two layers form, producing **violet ring** appear as liaison between the surface separations.

➤ Results:

Tube	Observation
Glucose	
Lactose	
Starch	

⚠ CAUTION

Concentrated sulfuric acid is extremely corrosive and can cause serious burns when not handled properly.




Experiment 2 : Benedict's test

➤ Objective:

To distinguish between the reducing and non-reducing sugars (to detect the presence of reducing sugar).

➤ Principle:

- The copper sulfate (CuSO_4) present in Benedict's solution reacts with electrons from the **aldehyde or ketone group of the reducing sugar** in **alkaline medium**.
- Reducing sugars are oxidized by the copper ion in solution to form a carboxylic acid and a **reddish precipitate** of copper oxide.
- The non-reducing sugars give negative result.




Experiment 2 : Benedict's test

➤ Method:

1. One ml of a sample solution is placed in a test tube.
2. Two ml of Benedict's reagent is added.
3. The solution is then heated in a boiling water bath for five minutes.
4. A positive test is indicated by: The formation of a reddish precipitate.

➤ Results:

Tube	Observation
Glucose	
Lactose	
Sucrose	




Experiment 3 : Barfoed's test

➤ Objective:

This test is performed to distinguish between reducing monosaccharides, reducing disaccharides and non reducing di- and polysaccharides.

➤ Principle:

- Barfoed's test used **copper (II)** ions in a slightly **acidic medium**.
- Reducing saccharides are oxidized by the copper ion in solution to form a carboxylic acid and a **reddish precipitate** of copper (I) oxide.
- **Different types of reducing sugars react at different rates** ➔ Reducing **monosaccharides** react **quickly** with Barfoed's reagent (acidic condition), but reducing **disaccharides** react very **slowly** or not at all.
- The non-reducing sugars give negative result.



Experiment 3 : Barfoed's test

➤ **Method:**

1. Place one ml of a sample solution in a test tube.
2. Add 3 ml of Barfoed's reagent (a solution of cupric acetate and acetic acid).
3. Heat the solution in a boiling water bath for 6 minutes (after the 3 min check the tubes).

➤ **Results:**

Tube	Observation
Glucose	
Lactose	
Starch	



Glucose (+)





Experiment 4 : Bial's test

➤ Objective:

To distinguish between pentose monosaccharide and hexose monosaccharide (to detect pentoses).

➤ Principle:

- Bial's reagent (orcinol, concentrated HCl as a dehydrating acid and ferric chloride as condensation reagent).
- The test reagent dehydrates **pentoses** to form furfural, furfural further reacts with orcinol and the iron ion present in the test reagent to produce a **bluish or green product**.
- While **hexoses** yield **muddy-brown to grey** condensation product.



Experiment 4 : Bial's test

➤ Method:

1. Put 2 ml of a sample solution in a test tube.
2. Add 2 ml of Bial's reagent to each tube.
3. Heat the tubes gently in hot water bath.
4. If the color is not obvious, more water can be added to the tube.

➤ Results:

Tube	Observation
Glucose	
Ribose	



Ribose
(+)



Glucose
(-)



Experiment 5 : Seliwanoff's test

➤ Objective:

To distinguish between aldoses and ketoses (to detect ketoses).

➤ Principle:

- Seliwanoff's test uses 6M HCl as dehydrating agent and resorcinol as condensation reagent.
- The test reagent dehydrates **ketohexoses** to form 5-hydroxymethylfurfural → 5-hydroxymethylfurfural further condenses with resorcinol present in the test reagent to produce a **cherry red** product within **two minutes**.
- Aldohexoses** react to form the same product, but do so **more slowly** giving **yellow to faint pink color**.




Experiment 5 : Seliwanoff's test

➤ **Method:**

1. One half ml of a sample solution is placed in a test tube.
2. Two ml of Seliwanoff's reagent (a solution of resorcinol and HCl) is added.
3. The solution is then heated in a boiling water bath for two minutes.

➤ **Results:**

Tube	Observation
Glucose	
Fructose	



Homework:

- Explain, although starch has free hemiacetal bond it gives negative Benedict test?**
- What is the difference between Benedict and Barfoed's reaction?**