

Quantitative amino acids estimation by Ninhydrin method

BCH303 [Practical]

Amino acid quantitation :

- Free amino acids, as well as amino acids released from macromolecules such as peptides, proteins or glycoproteins.
- Rapid and accurate.
- **Importance ?**

Amino acid analysis :

- Refers to the methodology used to determine the amino acid composition or content of proteins, peptides, and other pharmaceutical preparations.
- **Methods ?**

Ninhydrin:

- Detecting amino acids.
- Microgram amounts.
- Non-selectivity ?
- Other convenient reagents are available which can react with the alpha amino group to form colored or fluorescent derivatives. These include fluorescamine, dansyl chloride, dabsyl chloride, etc., **used in the detection of trace amounts of amino acids at the nanogram level.**

Quantitative estimation of amino acid using Ninhydrin reagent:

- Ruhemann's purple (RP) was discovered by Siegfried Ruhemann in 1910.
- Ruhemann's purple formed by the reaction at 570 nm is measured, whereas for **imino acids**, the absorption happens at 440 nm.
- The color intensity produced is proportional to the amino acid concentration.
- **So ?**

Standard curve :

- **Direct relationship between color and concentration → direct relationship between concentration and absorbance.**
- **The standard curve (also called calibration curve):** is a type of graph used as a quantitative research technique that shows the relationship between different **known concentrations** of a substance and the **absorbance at a specific wave length**.
- Is most commonly used to determine the concentration of a substance (**unknown**), using serial dilution of solutions (**Standard solutions**) of known concentrations.
- **How ?**

Standard curve :




Figure 1. A standard curve showing the relation between the absorbance of different concentrations of a substance.

Practical part

Aim:

- Determination of amino acids quantity using ninhydrin reaction.
- Understanding and constructing a standard curve.

Principle:

- At neutral pH, ninhydrin destroys each primary α -amino acid to form Ruhemann's purple.
- Maximum **absorption at about 570 nm.**
- Imino acids yields a yellow- orange product at neutral pH.
- The intensity of the color resulted is linearly proportional to the concentration of the amino acids present in the solution

Results:

Table 1. Concentration of standard amino acid solution and their absorbance at 570 nm.

Test tube	Amino acid concentration [μg/ml]	Absorbance at 570 nm
Blank		
A		
B		
C		
D		
E		
Unknown sample	_____	




Figure 1. standard curve of amino acid using ninhydrin method.