

Experiment (2) Separation of Main Proteins in Plasma and Serum

Blood contents





immunoglobulins

Albumin

It is the most abundant protein present in plasma.

It has many functions including:

☑ Maintenance of the blood osmotic pressure.

✓ Adjusting blood pH.

Act as a transporter, transporting free fatty acids, bilirubin, drugs, steroid hormones, calcium and copper in the blood.

Globulins

FIBRINOGEN

 Alpha (α) and Beta (β) globulins are transport proteins, but γ-globulins are part of the immune system.

- It is a glycoprotein (proteins that contain oligosaccharide chains)
- It is converted by thrombin into fibrin during blood coagulation



Serum Proteins

Total serum protein consists of two main fractions, albumin and globulin.

In normal people the A / G ratio is from 1.2 to 1.5.

Generally, the decrease in total protein is due to a decrease in albumin fraction and the increase is due to an increase in globulin components.

Dehydration is one condition in which the increase in total protein is due to increase in both albumin and globulin fractions because of haemoconcentration \rightarrow In this case the A / G ratio remains unaltered.



- A low serum albumin may be due to:
- A heavy loss of albumin in urine.
- Loss or malabsorption of protein from the digestive tract.
- Decreased formation by the liver due to defective liver.
- Increase catabolism of protein or due to insufficient intake of protein in diet.
- Globulin:
- Increase in globulin occurs most commonly in advanced liver disease, multiple myeloma and a number of chronic infections.

• Plasma protein can be separated from each other by :

Salting Out

*⊡***Ultracenterifugation**



☑ Chromatography

Principle of salting out



- When high concentrations of salt is added to the protein solution, the solubility decreases, and the protein precipitates.
 - This can be explained by the following:
 - The salt molecules compete with the protein molecules in binding with water, leading to dehydration.

The salt concentration needed for the protein to precipitate out of the solution differs from protein to protein.

Principles of identification tests

Biuret test

In the presence of peptides that contain at least two peptide bonds(i.e. it is not given by **dipeptides and free amino acids**), a copper(II) ion forms violet/blue-colored complexes in an alkaline solution.

Protein + Biuret reagent \rightarrow Blue Color

The intensity of the color is proportional to the number of peptide bonds and thus is a measure of the concentrations of proteins.

Heat coagulation

Protein + weak acid heating protein precipitate (cloudiness)