




Characterization of DNA by Spectrophotometric Assay and Melting Temperature (T_m)

After DNA isolation:



Characterization of extracted DNA by spectrophotometric assay:

1. DNA concentration.
2. DNA purity.


By measuring the absorption of ultraviolet light (UV).

- DNA has maximum absorption at **260nm**. **WHY?**



1st. UV for quantification of nucleic acid concentration :

- Is determined by measuring absorbance at **260 nm**. **WHY?**
- At 260 nm double-stranded DNA has specific absorption coefficient of $0.02 (\mu\text{g/ml})^{-1}\text{cm}^{-1}$.
- So:
→ **Concentration of DNA** = $(A_{260} / \epsilon L) \times \text{Dilution Factor (DF)}$. 



Beer-Lambert Law:

$$A = \epsilon cl$$



2nd. DNA purity:

1. To detect nucleic acid purity from proteins contamination:

→ Calculate A_{260}/A_{280} WHY?


- Highly purified DNA samples have a A_{260}/A_{280} nm ratio of (1.8-1.9).
- What if the ratio is below 1.8? What that means?
- What if the ratio is higher than 1.9? What that means?

2. To detect nucleic acid purity from carbohydrates, peptides, ethanol or any organic compounds:

→ Calculate A_{260}/A_{230} WHY?

- Purified DNA samples have a A_{260}/A_{230} nm ratio of (2-2.2).

DNA and protein absorption spectrum:





***What is the effect of the contaminants on DNA concentration?**

*** What if the samples contaminated by proteins or organic compound?**



3^{ed}. GC content:

- The two strands of a DNA molecule can be dissociated ("melted") into single strands by **heat or altered pH**, which breaks the hydrogen bonds between complementary bases (A=T and G≡C).
→ What that process called ?
- Hyperchromic and hypochromic effect.
- Melting temperature profile.
- The **melting temperature (T_m)** is the temperature at which **50% of the DNA is unpaired** (denatured).
→ Is the T_m same for all DNA molecules ? WHY?
→ What is the important of knowing T_m of DNA?

DNA melting curve:




FIGURE 4.4 DNA melting curve. A melting curve of DNA showing T_m (the melting temperature) and possible molecular conformations for various degrees of melting.



3rd. GC content:

- GC content can be calculated by generating T_m profile (DNA melting curve).

$$\%(G+C) = 2.44 (T_m - 69.3)$$

Relationship between T_m and GC%:





FIGURE 4.5 Effect of G-C content on DNA melting temperature. T_m increases with increasing percent of G + C.



What do you notice about the GC content in relation to T_m ?



Practical Part



Aim:

- Determination the concentration and purity of extracted DNA using UV spectrophotometer. [SEP]
- Determination of DNA melting temperature and GC content percentage.

Principle:

- DNA and proteins have a maximum absorbance at 260 and 280 respectively.
- dsDNA will be separated to ssDNA by heat (denaturation).
- O.D at 260 nm will increase during denaturation... Why?
- Temperature for midpoint of denaturation gives T_m .
- The DNA of each species has a specific denaturation curve.. Why?

Results:

1. Characterization of DNA by Spectrophotometric Assay:

| Wavelength (nm) | Absorbance of DNA | |
|-----------------|-------------------|-------|
| | Blood | Plant |
| 230 | | |
| 260 | | |
| 280 | | |

→ Find out the concentration of DNA using the following equation:

$$\text{Concentration of DNA} = (A_{260} / \epsilon L) \times \text{Dilution Factor}$$

→ Determine the purity of the DNA.

Results:

2. Melting Temperature of DNA: [L] [SEP]

| Temperature (°C) | DNA Absorbance at 260 nm | |
|-------------------|--------------------------|-------|
| | Blood | Plant |
| 25 | | |
| 50 | | |
| 60 | | |
| 70 | | |
| Boiling | | |

→ Plot the value of absorbance vs. temperature and calculate the T_m for sample DNA.

→ Find out the GC content of your sample using the following formula:

$$\text{Percent of G + C} = (T_m - 69.3) \times 2.44 [L] [SEP]$$

[L] [SEP]

Home Work:

• Watch the following videos :

- <https://www.youtube.com/watch?v=wXiiTW3pflM>
- https://www.youtube.com/watch?v=U2-5ukpKg_Q

→ What is agarose gel electrophoresis?