# Effect of various factors on protein solubility and structure

BCH303 [Practical]

### Proteins :

• Proteins are polymers of <u>amino acids.</u>



- Peptide bond.
- How peptide bond formed?

→ By removal of the elements of water (dehydration) from the  $\alpha$ -carboxyl group of one amino acid and the  $\alpha$ -amino group of another.



2

## **Proteins precipitation :**

- What is ?
- Proteins precipitation is widely used in downstream processing of biological products in order to concentrate proteins and purify them from various contaminants.
- Factors?
- The change of one of these factors will lead to protein precipitation and/ or denaturation.



### **Proteins denaturation:**

• **Denaturation** is a process in which the proteins **losing its quaternary structure, tertiary structure and secondary structure**, by application of some external factor or compound such as a strong acid or base, an organic solvent (e.g., alcohol or chloroform), or heat.



- No alteration on the molecule's primary structure.
- Solubility ?
- <u>Activity?</u>

# **Practical part**

## Tests of proteins



Effect of strong acids on protein solubility and structure.

2

3

Effect of salts of heavy metals on protein solubility and structure.

Effect of heat on protein solubility and structure.

# Experiment 1 : Effect of salt concentration on the protein solubility

#### **Objective:**

• To investigate the effect of different salt concentration on protein solubility.

#### **Principle:**



#### • Notes:

1.Each protein can be precipitated at <u>specific</u> salt concentration.

2.It is <u>reverse process</u>, the protein can again become soluble when we add water .

# Experiment 1 : Effect of salt concentration on the protein solubility

### Method:

1.Label one tube as A.

- 2. Add 2ml of albumin.
- 3. Add drops of **0.1M NaCl** solution, Concentrate your vision on the tube while adding.
- 4. Record your results.
- 5. In the same tube add few amounts of 100% solid  $(NH_4)_2SO_4$ , shake it well.
- 6. Record your results.
- 7. Compare between the two results.

Tube	Observation
Albumin + NaCl	
Albumin+100%saturate (NH <sub>4</sub> ) <sub>2</sub> SO <sub>4</sub>	

# Experiment 2 : Effect of strong acids on protein solubility and structure

#### **Objective:**

• To investigate the effects of strong acids on the protein solubility.

#### **Principle:**

- This test depend on affecting solubility of the protein as a function of changes in pH.
- In **highly acidic media**, the protein will be <u>positively charged</u>, which is attracted to the <u>acid</u> <u>anions</u> that cause them to precipitate.
- Applications:
- Detection of small amount of protein in urea sample.
- ➢ Stop the enzyme reaction.

# Experiment 2 : Effect of strong acids on protein solubility and structure

#### Method:

1.Label two tubes A and B.

2. In tube A: add 3 ml of conc. nitric acid (HNO<sub>3</sub>) CAREFULLY.

3. Then, Using a dropper add drops of albumin on the inner wall of tube A to form a layer up the acid.

4. Record your results.

5. In tube B: Add 3 ml of the albumin solution.

6. Then add 5-7 drops of TCA solution CAREFULLY.

7. Record your results.

Tube	Observation
Albumin + $HNO_3$	
Albumin+TCA	





# Experiment 3 : Effect of salts of heavy metals on protein solubility and structure

#### **Objective:**

• To identify the effect of heavy metal salt on protein.

#### **Principle:**

- Heavy metal salts usually contain Hg<sup>+2</sup>, Pb<sup>+2</sup>, Ag<sup>+1</sup> Tl<sup>+1</sup>, Cd<sup>+2</sup> and other metals with high atomic weights.
- Heavy metal salt will **neutralize the protein**.
- The protein will precipitate as insoluble metal protein salt.



#### • Applications:

➤ To eliminate the poisoning by palladium Pb++ ,.....mercury salts Hg++

# Experiment 3 : Effect of salts of heavy metals on protein solubility and structure

### Method:

1.Label two tubes A and B.

- 2. In tube A and B add 1 ml of Albumin sample.
- 3. In tube A: using a dropper add few drops of AgNO<sub>3</sub>.
- 4. Record your results.
- 5. In tube B: using a dropper add few drops of HgCl<sub>2</sub>.
- 6. Record your results.

Tube	Observation
Albumin $+$ AgNO <sub>3</sub>	
Albumin + $HgCl_2$	



# Experiment 4 : Effect of heat on protein solubility and structure

#### **Objective:**

• To investigate the effect of high temperature on protein structure.

#### **Principle:**

• Non-covalent bond can be broken by heating, leading to protein denaturation and the precipitation.



# Experiment 4 : Effect of heat on protein solubility and structure

#### Method:

- 1- Take 3 ml of protein Albumin.
- 2- Place it in a boiling water bath for 5-10 minutes
- 3-Remove aside to cool to room temperature.
- 4-Note the change.

Tube	Observation
Albumin + heating	



### Homework:

• From today lab, which factors lead to protein denaturation and which lead to precipitation? Differentiate between them regarding the protein activity.