

Detection and quantitative estimation of proteins by different methods

BCH303 [Practical]

Protein quantification:

- The accurate quantitation of protein content is a critical step in protein analysis.
- Importance of protein quantification ?
- Depending on the **accuracy required** and the **amount and purity** of the protein available:
 - ➔ different methods are appropriate for determining protein concentration.

Different methods of protein quantification:

- **Methods:**

1. Direct assay: measure the absorbance at 280 nm.
 2. Colorimetric/fluorescent and reagent-based protein assay: Protein is added to the reagent, producing a color change or increased fluorescence in proportion to the amount added.
- The most commonly used reagent-based techniques involve:
 - Biuret test.
 - Bradford test.
 - Bicinchoninic acid assay (BCA assay).
 - Lowry test.

Choosing the compatible method:

- **Best or ideal method ? WHY?**
- Each method has its **advantages** and **disadvantages**.
- How to choose the appropriate method?
 - ➔ Compatibility with the sample.
 - ➔ Availability.
 - ➔ Interfering substances .
 - ➔ Accuracy.
 - ➔ Sensitivity.
 - ➔ Time.
 - ➔

Choosing the compatible method:


| Method | Sensitivity | Time | Reagent | Interferences | Disadvantages and comments |
|---|-----------------------|----------------------|---|--|---|
| Biuret | Low 1-20 mg | Moderate 20-30min | Alkaline copper sulphate | Zwitterionic buffers, Some amino acids | Similar color with all proteins. Destructive to protein samples. |
| Lowry | High ~ 5 µg | Slow 40-60min | Cu ⁺² Folin– Ciocalteu | Ammonium sulphate, glycine, Zwitterionic, buffers, Mercaptans | Time-consuming. Color varies with proteins. Destructive to protein samples. |
| Bradford | High ~ 1 µg | Rapid 15 min | Coomassie Brilliant Blue G-250 | Strongly basic Buffers, detergents TritonX-100, SDS | Stable color, which varies with proteins. Reagent commercially available. Destruction to protein samples. Discoloration of glassware. |
| BCA | High ~ 1 µg | Slow 60 min | Cu ²⁺ , bicinchoninic acid | EDTA, DTT, Ammonium sulphate | Compatible with detergents. Reagents commercially available. Destructive to Protein samples. |
| Spectroph- otometric (A₂₈₀) | Moderate 50-100 µg | Rapid | - | Purines, pyrimidines, Nucleic acids | Useful for monitoring column eluent. Nucleic acid absorption can be corrected. None-destructive to protein samples. Varies with proteins. |

Criteria for choosing an assay:

- Therefore, successful use of protein assays involves selecting the method that is:
 - **Most compatible with the samples to be analysed, choosing an appropriate assay standard, and understanding and controlling the particular assumptions and limitations that remain.**


Determination of protein concentration:

- Protein concentration is determined by reference to a standard curve consisting of known concentrations of **a purified reference protein**.



- Next lab.**
- Typically, standard curves are constructed using at **least two replicates** for each point on the curve.

Determination of unknown concentration by standard curve:



Practical part

Experiment 1 : Qualitative detection of proteins by biuret test.

Objective:

- To detect the presence of a protein or peptides using biuret test.

Principle:

- In this reaction, peptide bonds in the proteins and peptides treated with an alkaline solution of dilute **copper sulphate** CuSO_4 (biuret reagent) forming a **purple coloured complex**.
- The colour density is proportional to the amount of proteins present.
- Two or more peptide bonds.
- Name?




Figure 1. The formation of biuret complex in biuret reaction.


Experiment 1 : Qualitative detection of proteins by biuret test.

Method:

1. Label three test tubes as **A and B**.
2. In tube **A**: add 1 ml of animal crude extract.
3. In tube **B**: add 1 ml of water.
4. Add 1 ml of biuret reagent to all tubes and mix well.

Results:

| Tube | Observation |
|----------------------|-------------|
| Animal crude extract | |
| Water | |



Blue color is the biuret reagent color

Experiment 2 : Quantitative estimation of proteins by biuret assay.

Objective:

- To determine the concentration of extracted protein by biuret assay.

Principle:

- Biuret method is based on copper ions Cu^{2+} binding to peptide bonds of protein under alkaline condition to give a violet colour that have a **maximum absorbance at 540 nm**.
- The intensity of the color, and hence the absorption at 540 nm, is directly proportional to the protein concentration, according to the Beer–Lambert law.




Figure 1. The formation of biuret complex in biuret reaction

From lower to higher concentration



There is a **linear relationship** between **purple color developed** and **concentration**.

Experiment 2 : Quantitative estimation of proteins by biuret assay.

Results:

Table 1. Concentration of standard BSA solution and their absorbance at 540 nm.

| Test tube | Protein concentration (g/L) [X- axis] | Absorbance at 540 nm [Y- axis] |
|---------------------------|--|-----------------------------------|
| Blank | | |
| A | | |
| B | | |
| C | | |
| D | | |
| E | | |
| F | | |
| G | | |
| Animal crude extract (D1) | _____ | |
| Animal crude extract (D2) | _____ | |
| Plant crude extract (D1) | _____ | |
| Plant crude extract (D2) | _____ | |




Figure 1. Standard curve of BSA using biuret method.

Experiment 3 : Quantitative estimation of proteins by Lowry assay.

Objective:

- To determine the concentration of extracted protein by Lowry assay.

Principle:

- Replaced by the more sensitive methods.
- The method is based on two chemical reactions.
- The resultant strong blue colour is partly dependent on the tyrosine and tryptophan content of the protein sample.




Figure 1. Series of reaction on Lowry method.