


Estimation of Serum Urea

-Urea:


- Urea is the highest **non-protein nitrogen** compound in the blood.
- Urea is the major excretory product of protein metabolism.
- It is formed by **urea cycle** in the **liver** from **free ammonia** generated during protein catabolism .
- Since historic assays for urea were based on measurement of nitrogen, the term **blood urea nitrogen (BUN)** has been used to refer to urea determination.



-Urea synthesis:

- **Protein metabolism** produces amino acids that can be oxidized.
- This results in the release of **ammonia** which is **converted to urea** (via urea cycle) and excreted as a waste product.
- Following synthesis in the **liver**, urea is carried out in the **blood to the kidney** which is readily **filtered from the plasma by glomerulus**.
- **Most** of the urea in the glomerular filtrate is excreted in the urine, and **some** urea is reabsorbed through the renal tubules.
- The **amount reabsorbed** depends on urine flow rate and **extent of hydration** (the amount of urea reabsorbed increases with dehydration.).
- **The concentration of urea in the plasma is determined by:**
 - ➔ Renal and liver function,
 - ➔ the protein content in diet,
 - ➔ and the rate of protein catabolism.

-Urea synthesis:



-Clinical Application:

- Measurement of urea used to in :
- Evaluate renal function.
- To assess hydration status.
- To determine nitrogen balance.
- To aid in the diagnosis of renal diseases.
- To verify adequacy of dialysis .
- Check a person's protein balance.

1-Plasma urea Concentration:

- Measurement of **Blood Urea Nitrogen (BUN)** alone is less useful in diagnosing kidney diseases because it's blood level is influenced by **dietary protein and hepatic function** (why?).
- But its diagnostic value improves with **serum creatinine values**.

	Type	Cause	Note
High urea (High urea concentration in plasma is called azotemia)	Pre-renal	<ul style="list-style-type: none">• Cognitive heart failure.• <u>Dehydration</u>.• High protein diet.• Increased protein catabolism.	<ul style="list-style-type: none">• Cognitive heart failure → reduced renal blood flow, less blood is delivered to kidney, then less urea is filtered.
	Renal	<ul style="list-style-type: none">• Renal failure .	
	Post-renal	<ul style="list-style-type: none">• Urinary tract obstruction.	
Low urea		<ul style="list-style-type: none">• Low protein intake.• Liver disease.• Pregnancy.	

2-Urine urea Concentration:

- The **Urine Urea Nitrogen test (UUN)** determines how much urea is in the urine to assess the amount of protein breakdown.
- The test can help determine **how well the kidneys are functioning**, and if the intake of **protein is too high or low**.
- Specimen: The urine urea nitrogen test is performed by collecting a **24-hour urine sample**.

	Cause
High urea in urine	<ul style="list-style-type: none">• Too much protein in the diet.• Too much protein breakdown in the body.
Low urea in urine	<ul style="list-style-type: none">• Malnutrition.• Too little protein in the diet.• Kidney issues.

Practical Part

-Objective:

- Estimation of blood urea nitrogen (BUN).

-Principle (of the kit used):

- The Reagent used contains: Urease, Glutamate Dehydrogenase, NADH, α -ketoglutaric acid, buffers and stabilizers .

1. **Reaction one:** Urea is hydrolyzed in the presence of urease enzyme and water to yield ammonia and carbon dioxide.



2. **Second reaction:** The ammonia reacts with α -ketoglutaric acid and reduced nicotinamide adenine dinucleotide (NADH) in the presence of glutamate dehydrogenase (GLDH) to yield glutamic acid and nicotinamide adenine dinucleotide (NAD).



- The amount of the urea in the sample is proportionally related to the reduced absorbance at 340 nm as a result of NADH oxidation to NAD.

-Reference Value:

SPECIMEN	UREA NITROGEN	UREA
Serum/Plasma	5-23 mg/dL	10-50 mg/dL
Urine 24 h	9-16g/24h	20-35 g/24 h

-Materials:

- BUN-ZYME Kit.

-Method:

	Standard	Serum
Reconstituted Reagent	3ml	3ml
Pre-warm at 37°C for 2 min. and add:		
Standard	0.025/25 μ l	-
Serum	-	0.025/25 μ l

- After exactly **30 seconds**, read and record absorbance A_1 against distilled water at **340 nm**.
- At exactly **60 seconds** after A_1 , read and record the absorbance A_2 and determine ΔA ($A_1 - A_2$).

-Calculations of the Results :

UREA NITROGEN (BUN)	UREA
SERUM OR PLASMA	
Urea Nitrogen (mg/dL) = $\frac{\Delta A (\text{Sample})}{\Delta A (\text{Standard})} \times 25$	Urea (mg/dL) = $\frac{\Delta A (\text{Sample})}{\Delta A (\text{Standard})} \times 53.57$

Homework:

- Determine the level of plasma Urea in the following cases and support your answer by some causes :
 - 1- liver disease.
 - 2- Diet high in protein.
 - 3- Dehydration.

-References:

- Clinical Chemistry: Techniques, Principles, Correlations (Bishop, Clinical Chemistry) Mar 31, 2009, by Michael L. Bishop MS MT (ASCP) CLS (NCA) and Edward P. Fody MD
- <http://www.nlm.nih.gov/medlineplus/ency/article/003605.htm>