(8): Separation of a mixture of phenols by thin layer chromatography (TLC)

The idea of the experiment:

Separation of a mixture of phenols (colorless compounds) by thin layer chromatography.

Observing the separated spots can be performed using:

- 1- UV light.
- 2- An iodine (I2) chamber. Iodine sublimes and will absorb to organic molecules in the vapor phase.

Materials and tools used:

Thin layer (a sheet of glass coated with silica gel). Phenols : phenol, catechol, pyrogallol. Unknown mixture, two mobile phases: 1- (hexane 5: 2ethylacetate), 2- (ethyl acetate 2: 5 dichloromethane), iodine.

Procedure:

- 1- Draw a line (in pencil not pen) across the bottom edge of the plate 1 cm up from the bottom.
- 2- Spot three spots along the line drawn on the plate.
- 3- Pour 10 ml of mobile phase in the jar and leave it few minutes to help to saturate the atmosphere with solvent vapor.
- 4- Put the plate inside the jar.
- 5- Remove the plate and mark the solvent front with a pencil.
- 6- Allow the plate to dry for a few minutes.
- 7- Calculate R_f for each substance.
- 8- Compare between R_f values of an unknown mixture and the known phenols.
- 9- Determine the components of an unknown mixture of phenols.
- 10- Repeat the same steps with another mobile phase and compare between R_f values of two mobile phases.