


Blood Biochemistry BCH 220 [Practical]

**Lab (2) Determination of Non-functional Plasma Enzymes
in Serum**



Objectives

- To determine the level of **lactate dehydrogenase (LDH)** in serum.
- To evaluate the presence of tissue damage.

Blood Enzymes

- Plasma, serum or **blood proteins**, are proteins present in blood plasma which have several functions.
- Some blood proteins also act as **enzymes**.
- **Enzymes** are biocatalysts that increase the rate of the chemical reaction.
- **Clinical enzymology** refers to measurement of enzyme activity in body fluids for the diagnosis and treatment of diseases.
- Most clinical enzyme measurements using **serum** or **plasma**, occasionally other fluids, such as urine and gut secretions are also investigated.
- The most commonly used body fluid for this purpose is **SERUM**. (Why?)

Differences Between Plasma Enzymes

Plasma Enzymes

1. Plasma-specific Enzymes (Functional)

Enzymes that are normally present in the plasma and perform their primary function in the blood.

2. Non-plasma specific Enzymes (Non functional)

Intracellular enzymes that are normally present in very small amount in blood and perform no known function in blood.


	Functional plasma enzymes	Non functional plasma enzymes
Their substrate	Always present in the blood	Absent from the blood
Site of synthesis	Liver	Different organs
Effect of diseases in its plasma levels	Decrease in liver diseases	Different enzymes increase in different organ diseases
Examples	Thrombin Plasmin	ALT LDH

 **Pause and Think** Which of these enzymes is a better diagnostic indicator ? Why?

Sources of Non functional Plasma Enzyme

1. **Cell damage** with the release of its content of enzymes into blood
2. **Block in the secretory pathway**
3. **Increase enzyme synthesis**

So estimation of the plasma concentration of these enzymes in blood is useful for the diagnosis of disease depending on their tissue origin.



Clinical Significance of Non-Functional Plasma Enzymes

Measurement of non-functional enzymes is important for:



- 1. Diagnosis of diseases.**
- 2. Prognosis of the disease:** following up of the treatment by measuring plasma enzymes before and after treatment.

Lactate Dehydrogenase (LDH)


- **Lactic acid dehydrogenase (LDH)** is an enzyme that helps produce energy.
- It is present in almost all of the tissues in the body and becomes elevated in response to cell damage
- LDH is most often measured to evaluate the presence of tissue damage (**diagnostic**).
- The enzyme LDH is in many body tissues, especially the **heart, liver, kidney, skeletal muscle, brain, blood cells and lungs**.

LDH Reaction

- LDH is a hydrogen transfer enzyme which catalyzes the **interconversion of pyruvate and lactate** with the mediation of **NAD⁺** as hydrogen acceptor, eventually **converting pyruvate to glucose**.
- The optimum pH for lactate pyruvate (L→ P) reaction is **8.8 – 9.8**, While for pyruvate to lactate (P→ L) is **7.7 – 7.8**.
- The enzyme is inhibited by **sulphydryl reagents** and **mercuric ions**.




Practical Part



Lactate Dehydrogenase Assay

Principle

LDH catalysis the following reaction:



- The rate of NAD^+ formation is indicated by **decrease in the absorbance at 340 nm** and it is directly proportional to serum LDH activity.

If:

- NADH is **reactant** → **decrease** the absorbance/min
- NADH is **product** → **increase** the absorbance/min

Method

	Tube
Sample (serum)	10 μ l
BUF (buffer/substrate)	1000 μ l
Mix and incubate at 37°C for 1-5 minutes	
SUB (substrate)	250 μ l
Mix, read the absorbance after 1 minute and at the same time start the stop watch. Read the absorbance again exactly after 1, 2 and 3 minutes	

Measure enzyme kinetics using UV-visible spectroscopy:

2) Applications → 2) Simple Kinetics → wave length (340 nm) → 1) Seconds → Duration (240 sec = 4 min) →

Intervals (60 sec= 1 min) → Print Data Table (off) → Press start (2 times) (*note: neglect the first reading*)

Results and Calculations

Time (min)		Absorbance at 340 nm
A ₁	1	
A ₂	2	
A ₃	3	
A ₄	4	

Calculations

- $\Delta A_1 = A_1 - A_2$ $\Delta A_2 = A_2 - A_3$ $\Delta A_3 = A_3 - A_4$
- $\Delta A/\text{min} = (\Delta A_1 + \Delta A_2 + \Delta A_3) / 3$
- $\text{LDH(U/L)} = \Delta A \times 20000$

Normal Values **225 to 450 U/L Adults**