


BCH 445- Biochemistry of Nutrition [Practical]
Estimation of inorganic phosphate in soft drinks



Phosphate in food

- Phosphate occurs **naturally** in the form of **organic esters** in many kinds of food, including meat, potatoes, bread, and milk.
- Phosphate also used as a food additive (inorganic phosphate) as a preservative, a flavor or color enhancer, extend shelf life, and retain moisture.



Phosphate

Figure 1. Chemical structure of phosphate

Soft drinks

- Soft drinks are complex mixtures containing a variety of substances such as coloring compounds, flavoring agents, acidifiers, sweeteners, preservatives, and caffeine.
- The most common **acidifier used in soft drinks** is phosphoric which gives a **tangy taste** in the mouth.
- Phosphoric acid can also acts as a **preservative, keeping the contents of the bottle fresh.**

Brand	Ingredients (as listed by the manufacturers)	Volume (mL)
A	Carbonated purified water, high-fructose syrup, sucrose, cane sugar, phosphoric acid , nature flavorings, caffeine	600

Figure 2. Ingredients of soft drinks

Coke vs Vinegar

Due to the use of phosphoric acid, cola **has similar acidity to vinegar** which no body can drink straight, but a ton of sugar, dyes and flavoring are added to **mask the acidity**.




Figure 3. pH comparison among pop-Cola drinks

Practical Part



Objective:

- Estimation of inorganic phosphate in soft drinks using ascorbic acid as a reducing agent.

Ammonium molybdate solution is prepared in **sulphuric acid**, **safety goggles** and **gloves** should be worn when handling the reagent

Principle

- Phosphoric acid is colorless, it **cannot be directly** determined using visible-light spectrophotometry, instead we'll **quantitatively** convert them into a colored substance, whose absorbance can be easily measured.
- Inorganic phosphate reacts with ammonium molybdate in an acid solution (ammonium molybdate prepared in sulphuric acid in this experiment) to form phosphomolybdic acid.
- Phosphomolybdic acid is then **reduced by a reducing agent (3% ascorbic acid)** to give **molybdenum blue a green/ blue color** that absorbs at 650nm .



Method

	Standard	Soft drink sample	Water	Ammonium molybdate	Ascorbic acid
Blank	-----	-----	2	0.5 ml	0.5 ml
3 ppm	2	-----	-----		
4.5 ppm	2	-----	-----		
6 ppm	2	-----	-----		
12 ppm	2	-----	-----		
15 ppm	2	-----	-----		
Sample (try different concentration)	-----	2	0		
		1	1		


Mix **thoroughly** after each addition → Allow to stand for 10 min
 (a deep blue/green color should develop) → Measure the absorbance at 650 nm.

Results

- Plot a graph between absorbance and concentration of phosphate in various standard solutions and obtain the calibrated curve.
- From the curve determine the amount of phosphate in the test solution.

Tube	Absorbance at 650 nm
Blank	
3 ppm	
4.5 ppm	
6 ppm	
12 ppm	
15 ppm	
Sample	

Estimation of inorganic phosphate in soft drink sample using ascorbic acid



Calculations

- **Inorganic phosphate concentration**= dilution factor x concentration from the curve = ----- ppm
- **Dilution factor**= final volume / aliquot volume

SD1=

SD2=

Note: the original sample was diluted 1:50

Homework

1. Enumerate 3 health effects of soft drinks.
2. How soft drinks consumption lead to osteoporosis?