


BCH 445- Biochemistry of Nutrition [Practical]
Lab (9) Determination of Calcium in Milk



Calcium

- **Calcium** is an important component of a healthy diet and a **mineral** necessary for life.
- **Calcium** is the most abundant mineral in the human body, making up 1.5 to 2% of the total body weight.
- It is a mineral that people need to build and maintain strong bones and teeth.
- It is also very important for other physical functions, such as **muscle control and blood circulation**.
- If we do not have enough calcium in our diets to keep our bodies functioning, calcium is removed from where it is stored in our bones.
- Overtime, this causes our bones to grow weaker and may lead to **osteoporosis** (a disorder in which bones become very fragile).



Milk and Calcium

- **Milk** is a heterogeneous mixture of **proteins, sugar, fat, vitamins** and **minerals**.
- Milk and milk products are some of the **natural sources of calcium**.
- The calcium concentration in bovine milk is about **1g/L**.
- Cow's milk has **good bioavailability** of calcium (21 to 45%).
- Milk is an excellent source of dietary calcium for those whose bodies tolerate it because it has a high concentration of calcium and the calcium in milk is excellently absorbed.
- In the United States, approximately 72% of calcium intakes come from dairy products and foods with added dairy ingredients.



Practical Part

Objective:

- Determination of Calcium in milk sample.

Principle

- In this experiment, The determination of calcium in milk is based on a **complexometric titration** of calcium with an aqueous solution of the disodium salt of EDTA at high pH value=12.
- **Complexometric titration** is a type of titration based on complex formation between the analyte and titrant.
- Such compounds (**EDTA**) are capable of forming **chelate complex** with many cations (**metal ion**) in which the cation is bound in a ring structure.
- The ring results from the formation of a salt-like bond between the cation and the carboxyl groups together with a coordinate bond through the lone pair of electrons of the nitrogen atom.




Figure 1. EDTA complex with metal ion

Principle

- The common form of the agent is **disodium salt $\text{Na}_2\text{H}_2\text{EDTA}$** .
- It is colorless and can be weighed and dissolve in water to form a stable solution.
- At high pH (>10) the remaining protons leave EDTA forming EDTA^{4-} anion.








Figure 2. Protonated and deprotonate form of EDTA



EDTA at low pH
(fully protonated)



EDTA at high pH
(Deprotonated)





Figure 3. EDTA reaction with calcium ion at a high pH


Solochrome indicator

- The Solochrome dark blue indicator is a suitable indicator in this case.
- The dye **itself has a blue color**.
- This blue dye also forms a complex with the calcium ions changing color from **blue** to **pink/red** in the process, but the dye–metal ion complex is less stable than the EDTA–metal ion complex.
- As a result, when the calcium ion–dye complex is titrated with EDTA the Ca^{2+} ions react to form a stronger complex with the EDTA changing the dye color to **blue**.



At the beginning of reaction

At the end of reaction



Excess Ca^{2+} ions present to complex with indicator

Ca^{2+} ions almost all complexed by EDTA

All Ca^{2+} ions complexed by EDTA, indicator completely un-complexed

Figure 4. Color changes for calcium-EDTA titration

How to determine calcium in the presence of Mg?

- This method for determining Ca^{2+} concentration in the presence of Mg^{2+} relies on the fact that the pH of the solution is sufficiently high ((The pH will be approximately 12.5 due to the addition of **concentrated NaOH solution**)) to ensure that all magnesium ions precipitate as magnesium hydroxide before the indicator is added.
- In this condition, **magnesium ions** are precipitated as hydroxide and do not interfere with the determination of calcium.




Figure 5. Calcium and magnesium complex with EDTA

Method

1. Combine 10mL of **sample**, 40mL distilled water, and 4mL of 8M sodium hydroxide solution into an Erlenmeyer flask and allow solution to stand for about 5 minutes with occasional swirling.
2. A small of magnesium hydroxide may precipitate during this time. Do not add the indicator until you have given this precipitate a chance to form.
3. Then add 6 drops of the **Solochrome dark blue** solution.
4. After that start to titrate with EDTA solution.
5. Repeat titration for three trials.

Results

	EDTA volume (ml)
1	
2	
3	
Average	

$$\text{Molarity (M)} = \frac{\text{number of moles}}{\text{liters of solution}}$$

$$\text{Number of moles (n)} = \frac{\text{Weight of substance}}{\text{molecular weight}}$$

Calculation

1. Calculate the moles of EDTA required to complex the Ca^{2+} ions in the sample:

→ Number of moles (for EDTA) = Molarity of EDTA x volume of EDTA in L

Note : Ratio $\text{Ca}^{2+}:\text{EDTA} = 1 : 1$ (i.e moles of EDTA = moles of Ca^{2+})

Molarity of EDTA = 0.03408 M

1. Calculate weight of Ca^{2+} :

→ Weight of Ca^{2+} = Number of moles x molecular weight (40.78)

• % of Ca^{2+} = (weight of Ca^{2+} / weight of sample) x 100

OR

- $\text{Ca}^{2+} \% = \frac{[\text{Molarity of EDTA X vol. of EDTA (in liter)] \times 40.78}{\text{Weight of the sample}} \times 100$

* 40.78 g/ mol is the molecular weight of Ca^{2+}

Additional resources

- <https://www.youtube.com/watch?v=AiWFA-Bvg60>
- <https://www.youtube.com/watch?app=desktop&v=hTy9JB1lUVg>