

Lab sheet #7

-Arginase Activity in Liver-

A. Method

1. Set two test tubes as 'test' and 'standard'
2. Add 1 ml of working reagent to each test tubes
3. Pre-incubate for 3 min. at 37°C
4. Add 0.01 ml of liver sample and standard to the 'test' and 'standard' respectively
5. After exactly 30 seconds, read and record absorbance A₁ against distilled water at 340nm.
6. At exactly 60 seconds after A₁, record the absorbance A₂ and determine ΔA (A₁ - A₂).

B. Results:

	Time (seconds)	Absorbance at 340 nm
Standard	30 (A ₁)	
	60 (A ₂)	
	ΔA (A ₁ - A ₂) =	
Serum	30 (A ₁)	
	60 (A ₂)	
	ΔA (A ₁ - A ₂) =	

C. Calculations:

$$A \text{ Urea (mmol/L)} = \frac{\Delta A \text{ sample}}{\Delta A \text{ standard}} \times \text{Std. Conc. (30)}$$

1. Urea produced in 1 minute = **A** mmol /min /0.01ml of liver extract.
2. Urea in 1ml of diluted liver extract = _____ x 100 = _____ mmol/min/ml of diluted liver extract.
3. Urea concentration in 1ml of undiluted liver extract = _____ x 3 x 5
4. Urea concentration in micromoles = _____ x 1000 = _____ micromoles/min/ml.
5. Total activity present in liver = _____ x total volume of liver extract (80ml) = _____ micromole.
6. **Arginase activity as (μmol/g) of liver** = total activity in liver/ wt of liver (20g).
_____ μmol/g