Lab sheet #6

Glutathione estimation in serum

Method:

- 1- Prepare 5 serial dilutions from GSH stock (10mg/dl) in concentrations (2, 4, 6, 8 and 10 mg/dl) with a final volume of 3 ml diluted with phosphate solution
- 2- Add 3 ml of phosphate solution to a separate tube and use it as a blank
- 3- Add 0.5 ml of the previously prepared tubes, then add 2, 0.25 and 0.25 ml of phosphate solution, DTNB and water respectively
- 4- In the sample tube add 0.3 ml of serum, then add 2, 0.25 and 0.45 ml of phosphate solution,
 DTNB and water respectively
- 5- Incubate all tubes for 10 min at 37°C
- 6- Read absorbance at 412 nm
- 7- Plot the standard curve and determine glutathione concentration from the graph expressed as μ g/ml

Result:

Tubes	Abs. at 412 nm
Blank	
1	
2	
3	
4	
5	
Sample	