



# PHYSICAL PROPERTIES AND DETECTION OF NORMAL CONSTITUENTS OF URINE

# Urinary system:

- The **kidneys, ureters, bladder** and **urethra** make up the urinary system
- The kidneys form **urine**, which passes through the ureters to the bladder for storage prior to excretion
- Waste products are excreted selectively, electrolyte levels are controlled and pH (acid-base balance) is maintained by excretion of hydrogen ions
- The composition of urine reflects exchange of substance between the nephron and blood in the renal capillaries




Figure 1. urinary system





Figure 2. nephron structure

# Urine Formation:

There are three processes involved in the formation of urine:

1. **Filtration**
2. **Selective reabsorption**
3. **Secretion**



# The three processes of urine formation





Figure 3. schematic representation of urine formation processes

# Step one: Filtration:

- This takes place through the **semipermeable membrane of glomerulus and glomerular capsule (Bowman's capsule)**
- **Water and small molecules** move from the glomerulus to the **inside** of the glomerular capsule
- Molecules which have molecular weight **more than** 70,000 Dalton can not pass the glomerulus
- **Blood cells, plasma proteins** and other large molecules are too large to filtrate (not filtrated)
- Inside the glomerular capsule now contains **glomerular filtrate** which is very similar in composition of plasma **except of plasma proteins and blood cells**
- (Non-selective filtration occurs).




Figure 4. Filtration

## Step two: Selective reabsorption:

- Is the process of restoration **water** and **some solutes** from the tubular fluid and returning them to the blood.
- **Reabsorption** is the movement of water and solutes from the tubule **back** into the blood.
- Nutrients such as glucose and amino acids **return** to the peritubular capillaries almost exclusively at the proximal convoluted tubule.

💡 Why is this step referred to as a selective reabsorption?




Figure 5. Selective reabsorption

## Step three: Secretion:

- Is a second way by which substances are removed from blood and added to the tubular fluid.
- Is a process in which the renal tubule extracts chemicals from the capillary blood and secretes them into the tubular fluid.
- Hydrogen ions (H<sup>+</sup>), creatinine, and drugs such as penicillin are some of the substances moved by active transport from blood into the kidney tubule.
- Tubular secretion is now known to occur along the length of the kidney tubule.








Figure 6. Secretion

In the end, **urine contains:**

1. Substances that have undergone glomerular filtration (**step one**) but have not been reabsorbed (**step two**).
2. Substances that have undergone tubular secretion (step three).



# Glomerular filtrate VS Urine

Constituent	Daily Excretion	
	Glomerular Filtrate	Urine
Water	130,000 ml 	1500 ml
Sodium	20,000 mmol 	150 ml
Albumin	4 g (60 $\mu$ mol) 	0.04 g (6 $\mu$ mol)
Urea	900 mmol 	400 mmol

100

# Composition of Normal Urine:

- Water 96%
- Urea 2%
- Uric acid
- Creatinine
- Oxalate
- Ammonia
- Sodium
- Potassium
- Chloride
- Phosphate
- Sulphate

2%





Figure 7. Composition of Normal Urine

# Urinalysis:

- **Urinalysis (UA)** is one of the most frequently ordered tests.
- Two unique characteristics of urine specimens:
  1. Urine is readily available and easily collected specimen.
  2. Urine contains information about many of the body's major metabolic functions, and this information can be obtained by simple laboratory tests.

## Laboratory testing for routine urinalysis (**types of testing**):

- First, the **physical characteristics** of the urine are noted and recorded.
- Second, a series of **chemical tests** is run. A chemically impregnated **dipstick** can be used for many of these tests.
- Third, the urine sediment is examined under **microscopic** to identify components.



Patient Name: \_\_\_\_\_

Age: \_\_\_\_\_  M  F

Physician's Name: \_\_\_\_\_

Collection Date: \_\_\_\_\_ Test Date: \_\_\_\_\_ Tester's Initials: \_\_\_\_\_

**Physical Characteristics:**

Color:  colorless  yellow  amber  other  orange  green  red

Appearance:  clear  hazy  cloudy  turbid

**Chemical Measurements:** (circle one)

urobilinogen (mg/dL)	normal	2	4	8			
glucose (mg/dL)	neg	50	100	250	500	1000	
ketone (mg/dL)	neg	trace/5	+ /15	++ /40	+++ /80	++++ /160	
bilirubin	neg		+	++	+++		
protein (mg/dL)	neg	trace		+ /30	++ /100	+++ /300	++++ /2000
nitrite	neg	pos (any pink color is considered positive)					
leukocytes	neg	trace	+	++	+++		
blood	neg	trace	moderate	trace	+ /small	++ /mod	+++ /large
		Non-Hemolyzed		Hemolyzed			
pH	5	6	6.5	7	8	9	
specific gravity	1.000	1.005	1.010	1.015	1.020	1.025	1.030

**Microscopic Examination:**

WBC \_\_\_\_\_ /HPF Crystals \_\_\_\_\_ Parasites \_\_\_\_\_

RBC \_\_\_\_\_ /HPF Bacteria \_\_\_\_\_ Spermatozoa \_\_\_\_\_

Casts \_\_\_\_\_ /LPF Yeast \_\_\_\_\_ Artifacts \_\_\_\_\_

Epithelial Cells \_\_\_\_\_ /HPF Other \_\_\_\_\_

Physical characteristics

Chemical tests

Microscopic examination

Figure 8. Urine Analysis Report

# Urine dipstick / Urine test strips:

- The **test strips** consist of a ribbon made of absorbent microfiber cellulose pads attached to it.
- Each pad contains the dried **reagents** needed for a **specific test** that react with the compounds present in urine producing a **characteristic colour**.
- The depth of color produced relates to the concentration of the substance in the urine.
- It provides quick **Semi-quantitative determinations** of pH, protein, glucose, ketones, bilirubin, hemoglobin (blood), nitrite, leukocyte, urobilinogen, and specific gravity.
- Color changes then matched to the control chart at the correct time after each stick is dipped into the urine specimen.




Figure 9. Urine dipstick

# Urine dipstick / Urine test strips:

## How to test your urine(visual read)?




**A**

**Prepare some fresh urine sample.**



**B**

**Dip the dry strip into the urine.**



**C**

**Absorb the excess urine with absorbent paper.**



**D**

**Contrast color chart, close to which color?**

Figure 10. Urine dipstick usage guide

How to use urine test strips : [https://youtu.be/g\\_4PiTezF6s](https://youtu.be/g_4PiTezF6s)

# Simple Examination of the Urine

## Physical Examination

Volume, Specific gravity,  
Color, Appearance, Odor, pH

## Chemical Examination

### Organic

Uric acid, Creatinine

### Inorganic

Chloride, Phosphate,  
Bicarbonate, Sulphate,  
ammonia.

# Physical Examinations:

## 1. Volume:

- The daily output of urine on an average diet and normal fluid intake is between **800-2500 ml** with an average of **1500 ml/day**.
- Effected by : 1) **Physiological factors** 2) **Pathological factors**

### Polyuria

- More than 2500 ml/day.
- Diabetes mellitus.
- Chronic renal insufficiency.

### Oliguria

- Below 500 ml /day.
- In case of deficient intake of water
- Excessive loss of fluids by other routs like hemorrhage or as diarrhea and vomiting.

### Anuria


- Below 100 ml /day.
- Stones or tumors in the urinary tract can also cause it by creating an obstruction to urinary flow.



# Physical Examinations:

## 2. Color

- Normally, Urine is **amber** in color due to the presence of urobilin (urochrome).
- **Pale** urine has a **low** specific gravity, a **dark** line has a **high** specific gravity ( a **direct relationship** between the colour and the specific gravity).
- **Coloured urines** occur in certain diseases or metabolic disorders, and after the administration of many drugs



The more that your body is dehydrated , the more the urine is dark

Figure 11. Urine color chart

## Change in the color or odour is a warning sign!




Figure 12. Urine color chart

# Physical Examinations:

## 3- Appearance:

- Normal urine is clear.
  - **Urine clarity is typically classified as:** clear, mildly cloudy, cloudy, or turbid.
- ➔ **Note:** cloudy or turbid urine can indicate dehydration, urinary tract infection or presence of RBCs, WBCs, epithelial cells or bacteria.



Figure 13. Urine Turbidity Scale

# Physical Examinations:

## 4- Odour:

- Normally urine smells **aromatic** due to the presence of volatile organic acids.
- The urine of patients with **diabetes ketoacidosis** may have a fruity (acetone) odour because of ketosis.
- Urine which is infected with **Gram-negative organisms** often has a distinctive unpleasant smell (fish-like odour).
- Certain drugs impart a typical odour.

## 5- pH:


- On a normal mixed diet the urine is **usually acid**, generally varying in pH between **5.5 and 8.0**, with a **mean of 6 in 24 hours**.

### Acidic Urine:

- Diabetic ketosis, urinary tract infection, diarrhea and starvation.

### Alkaline Urine:

- A **vegetarian diet** which causes a tendency to alkalosis.
- It may also be grossly increased by **bacterial infection** of the urinary tract.



# Physical Examinations:

## 6- Specific gravity (SG):

- The normal specific gravity (correctly called relative density) of a pooled 24 hour urine sample is between **1.010 and 1.025**.
  - There are **direct relationship** between **concentration** of substance in urine (concentration of urine) and **SG**.
- The concentration of urine is **highest** in the a morning specimen (**overnight urine**) and is **lowest** in a specimen passed an **hour after much fluid has been taken**.




Figure 14. Urine dipstick analysis

# Practical Part


---

# Objectives:

1. Simple **physical** examination of urine.
2. To detect some of the normal organic constituents of urine (**Qualitative**).
3. To detect some of the normal inorganic ions present in urine (**Qualitative**).

- **Note:**

All the examination in 24 hour collection of urine.



# Physical Examinations:

## Method:

### 1- Volume:

Measure the volume of the 24 hour collection of normal urine.

### 2- Odour:

State whether it is normal urine like ammonical, or not.

### 3- Colour:

Visually examine its colour.

### 4- Appearance:


State whether it is clear, cloudy or whether deposits or precipitates are present.

### 5- pH:

Record the pH of the sample by test strips.

### 6- Specific gravity:

Record the specific gravity of the sample by test strips.





# Chemical Examinations:

**Principle:** Each test based on the chemical properties of the substance.

## 1- Organic:

### A. Uric acid:


Method:

1. To 2 ml of urine add 1 ml of **Benedict reagen**.
2. Then heated in a boiling water bath for three minutes .
3. Changes to the **white precipitate** indicates the presence of uric acid.

### B. Creatinine:

Method:

1. To about 5 ml of urine add a few drops of a saturated solution of picric acid.
2. On rendering the solution alkaline with a few drops of 10% sodium hydroxide solution, a deep **red color** or orange due to creatinine picrate appears.
3. On acidification, with 2N HCl, the color changes to **yellow**.



# Chemical Examinations:

## 2- Inorganic

### A. Chloride:

-Method:

Add 5 ml of Urine +5 drops of 2N nitric acid+ 2N silver nitrate solution

→ A white precipitate of **silver chloride** is formed (Silver chloride is precipitated in the presence of nitric acid and silver nitrate).

### B. Phosphate:

-Method:

Add 5 ml of urine +5ml nitric acid+4 ml of ammonium molybdate, then heat in water bath.


→ A yellow crystalline precipitate of **ammonium phospho-molybdate** appears.

### C. Bicarbonate:

-Method:

Add 4 drops of concentrate hydrochloric +5 ml of urine.

→ A slight effervescence occurs due to **CO<sub>2</sub> evolution**. Test the gas evolved with lime water.



AM

# Chemical Examinations:

## 2- Inorganic

### D. Sulphate:

-Method:

To Acidify add 10 ml of urine with 1ml dilute hydrochloric acid + 4 drops of 5% barium chloride solution


→ A white precipitate sulphate is precipitated as of **barium sulphate** is formed.

### E. Ammonia:

-Method:

Add 1 ml of 10% sodium hydroxide solution +5 ml or urine, then heat in water bath.

→ The evolved **ammonia** may be detected its occur in confirmed by turning the moist red litmus paper to blue.



AM

# Summery:

Physical examination	
The normal constituent of 24 hour urine	
<b>Volume</b>	800-2500 ml with an average of 1500 ml
<b>Color</b>	Amber in color
<b>Appearance</b>	Clear
<b>Odour</b>	Urine like
<b>pH</b>	5.5 - 8.0, with a mean of 6
<b>Specific gravity</b>	1.010 - 1.025
Chemical examination	
Chemical	Positive result
<b>Uric acid</b>	White precipitate
<b>Creatinine</b>	Deep orange color
<b>Chloride</b>	White precipitate
<b>Phosphate</b>	Yellow precipitate
<b>Bicarbonate</b>	CO <sub>2</sub> bubble appeared
<b>Sulphate</b>	White precipitate
<b>Ammonia</b>	Litmus paper turns to blue