


# DETECTION AND ESTIMATION OF SOME ABNORMAL CONSTITUENTS IN URINE



# Abnormal urine

## Physical Examination

Odor

Volume

Color

pH

Specific gravity

## Chemical Examination

Blood (RBC)

Leukocyte (WBC)

Ascorbic acid

Glucose

Nitrite

Ketone bodies


Proteins

Amino acids

Bilirubin


Uroblinogen

# 1- Physical Examination (abnormal):




## 2- Chemical Examination

- The following are some abnormal constituents that **not normally** found in **detectable** amount:

| Positive in Urine                 | Cause  | Notes   |
|-----------------------------------|--|---|
| <b>Blood (RBC)</b><br>(hematuria) | <ul style="list-style-type: none"><li>Bleeding because of damage to kidney or genitourinary system, eg: <b>Renal Calculi, Renal Tumor, Trauma to kidneys.</b></li><li>Urinary tract infection.</li><li>Malignant hypertension.</li></ul> | <ul style="list-style-type: none"><li>Any pink, red or brown urine must be considered as bloody until proved otherwise.</li></ul>  |
| <b>Hemoglobinuria</b>             | <ul style="list-style-type: none"><li>Intravascular hemolysis due to hemolytic anemia.</li></ul>   |   |
| <b>Leukocyte (WBC)</b>            | <ul style="list-style-type: none"><li>Urinary tract infection bacteria.</li></ul>  | <ul style="list-style-type: none"><li>Urine with positive results from the dipstick should be examined microscopically for WBCs and bacteria.</li></ul>   |
| <b>Ascorbic acid</b>              | <ul style="list-style-type: none"><li>Large urinary concentrations arise from therapeutic doses of vitamin C.</li></ul>  |   |

## 2- Chemical Examination

| Positive in Urine                   | Cause   | Notes  |
|-------------------------------------|---|--|
| <b>Glucose</b><br>(Glycosuria)      | <ul style="list-style-type: none"><li>Blood glucose level exceeds the reabsorption capacity of the tubules, eg, <b>Diabetes mellitus</b>.</li><li>Defect in the tubular reabsorption eg. <b>fanconi syndrome</b>.</li></ul> | Normally, glucose is present in the glomerular filtrate and reabsorbed by the proximal tubules. ( <i>see next slide</i> )  |
| <b>Ketone bodies</b><br>(ketonuria) | <ul style="list-style-type: none"><li>Occurs whenever increased amounts of fat are metabolized eg, <b>Diabetes mellitus, starvation and altered carbohydrate metabolism</b>.</li></ul>                                      | Urine may have a fruity (acetone) smell.   |
| <b>Nitrite</b>                      | <ul style="list-style-type: none"><li>Urinary tract infection.</li></ul>  | Bacteria that can reduce the nitrate to nitrite.<br> <pre>graph TD; A[Nitrate (NO3-)] -- Nitrate reductase --&gt; B[Nitrite (NO2-)]</pre> |

# Note:

Glucose level exceeds the reabsorption capacity in diabetes patients:




Figure 1. Glucose reabsorption by the kidney


💡 How are high glucose levels related to polyuria?

## 2- Chemical Examination

| Positive in Urine                    | Cause  | Notes   |
|--------------------------------------|--|---|
| <b>Bilirubin</b>                     | <ul style="list-style-type: none"> <li>Elevated amount of bilirubin in the blood stream, eg, <b>Bile duct obstruction.</b></li> </ul>  | <ul style="list-style-type: none"> <li>The urine may be dark with a <u>yellowish foam</u> if much is present.</li> </ul>        |
| <b>Urobinogen</b>                    | <ul style="list-style-type: none"> <li>Increased production eg, <b>hemolytic anemia.</b></li> </ul>  | <ul style="list-style-type: none"> <li>Its presence does <u>not</u> give a colored foam (urobilinogen is colorless).</li> </ul> |
| <b>Amino acid</b><br>(aminoaciduria) | <ul style="list-style-type: none"> <li>Blood amino acid level exceeds the reabsorption capacity of the tubules eg, <b>Phenylketonuria, Alkaptonuria.</b></li> <li>Defect in the tubular reabsorption eg, <b>fanconi syndrome, cystinuria.</b></li> </ul> |   |
| <b>Protein</b>                       | <ul style="list-style-type: none"> <li>Acute infection.</li> <li>Primary kidney disease.</li> <li>Secondary kidney disease.</li> </ul>   |   |

## Common Correlations in Urinalysis:

| <b>Microscopic Elements</b> | <b>Physical Examination</b>   | <b>Dipstick Measurement</b>      |
|-----------------------------|-------------------------------|----------------------------------|
| Red blood cells             | Turbidity, red to brown color | Blood                            |
| White blood cells           | Turbidity                     | Protein<br>Nitrite<br>Leukocytes |
| Epithelial cast cells       | Turbidity                     | Protein                          |
| Bacteria                    | Turbidity, odor               | pH<br>Nitrite<br>Leukocytes      |
| Crystals                    | Turbidity, odor               | pH                               |





## Test strip (dipstick):


- Normally, substances such as nitrite, proteins, glucose, ketone bodies, bilirubin, urobilinogen and blood are present in **very small quantities** that is **not** capable of detection by this method.

→ But present in detectable amount are **not normal**.

**(False positive and false negative are common when using dipstick)**



|                    | False-positive            | False-negative         |
|--------------------|---------------------------|------------------------|
| <b>Protein</b>     | Alkaline Urine<br>Ammonia | Dilute Urine           |
| <b>Glucose</b>     | Strong oxidizing agent    | Ascorbic acid          |
| <b>Blood</b>       | Oxidizing contaminants    | High ascorbic acid     |
| <b>Bilirubin</b>   | Certain drugs             | Ascorbic acid, nitrate |
| <b>Uroblinogen</b> | Alkaline Urine            | Nitrite, formaline     |
| <b>Nitrite</b>     | Pigmented urine           | Ascorbic acid          |



## Notes in using test strip:

- Reagent strips should be stored in their original container.
- The lid should be kept tightly closed.
- Strips should **not** be used if expired or discolored.
- Strips should **not** be exposed to sunlight, moisture, heat, or cold.
- Specific reagents should be read at the appropriate time after dipping in urine, as recommended by the manufacturer.
- The strip should **not** be dipped for more than a second in the urine, and excess urine should be blotted off on the edge of absorbent paper to prevent mixing of reagents.

# Types of urine specimens:

- Type of specimen and collection procedure are determined by physician and depend on the tests to be performed.
- There are basically four types of urine specimens:

| Sample type           | Sampling  | Purpose  |
|-----------------------|---|--|
| Random specimen       | No specific time<br>most common, taken<br>anytime of day  | Routine screening, chemical  |
| Morning sample        | First urine in the morning,<br>most concentrated  | Pregnancy test, microscopic<br>test                                |
| Clean catch midstream | Discard first few ml, collect<br>the rest   | Culture  |
| 24 hours              | All the urine passed during<br>the day and night and next<br>day 1 <sup>st</sup> sample is collected. | used for quantitative and<br>qualitative analysis of<br>substances |

- **Note:** 24h sample is necessary for accurate quantitative results.

# Practical Part

---

# Objectives:

1. The semi-quantitative detection of some abnormal constituents using **test-strips**.
2. The detection of amino acids in a urine sample using **ninhydrin**.
3. The effect of the type of urine collection in the detection of urine constituents.



# 1- Detection of some abnormal constituents of urine using test strip:

## -Method:

1. You will have one urine samples.
2. You have to fill the following table and then the probable diagnosis.

| Test                             | Sample 1 |
|----------------------------------|----------|
| Volume                           | 3000 ml  |
| Color                            |          |
| Odor                             |          |
| pH                               |          |
| Specific gravity                 |          |
| Protein                          |          |
| Blood                            |          |
| Bilirubin                        |          |
| Uroblinogen                      |          |
| Glucose                          |          |
| Ketone                           |          |
| Nitrite                          |          |
| Leukocyte                        |          |
| Clinical Diagnosis for sample 1: |          |

## 2- Detection of amino acid using ninhydrin:


### Principle :

Ninhydrin reacts with **all amino acids** except proline and hydroxyproline at pH 3-4 to give a **purple colored** compound.

→ **Proline** will give a **yellow color**.

1. Initially, the amino acid is oxidized to an aldehyde containing one carbon atom less, together with the release of **ammonia** and **carbon dioxide**.

2. Then the **ammonia**, **ninhydrin** and the reaction product **hydrindantin** react to form the **purple product**.




## Method:

- As standards, use proline and glycine as the following table:

| Solution     | Volume (ml) |
|--------------|-------------|
| Glycine      | 1           |
| Proline      | 1           |
| Urine Sample | 1           |

- Add a 1 ml of ninhydrin solution to each test-urine.
- Boil the contents of each test tube for 2 minutes.
- Record your observations.






### 3- The effect of the type of urine collection on the detection of Urine constituents:

#### -Method:

1. You have two samples, one is random urine sample, the other is 24-hour urine sample from the **same patient**.
2. Compare between the two samples in the presence of the proteins using the test strip.

| Test Parameter      | 24 hour Urine sample | Random urine Sample |
|---------------------|----------------------|---------------------|
| Protein<br>(+ or -) |                      |                     |



## References:

1. A Manual of Laboratory and Diagnostic Tests 9th edition (January , 2014), Frances T Fischbach RN, BSN, MSN By Lippincott Williams & Wilkins Publishers
2. The Abnormal Urinalysis, Hiren P. Patel, MD *Pediatr Clin N Am* 53 (2006) 325 – 337
3. Clinical Biochemistry, An Illustrated Colour Text 4th edition, Allan Gaw, Michael J. Murphy, Robert A. Cowan, Denis St. J. O'Reilly, Michael J. Stewart, James Shepherd
4. BCH 472 BCH practical note