


QUANTITATIVE PROTEIN ESTIMATION OF URINE



Introduction

- In a **healthy renal and urinary tract system**, the urine contains no protein or only trace amounts.
- The presence of **increased amounts of protein** in the urine can be an important indicator of renal disease. It may be the first sign of a serious problem and may appear before any other clinical symptoms.
- However, there are other **physiologic conditions** (eg, exercise, fever) that can lead to increased protein excretion in urine.
- Also, there are some renal disorders in which proteinuria is **absent**.




Figure 1. Illustration of proteinuria

Proteinuria



- Proteinuria means **protein in urine**.
- Protein in normal urine should be **less than 150 mg/L**.
- **Proteinuria** is defined as urinary protein excretion of greater than 150 mg per day (per one liter).
- High concentrations of protein cause frothy or sudsy urine.
- Although proteinuria is usually benign, the condition can be a marker for a serious underlying renal disease or systemic disorder.

Notes:

1. **Dipsticks** (is the most common **initial** screening test for proteinuria) can only detect around 150 mg/L of albumin.

→ The dipstick will not detect non-albumin proteins.

1. /L = /24-hour = /day [**for normal, because the average normal output per day is 1000 ml or 1L**]

Types of Proteinuria as related to the cause:

Some of the causes of proteinuria are:

1. **Primary** kidney disease called **nephritis**.
2. **Secondary** kidney disease due to:
 - Diabetes
 - Hypertension
 - Cancer
 - Pregnancy
 - Medications
3. **Other:**
 - Acute infection

Mechanisms of proteinuria

Type

Glomerular proteinuria (Filtration)

- **Cause:** due to glomerular disease and abnormal permeability of the glomerular capillaries to protein lead to increased filtration of normal plasma protein and because **albumin** has the highest concentration in the plasma it is called **abuminuria** eg. Malignant hypertension, Glomerulonephritis.
- *The most common kind of proteinuria.*

Tubular proteinuria (Reabsorption)

- **Cause:** Defect in the *reabsorption* process eg, **Fanconi Syndrome**.
- **Signs:** Low molecular weight protein that is found in urine eg, beta-2 microglobulin.
- **Note:** The amount of proteinuria is < 2 g and dipstick may be negative.

Secretory proteinuria (Secretion)

- **Cause:** Over-secretion of certain proteins in the tubules, most notably the over secretion of Tamm-Horsfall proteins eg, in **interstitial nephritis**.

Mechanisms of proteinuria

Type	
Overflow proteinuria	<ul style="list-style-type: none">▪ Cause: Commonly associated with <u>increased production</u> of abnormal low molecular weight proteins (eg, light chains in multiple myeloma) i.e overflow of high plasma proteins, that exceeds the reabsorption capacity of the tubules.▪ Signe: concentrations of low molecular weight protein found in urine.
Functional proteinuria	<ul style="list-style-type: none">▪ Cause: Occurs when <u>increased renal blood flow</u> delivers increased amounts of protein to the nephron, eg. exercise, fever, high-output heart failure, resulting in increased protein in the urine.▪ Note: usually < 1 g/day and reverses when renal blood flow returns to normal.

Mechanisms of proteinuria




Figure 2. Mechanisms of proteinuria

Types of Proteinuria as Related to Quantity:

Protein amount per 24-hour	Type of proteinuria
0.15 to 2.0 g	<ul style="list-style-type: none">▪ Tubular proteinuria.▪ Overflow proteinuria (<i>an increased proportion of low molecular weight proteins</i>).
2.0 to 4.0 g	<ul style="list-style-type: none">▪ Usually glomerular.
> 4.0 g	<ul style="list-style-type: none">▪ Always glomerular (<i>mainly albumin</i>).

Protein determination

- The quantitative estimation of the daily excretion of protein is of value to the clinician in order to give a general idea of the type of renal disease (HOW?) its severity and to monitor the results of treatment given.
- The protein content can be determined by numerous methods, in this lab **turbidimetric method** will be used.
- Determination of total protein by measurement of **protein turbidity produce by mixed with an anionic organic acid** such as sulfosalicylic acid , TCA , or benzethonium chloride.
- **Sulphosalicylic acid** is used in this experiment to precipitate the protein in a 24 hour sample of urine.
- The turbidity is proportional to the concentration of the protein, and may be measured with a spectrophotometer at 500 nm.

Practical Part

Sulfosalicylic acid (SSA) test:

- **The sulfosalicylic acid (SSA)** turbidity test quantitatively screens for proteinuria.
- The advantage of this easily performed test is its greater sensitivity for proteins such as Bence Jones.
- The SSA reaction will detect globulin and Bence-Jones proteins, in addition to albumin (*although it is more sensitive to albumin*).




Figure 3. Sulfosalicylic acid structure




Figure 4. Protein turbidity test results

Sulfosalicylic acid (SSA) test:

Objective:


- Quantitative estimation of protein in urine by turbidimetric methods using sulfosalicylic acid.

Principle:

- Sulfosalicylic acid is an **anion(-)** which neutralizes the protein **cations(+)** leading to its precipitation (pH in highly acidic media, the protein will be positively charged, which is attracted to the acid anions that cause them to precipitate).
 - Then the radiation of a wavelength which is not absorbed by the solution is made to pass through the suspension and the apparent absorption will be solely because of the scattering by the particles.
- (The higher protein concentration , the lower transmittance value).



Increased concentration
Decreased transmission



Method:

1-Set up a series of test tube as follows, label from 1- 7:

Tube	Protein Stock Solution(140 mg/dl)	0.85% Saline
1	4.5	1.5
2	3	3
3	2.4	3.6
4	1.5	4.5
5	0.9	5.1
6	0.3	5.7
7(Blank)	0	6
Urine Sample	-	-

2-Set another 8 test tube labelled 1-7 and pipette in each one add 8 ml of sulfosalicylic acid :

Tube	Sulfosalicylic acid
1	8 ml
2	8 ml
3	8 ml
4	8 ml
5	8 ml
6	8 ml
7(Blank)	8 ml
Urine Sample	8 ml


- 3- Into tube 1 pipette 2 ml of protein solution 1 (that you prepared before), into tube 2 pipette 2 ml of protein solution 2 etc. For the Urine Sample pipette 2ml of the Sample.
- 4-Mix contents of each tube well and allow standing for 5 minutes.
- 5-Using solution 7 (Blank) to set transmittance 100% at 500nm.
- 6-Then use solutions from 1-6, to recorded respective transmittance of each suspension.
- 7-Record your results.

Tube	Transmittance at 500 nm	Protein concentration mg/dl
7(Blank)	100 %	
1		
2		
3		
4		
5		
6		
Urine Sample		

Results:

- Plot Transmittance against Protein concentration (mg/dl).
- Read the Protein concentration of Urine Sample from the standard curve.
- Compare the result you got with the normal range of protein excretion in 24 h urine specimen if you know that the protein excretion in healthy sample (0- less than 0.150g/24 h).

Note: Assuming that the 24 hour urine sample for the patient = 1000 ml



References:

- The Washington Manual™: Nephrology by Steven Cheng, Anitha Vijayan, Katherine E. Henderson, and Thomas M. De Fer 3rd edition published by Wolters Kluwer (<https://www.inkling.com/read/washington-manual-nephrology-cheng-vijayan-3rd/chapter-4/approach-to-proteinuria>).
- A Manual of Laboratory and Diagnostic Tests 9th edition (January , 2014), Frances T Fischbach RN, BSN, MSN By Lippincott Williams & Wilkins Publishers.
- Clinical Biochemistry, An Illustrated Colour Text 4th edition, Allan Gaw, Michael J. Murphy, Robert A. Cowan, Denis St. J. O'Reilly, Michael J. Stewart, James Shepherd.
- Zandi-Nejad K, Eddy AA, Glasscock RJ, Brenner BM. Why is proteinuria an ominous biomarker of progressive kidney disease? *Kidney Int Suppl.* 2004(92):S76–S89.
- Eric P Cohen, MD Professor, Department of Medicine, Division of Nephrology, Medical College of Wisconsin; Nephrology Section Chief, Zablocki Veterans Affairs Hospital (Nephrotic Syndrome) <http://emedicine.medscape.com/article/244631-overview>
- BCH 472 practical note