


ESTIMATION OF SERUM CREATININE, URINE CREATININE AND CREATININE CLEARANCE



Kidney functions

The kidneys serve three essential functions:

1. **They function as filters**, removing metabolic products and toxins from the blood and excreting them through the urine.
2. **They regulate the body's** fluid status, electrolyte balance, and acid-base balance.
3. **The kidneys produce or activate hormones** that are involved in erythropoiesis, Ca^{2+} metabolism, and the regulation of blood pressure and blood flow.



Renal function tests (RFT)

- Are used to detect the presence of **renal diseases** and assess their progress.
- The most widely used test is to measure the **glomerular filtration rate (GFR)**, that is, the rate of filtrate formation by the kidneys.

Glomerular Filtration Rate (GFR)

- Under **normal** conditions, approximately 625 mL of plasma flow through the kidneys each minute and the volume of plasma **filtered** is 125 ml/ min which is called the **glomerular filtration rate**.
- Glomerular filtration rate (GFR): is the volume of plasma filtered by the kidneys per unit of time.
- GFR is an **important** and the **best** overall measurement in the evaluation of kidney function.

Measuring the GFR

- Accurate measurement of the GFR by clearance tests requires determination of the **concentration**, in **plasma** and **urine**, of a substance is known to be **completely filtered** from the plasma at the glomerulus.
- This substance must not be reabsorbed nor secreted by renal tubules, broken down, or accumulated by the tubules and must remain at a constant concentration in the **plasma** throughout the period of urine collection.
- **Clearance is given by:**

$$\text{Clearance (ml/min)} = U.V / P$$

→ Where:

U= concentration of any substance in urine.

P= concentration of the same substance in plasma.

V= volume of urine (ml/min).

Substances used for Measuring GFR (clearance test)

| | Inulin Clearance | Creatinine Clearance | Urea Clearance |
|----------------------|---|--|--|
| Source | Non-toxic fructose polymer. | End-product of skeletal muscle creatine metabolism. | End-product of protein metabolism. |
| Advantages | Not reabsorbed or secreted. | An <u>endogenous</u> product of muscle metabolism; near constant production. | An <u>endogenous</u> product of protein. |
| Disadvantages | Not made by body; must be injected (exogenous). | Small amount is secreted. | -Partially reabsorbed. -Synthesis varies with diet. |

Creatinine

- **Creatinine** is derived from “**creatine**” which is synthesized in the liver, kidney and pancreas it moves through the circulation and is taken up entirely by muscles.
- Creatinine is a substance that, in health, is easily excreted by the kidney.
- It is the by-product of muscle energy metabolism and is produced at a **constant rate according to the muscle mass of the individual**.
- Endogenous creatinine production is **constant** as long as the muscle mass remains constant.




Figure 1. Creatinine reaction

Creatinine clearance

- A measure of the amount of creatinine eliminated (filtered) from the blood by the kidneys.
- Creatinine is cleared from the body fluids **almost entirely** by glomerular filtration (*small amount is secreted by kidney tubules*).
- Therefore, the clearance of creatinine can be used to assess GFR.
- Because measurement of creatinine clearance does not require intravenous infusion into the patient, this method is much more widely used than inulin clearance for estimating GFR clinically.
- **Tubules to variable degree secrete creatinine**, which by itself, would lead to an ~20% overestimate of GFR in humans.

Serum Creatinine

- **High plasma creatinine:** plasma creatinine tends to be higher in subjects with a **large muscle mass**.

Other **non-renal** causes of increased plasma creatinine include the following:

- A high meat intake can cause a temporary increase.
- Transient, small increases may occur after vigorous exercise.

→ If non-renal cause can be excluded, an increased plasma creatinine indicates a fall in GFR (renal disease).

Urine Creatinine

Decreased urine creatinine is found in:

- Advanced renal disease.
- Renal stenosis.

Increased urine creatinine is found in:

- Diabetes mellitus.

Clinical Implications of creatinine clearance

Decreased creatinine clearance is found in any condition that decreases renal blood flow:

- Impaired kidney function.
- Shock, dehydration.
- Haemorrhage.
- Hypothyroidism

Increased creatinine clearance is found in:

- Pregnancy.
- Hyperthyroidism

Reference Values

- Urine creatinine :1- 2 g/24h
- Serum creatinine: 0.6–1.2 mg/dL
- Creatinine clearance: 100-130 ml/min/1.73m²

Note:

What 1.73 m² means?

- Kidney function is proportional to kidney size, which is proportional to body surface area. A of 1.73 m² is the normal mean value for young adults.
- Adjustment for body surface area is necessary when comparing a patient's estimated GFR to normal values or to the GFR criterion for the diagnosis of CKD, and to levels defining the stages of CKD.

Chart 1. Chronic kidney disease staging

| Stage | Description | GF (ml/min/1.73m ²) |
|-------|---|---------------------------------|
| I | Kidney lesion with normal or increased GF | ≥ 90 |
| II | Kidney lesion with mild GF decrease | 60-89 |
| III | Kidney lesion with moderate GF decrease | 30-59 |
| IV | Kidney lesion with marked GF decrease | 15-29 |
| V | Functional kidney failure or undergoing SRT | < 15 |

SRT- substitutive renal therapy. Source: National Kidney Foundation, 2002.

Practical Part


Objectives:

- To estimate creatinine in serum and urine.
- To calculate creatinine clearance value.

Principle:

(Jaffe's method):

Colorimetric estimation of creatinine using the alkaline picrate method:



Absorbance at 520nm

Method

1- Set up a series of 8 test tube as follows:

| Chemical | Standard (3mg/dl) | | Test (serum) | | Test (urine) | | Blank |
|--------------|-------------------|--------|--------------|--------|--------------|--------|-------|
| | (A) | (B) | (C) | (D) | (E) | (F) | |
| Water | 1.5 ml | 1.5 ml | 1.5 ml | 1.5 ml | 1.5 ml | 1.5 ml | 2 ml |
| Standard | 0.5 ml | 0.5 ml | - | - | - | - | - |
| Serum Sample | - | - | 0.5 ml | 0.5 ml | - | - | - |
| Urine Sample | - | - | - | - | 0.5 ml | 0.5 ml | - |
| Picric acid | 6 ml | 6 ml | 6 ml | 6 ml | 6 ml | 6 ml | 6 ml |

2- Immerse the Tubes carefully in the boiling water bath for 40 seconds.

3- Pipette 0.6 ml of NaOH to all tube.

4- Let the tubes stand for 20 min.

Results:

| Tube | Standard | | Test (Serum) | | Test (Urine) | |
|---------------------------------|----------|-----|--------------|-----|--------------|-----|
| | (A) | (B) | (C) | (D) | (E) | (F) |
| Absorbance at 520 nm | | | | | | |
| Average (Mean of Absorbance) | | | | | | |

Calculation

Patient information: 24h urine volume = 100ml, gender: women, body surface: 1.6m², DF=100.

1- Serum creatinine=

(Mean Absorbance of serum test ÷ Mean Absorbance of Standard) X concentration of standard = mg / dl

2- Urine creatinine=

(Mean Absorbance of Urine test ÷ Mean Absorbance of Standard) X concentration of standard X **DF**= mg / dl

(To compare with normal range, convert from mg/dl to g/24 h) **HOW?**

3- Creatinine Clearance:

=U.V/ P

= [(*Urinary creatinine (mg/dl)*) / (*plasmac creatinine (mg/dl)*)] x *Urine volume(ml/min)* = B

B -----> 1.6 m² (person surface area)

? -----> 1.73 m²

Corrected for surface area= ml/min/1.73 m²

Example:

Find the Creatinine Clearance if you know that the Urine creatinine $U = 488 \text{ mg/dl}$, Serum creatinine $P = 2.32 \text{ mg/dl}$, Volume of urine in $24 \text{ h} = 100 \text{ ml}$ and A "surface area" $= 1.6 \text{ m}^2$?

$$\begin{aligned}\rightarrow \text{Creatinine Clearance} &= U.V / P \\ &= (488 \text{ mg/dl} \div 2.32 \text{ mg/dl}) \times (100 \div 1440^*) = 14.6 \text{ ml/min}\end{aligned}$$

14.6 ml/ min in 1.6 m², find the creatinine clearance for 1.73 m² surface area :

$$= (14.6 \times 1.73) \div 1.6 = \underline{15.8 \text{ ml/min/1.73m}^2}$$

-----OR-----

$$\begin{aligned}\rightarrow \text{Creatinine Clearance} &= (U \times V \times 1.73) / (P \times 1440 \times A) \\ &= (488 \text{ mg/dl} \times 100 \times 1.73) / (2.32 \times 1440 \times 1.6) \\ &= \underline{15.8 \text{ ml/min /1.73m}^2}\end{aligned}$$

* To convert 24 hour to min (24x60 = 1440)

Discussion

- Comment on the concentration of creatinine in serum and urine.
- Comment on the value of Creatinine Clearance.

Homework

A man aged 40 years with surface area of 1.8m^2 has a serum creatinine of 6 mg/dl . A 24 h urine of 2200 ml is collected and the urine creatinine concentration found to be 450 mg/dl .

Calculate the Creatinine Clearance.

References

- Kidney function tests: MedlinePlus Medical Encyclopedia. (n.d.). Retrieved from <https://medlineplus.gov/ency/article/003435.htm>
- (N.d.). Retrieved from [https://www.mayoclinicproceedings.org/article/S0025-6196\(11\)60609-5/fulltext](https://www.mayoclinicproceedings.org/article/S0025-6196(11)60609-5/fulltext)
- Lecture Notes: Clinical Biochemistry Geoffrey Beckett, Simon W. Walker, Peter Rae.
- A Manual of Laboratory and Diagnostic Tests ,By Frances T Fischbach RN, BSN, MSN By Lippincott Williams & Wilkins Publishers.