


ESTIMATION OF URIC ACID IN SERUM



Uric acid production

1. **Uric acid** is the product of **catabolism of the purine** (adenosine and guanine) that result from the break down of ingested nucleic acid (*exogenous*) or from tissue destruction (*endogenous*).
2. **Uric acid** is transported by the plasma from the **liver** to the **kidney**, where it is filtered and where about **70% is excreted**. The remainder of uric acid is excreted into the **GI tract** and degraded.




Figure 1. Chemical structure of Uric acid




Figure 2. uric acid production and excretion

Uric acid excretion




Figure 3. Uric acid excretion

Purine metabolism




Figure 4. Pathways of purine metabolism, Source: sciencedirect.com

Clinical Application

1. Uric acid is measured to assess inherited disorders of **purine metabolism**.
2. To conform diagnosis and monitor treatment of **gout**. ??
3. To assist the diagnosis of **renal calculi** (uric acid kidney stones).
4. To detect **kidney dysfunction**.
5. Evaluation of **leukemia**.

Serum Uric Acid

Case	Cause
Increased uric acid serum (Elevated uric acid levels (hyperuricemia))	<ul style="list-style-type: none">• Gout (the amount of increase is <u>not</u> directly related to the severity of the disease).• Renal diseases and renal failure (decreased excretion of uric acid)• Leukemia, multiple myeloma, lymphoma.• Lesch-Nyhan syndrome (rare hereditary gout result from an enzyme defiance hypoxanthine-guanine phosphoribosyltransferase (HGPRT)).
Decreased uric acid serum (hypouricemia)	<ul style="list-style-type: none">• Liver disease (Decreased Production).• Fanconi syndrome (Increased excretion).

Urine Uric acid

- This test evaluates uric acid metabolism in gout and renal calculus formation.
- The uric acid urine test measured in a sample of urine collected over 24 hours.
- A **high level of uric acid** in the urine means that the patient is **more** likely to develop **uric acid kidney stones.**

Case	Cause
Increased urine uric acid (uricosuria)	Tubular reabsorption defect (Fanconi syndrome)
	Multiple myeloma, lymphoma
	Lesch-Nayhan syndrome
Decreased urine uric acid	Kidney disease

Gout

- Excess **monosodium urate crystallizes** and deposits (*needle like crystals*) in the joints, soft tissues, and organs.
- This will lead to **inflammation** of tissues → This inflammation is responsible for the crisis symptoms acute gouty arthritis.

Notes:

- Hyperuricemia **does not** always lead to gout. Less than 20% of cases develop into arthritic gout disease.
- Uric acid level is just one of several criteria necessary for diagnosis.
- Blood test results can be misleading, though → Some people have high uric acid levels, but never experience gout, and some people have signs and symptoms of gout, but don't have unusual levels of uric acid in their blood.




Figure 5. Chemical structure of monosodium urate




Figure 6. monosodium urate crystals

Practical Part

Objectives:

- To estimate the amount of uric acid in blood by using uric acid liquicolor kit

Principle (of the used kit):

1. Kit contains:

The enzyme reagent used includes: buffer, uricase, peroxidase, 4-Aminophenazone and 3,5-Dichloro-2-hydroxybenzene-sulfonic acid (DCHBS) .

1. Uric acid in the sample oxidized by **uricase** to allantoin and hydrogen peroxide.



2. Hydrogen peroxide reacts with 3,5-Dichloro-2-hydroxybenzene-sulfonic acid (DCHBS) and 4-aminophenazone (PAP) in the presence of **peroxidase** (Hydrogen peroxide oxidoreductase) to yield a red-violet quinoneimine dye (**chromogen**).



- The intensity of the dye is measured at 520nm and it is **directly proportional** to the concentration of uric acid present in the sample.

Materials:

Uric acid liquid (QCA) kit.

Method:

	Blank	Standard	Test
Buffer	1ml	1ml	1ml
Standard	—	0.02 ml	—
Sample	—	—	0.02 ml

1. Water bath at 37°C for 5 min.
2. Read absorbance at 520 nm.

Calculations of the Results:

$$\text{-Serum uric acid} = \frac{A (\text{Sample})}{A (\text{Standard})} \times \text{Std. Conc. (8 mg/dl)}$$

Reference value in serum:

Serum: 2.5–7.0 mg/dL

Urine: 250 –750 mg/24 hour

Homework

- Why does uric acid increase in leukemia?
- Why does Lesch-Nyhan syndrome cause hyperuricemia?
- What is the mechanism of action of allopurinol in managing gout?

References

- Bobulescu, I. A., & Moe, O. W. (2012). Renal Transport of Uric Acid: Evolving Concepts and Uncertainties. *Advances in Chronic Kidney Disease*, 19(6), 358–371. <http://doi.org/10.1053/j.ackd.2012.07.009>.
- A Manual of Laboratory and Diagnostic Tests 9th edition (January , 2014), Frances T Fischbach RN, BSN, MSN By Lippincott Williams & Wilkins Publishers.