


# DETERMINATION OF PLASMA AMYLASE



# Amylase

- Amylase is an enzyme that catalyzes the breakdown of starch and glycogen by hydrolysis of internal  **$\alpha$ -1,4-glycoside bonds** into smaller carbohydrate groups (maltose, oligosaccharides, glucose).
- It is produced in the **salivary glands, pancreas, liver, and fallopian tubes** and is **normally** excreted in **small** amounts in the urine.




Figure 1. Chemical reaction of amylase

# Amylase main sources

- Among healthy individuals, the **pancreas** and the **salivary glands** account for almost all serum amylase, 40-45% from the pancreas and 55-60% from the salivary glands.
- Electrophoresis shows that serum amylase is of **2 main types**:
  - P-type amylase from the **pancreas**
  - S-type amylase from the **salivary glands**




Figure 2. The anatomy of the pancreas




Figure 3. The locations of salivary gland.  
Source: my.clevelandclinic.org

# Amylase in Serum and Urine

- This test of **blood** and **urine** is most often **used to distinguish acute pancreatitis and other causes of abdominal pain** that require immediate surgery.
- If the pancreas or salivary glands are inflamed, much more of the enzyme enters the blood and, consequently, more amylase is excreted in the urine.
- Serum and urine amylase measurement in addition to other laboratory tests, **amylase clearance**, **amylase isozyme** , and measurement of **serum lipase levels**, increase the specificity of amylase measurement in the diagnosis of acute pancreatitis.

# Pancreas Function Test

- Blood levels of the pancreatic enzymes **amylase** and **lipase** are measured.
- This test used to diagnose and monitor treatment of **acute pancreatitis**.
- **lipase** test has become a much more **sensitive** and **specific biomarker** in diagnosing acute pancreatitis.




Figure 4. Pathology of acute pancreatitis. *Source: my.clevelandclinic.org*

## **Range of expected values of amylase:**

**Serum :** 16 -108 U/L

**Urine:** 0 - 14 U/hour

## **Increased plasma amylase (hyperamylasaemia):**

1. Salivary gland inflammation
2. Acute pancreatitis
3. Pancreatic cancer
4. Obstruction of pancreatic duct

## **Decreased plasma amylase:**

1. Pancreatic insufficiency
2. Liver disease
3. Kidney disease
4. Cystic fibrosis
5. Pregnancy

# Practical Part

---

## Objective:

To estimate the concentration of amylase in serum.

## Principle (of the used kit):

1- **Amylase** hydrolyzed p-nitrophenyl D-maltoheptoside (**PNPG7**) to P-nitrophenylmaltotriose (**PNPG3**) and **maltotetrose**:



2- **Glucoamylase** hydrolyzes **PNPG3** to P-nitrophenylglycosie (**PNPG1**) and **glucose**:



3-Then **PNPG1** is hydrolyzed by **glycosidase** to **glucose** and **P-nitrophenol** which produce a **yellow color** which absorb at 405nm, the rate of **increase** in Ab is measured at 405 nm and is proportional to the amylase activity in the sample:





## Materials:

Amylase (color/kinetic) kit (UDI).

2) Applications → 2) Simple Kinetics → wave-length (405 nm) → 1) Seconds  
→ Duration (120 sec = 2 min) → Intervals (30 sec) → Print Data Table (off)  
→ Press start (2 times)

## Method:

Chemicals	Sample
Amylase substrate	1.0 ml
Pre-warm at 37°C for 5 minutes and add:	
Sample1	0.025 ml

1. Mix and incubate at 37°C for **90 seconds** and read the absorbance at **405 nm** against distilled water.
2. Continue readings every **30 seconds for 2 minutes** and determine  $\Delta A/\text{min}$ .

# Results

<b>Time (Seconds)</b>	<b>Absorbance at 405 nm</b>
<b>0</b>	
<b>30</b>	
<b>60</b>	
<b>90</b>	
<b>120</b>	

# Calculations

-Amylase Activity in TEST (U/L)=  $\Delta A/\text{min}$  x 4824

$$\Delta A/\text{Min} = (\Delta A1 + \Delta A2) \div 2$$

$$\rightarrow \Delta A1 = (A60s - A30s) + (A30s - A0s)$$

$$\rightarrow \Delta A2 = (A120s - A90s) + (A90s - A60s)$$

## References

- Fischbach FT, Dunning MB. A Manual of Laboratory and Diagnostic Tests. Lippincott Williams & Wilkins, 2009 .p. 419-420.
- Ismail OZ, Bhayana V. Lipase or amylase for the diagnosis of acute pancreatitis? Clin Biochem. 2017 Dec;50(18):1275-1280. doi: 10.1016/j.clinbiochem.2017.07.003. Epub 2017 Jul 16. PMID: 28720341.
- Amylase Test: MedlinePlus Medical Test. (2025, Mar 3). Retrieved from <https://medlineplus.gov/lab-tests/amylase-test/>
- BCH472 practical note.