

ESTIMATION OF SERUM UREA



Urea

- Urea is the highest **non-protein nitrogen** compound in the blood.
- Urea is the major excretory product of protein metabolism.
- It is formed by **urea cycle** in the liver from **free ammonia** generated during protein catabolism.
- Since historic assays for urea were based on measurement of nitrogen, the term **blood urea nitrogen (BUN)** has been used to refer to urea determination.




Figure 1. Chemical structure of urea

Urea synthesis

- **Protein metabolism** produces amino acids that can be oxidized.
- This results in the release of **ammonia** which is **converted to urea** (via *urea cycle*) and excreted as a waste product.
- Following synthesis in the **liver**, urea is carried out in the **blood to the kidney** which is readily **filtered from the plasma by glomerulus**.
- Most of the urea in the glomerular filtrate is excreted in the urine, and **some** urea is reabsorbed through the renal tubules.
- The **amount reabsorbed** depends on **urine flow rate** and **extent of hydration** (the amount of urea reabsorbed increases with dehydration).

Dehydration → ↓ Urine Flow → Concentrated Urine ↑ Urea reabsorption

The concentration of urea in the plasma is determined by:

1. **Renal and liver function**
2. **The protein content in diet**
3. **The rate of protein catabolism**

Urea synthesis




Figure 2. Urea cycle

Clinical Application

Measurement of urea used to in:

- Evaluate renal function
- To assess hydration status
- To determine nitrogen balance
- To aid in the diagnosis of renal diseases
- To verify adequacy of dialysis
- Check a person's protein balance

Plasma urea Concentration

1. Measurement of Blood Urea Nitrogen (BUN) alone is less useful in diagnosing kidney diseases because it's blood level is influenced by **dietary protein and hepatic function** (*why?*).
2. But its diagnostic value improves with **serum creatinine values**.

	Type	Cause	Note
High urea (High urea concentration in plasma is called uremia) <i>(Uremia vs Azotemia)</i>	Pre-renal	<ul style="list-style-type: none"> ▪ Cognitive heart failure ▪ <u>Dehydration</u> ▪ High protein diet ▪ Increased protein catabolism 	<ul style="list-style-type: none"> ▪ Cognitive heart failure → reduced renal blood flow, less blood is delivered to kidney , then less urea is filtered
	Renal	<ul style="list-style-type: none"> ▪ Renal failure 	
	Post-renal	<ul style="list-style-type: none"> ▪ Urinary tract obstruction 	
Low urea		<ul style="list-style-type: none"> ▪ Low protein intake ▪ Liver disease ▪ Pregnancy 	

Urine urea Concentration

- The Urine Urea Nitrogen test (UUN) determines how much urea is in the urine to assess the amount of protein breakdown.
- The test can help determine **how well the kidneys are functioning**, and if the intake of **protein is too high or low**.
- **Specimen:** The urine urea nitrogen test is performed by collecting a **24-hour urine sample**.

	Cause
High urea in urine	<ul style="list-style-type: none">▪ Too much protein in the diet.▪ Too much protein breakdown in the body.
Low urea in urine	<ul style="list-style-type: none">▪ Malnutrition.▪ Too little protein in the diet.▪ Kidney issues.

Reference Value

SPECIMEN	UREA NITROGEN	UREA
Serum/Plasma	5-23 mg/dL	10-50 mg/dL
Urine 24 h	9-16g/24h	20-35 g/24 h

💡 Urea nitrogen VS urea ?




Figure 3. Chemical structure of urea

Practical Part

Objectives:

- Estimation of blood urea nitrogen (BUN).

Principle (of the used kit):

The Reagent used contains: **Urease, Glutamate Dehydrogenase, NADH, 2-oxoglutarate, buffers and stabilizers .**

1. **Reaction one:** Urea is hydrolysed in the presence of urease enzyme and water to yield ammonia and carbon dioxide.



2. **Second reaction:** The ammonia reacts with 2-oxoglutarate and reduced nicotinamide adenine dinucleotide (NADH) in the presence of glutamate dehydrogenase (GLDH) to yield glutamate and nicotinamide adenine dinucleotide (NAD).



The amount of the **urea** in the sample is proportionally related to the reduced absorbance at 340 nm as a result of **NADH oxidation** to NAD.

Materials:

Stanbio urea nitrogen (BUN) liquid-UV procedure

Method:

	Standard	Serum
Working reagent	1 ml	1 ml
Pre-warm at 37°C for 3 min. and add:		
Standard	0.01/10 μ l	-
Serum	-	0.01/10 μ l

1. After exactly **30 seconds**, read and record absorbance A_1 against distilled water at **340 nm**.
2. At exactly **60 seconds** after A_1 , read and record the absorbance A_2 and determine ΔA ($A_1 - A_2$).

Calculations of the Results:

$$\text{-Serum BUN (mg/dL)} = \frac{\Delta A (\text{Sample})}{\Delta A (\text{Standard})} \times \text{Std. Conc. (30mg/dL)}$$

$$\text{-Serum urea (mg/dL)} = \text{BUN} \times 2.14$$

References

- Clinical Chemistry: Techniques, Principles, Correlations (Bishop, Clinical Chemistry) Mar 31, 2009, by Michael L. Bishop MS MT (ASCP) CLS (NCA) and Edward P. Fody MD
- Urea nitrogen urine test: MedlinePlus Medical Encyclopedia. (2025, Feb 13). Retrieved from <https://medlineplus.gov/ency/article/003605.htm>