

CLS 281

Basic Biochemistry and Biomolecules


جامعة  
الملك سعود  
King Saud University




Experiment 4

# General Color Tests for Carbohydrates

# Carbohydrate Classification



# Aldose vs Ketose




**Aldose** has a **carbonyl group** at the end of the carbon chain.  
**Ketose** has a **carbonyl group** in the middle of the carbon chain.

Two monosaccharides are joined together by a

# Glycosidic Bond

formed by the loss of water molecule.




# General Color Tests for Carbohydrates

1. Molisch Test
2. Anthrone Test

# 01 Molisch Test

- It is a general test used to detect the presence of carbohydrates in a given sample.
- **Molisch's Reagent**
  - concentrated sulphuric acid H<sub>2</sub>SO<sub>4</sub>
  - Phenol-type molecules (such as α-naphthol or thymol).
    - α-Naphthol.
    - Thymol:
      - It can be used as a reagent instead of α-Naphthol.
      - Thymol is more stable than α-Naphthol and can be applied to insoluble carbohydrates like cellulose or wood.



α-naphthol


# 01 Molisch Test Principle

- **Principle**


Monosaccharides + H<sub>2</sub>SO<sub>4</sub>  $\xrightarrow{\text{Heat}}$  furfural (of furfural derivatives) + 3 H<sub>2</sub>O

furfural (of furfural derivatives) +  $\alpha$ -Naphthol (2 phenol groups)  $\rightarrow$  **purple ring** at the interface.

- Carbohydrate undergoes **dehydration** upon the introduction of concentrated sulphuric acid, resulting in the formation of an **aldehyde**.
- This aldehyde undergoes **condensation** along with two  $\alpha$ -naphthol, resulting in the formation of a **purple or reddish-purple-colored complex**.
- The **purple ring color** is due to condensation products of furfural or its derivatives with  $\alpha$ -Naphthol.



# Furfural



- Furfural is an organic compound derived from a variety of agricultural byproducts, including corncobs, oat, and wheat bran.
- The name furfural comes from the Latin word furfur, meaning bran, referring to its usual source.
- It is not a carbohydrate.



## 01

# Molisch Test

- It is effective for any compound which can be dehydrated to furfural or substituted furfural (such as hydroxymethyl furfural) by concentrated sulfuric acid.
- If the carbohydrate is an oligosaccharide (e.g., disaccharide, trisaccharide ...etc.) or a polysaccharide, the hydrolysis of the carbohydrate acetal linkage occurs simultaneously with the dehydration reaction (in polysaccharide, the color develops slower).

## 01

# Molisch Test Procedure

Note: swirl the samples and reagent bottles before use to aspirate homogenous solution.

Steps	Tube No.	Tube 1	Tube 2	Tube 3
1	Sample	2 ml of 0.5% starch	2 ml of 0.5% sucrose	2 ml of 0.5% glucose
2	$\alpha$ -Naphthol Reagent	2 drops	2 drops	2 drops
3	Mix			
4	Con. H <sub>2</sub> SO <sub>4</sub>  Reagent	2 ml	2 ml	2 ml




- Incline the test tube slowly and carefully add 2 ml of concentrated sulfuric acid (H<sub>2</sub>SO<sub>4</sub>) down the side of the tube to form a layer below the sugar solution.

Positive result >>> purple ring at the interface indicate the presence of carbohydrate.

## 01


# Molisch Test Result

- All carbohydrates are detected using this test.
- Monosaccharides give a rapid positive test.
- Disaccharides and polysaccharides react slower.
  
- A negative result by this reaction is very good evidence of the absence of carbohydrates, but a positive test is an indication of the probable presence of carbohydrates.



## 02 Anthrone Test


- It is another general test for carbohydrates.
- **Aim:** To detect the presence of carbohydrates in a given solution.
- **Reagent**
  - Concentrated sulphuric acid  $H_2SO_4$
  - Phenol-type molecules (anthrone)
- Based on similar reactions in which anthrone is used instead of  $\alpha$ -Naphthol in the Molisch test to form a green to blue-green color product.



Anthrone

## 02

# Anthrone Test Principle



1. Hydrolysis of sugar to monosaccharides.
2. Dehydration of monosaccharides by H<sub>2</sub>SO<sub>4</sub> producing furfural (of furfural derivatives).
3. The reaction of furfural with Anthrone reagent resulted in a **Blue-green complex**.

## 02

# Anthrone Test Procedure

Steps	Tube No.	Tube 1	Tube 2	Tube 3	Tube 4	Tube 5
1	Sample	Small piece of Filter paper	1 drop 0.5%starch	1 drop of 0.5%sucrose	1 drop 0.5%glucose	1 drop Blank
2	Water As a Diluent	1ml	1ml	1ml	1ml	1ml
3	Anthrone reagent	3ml	3ml	3ml	3ml	3ml

4. Heat in boiling water bath for 3 mins.

# 02

## Anthrone Test

- It is very sensitive.
- It will give a positive reaction with filter paper (cellulose).

### Uses:

- It can be used for the quantitative determination of glycogen and sugar in the blood.
- It can be used as a qualitative test since different sugars dehydrate at different rates and produce a variety of colors.

# Summary of General Color Tests for Carbohydrates

Test	Detect	Reagent	Principle	Positive Result	Negative Result	Note
<b>Molisch Test</b>	Detect the presence of <b>carbohydrates</b> in a given sample.	<ul style="list-style-type: none"> <li>- concentrated <u>sulphuric acid</u> H<sub>2</sub>SO<sub>4</sub>.</li> <li>- Phenol-type molecules (<b>α-Naphthol</b>).</li> </ul>	<p>1- Hydrolysis of sugar to monosaccharides.</p> <p>2-Monosaccharides + H<sub>2</sub>SO<sub>4</sub> → furfural (of furfural derivatives) + 3 H<sub>2</sub>O (dehydration reaction)</p> <p>3- furfural (of furfural derivatives) + α-Naphthol (2 phenol groups) → <b>purple ring</b> at the interface. (condensation reaction)</p>	<b>purple ring color</b>	Colorless	<p><u>Monosaccharides</u> give a <b>rapid</b> positive test.</p> <p><u>Disaccharides and polysaccharides</u> react <b>slower</b>.</p>
<b>Anthrone Test</b>	Detect the presence of <b>carbohydrates</b> in a given sample.	<ul style="list-style-type: none"> <li>- concentrated <u>sulphuric acid</u> H<sub>2</sub>SO<sub>4</sub>.</li> <li>- Phenol-type molecules (<b>anthrone</b>)</li> </ul>	<p>1- Hydrolysis of sugar to monosaccharides.</p> <p>2-Monosaccharides + H<sub>2</sub>SO<sub>4</sub> → furfural (of furfural derivatives) + 3 H<sub>2</sub>O (dehydration reaction)</p> <p>3-The reaction of furfural with Anthrone reagent resulted in a <b>Blue-green complex</b>.</p>	<b>Blue-green complex.</b>	Reagent color	<p>It is very sensitive.</p> <p>It will give a positive reaction with filter paper (cellulose).</p>



# Guideline for writing the lab report

## Total: 5 marks

All the following information should be included in your report:

- a) Course # (CLS 281)
- b) Experiment title
- c) Date of the experiment
- d) Student's names and university ID#
- e) Section #

The lab report is broken down into 6 sections:

1. Experiment **title**
2. The **aim** of the experiment (objective, or what the test detects specifically) (1 mark)
3. **Principle** (chemical reaction) (1 mark)
4. **Methodology** (written in **steps**, **not in tables**)
5. **Result** (1 mark)
6. **Interpretation or Comment** (2 mark)

Deadline: Next lab    Submission: via email