


## Designing PCR Primers

**A primer** is a short synthetic oligonucleotide which is designed to have a reverse sequence that is complement to a region of interest in a target DNA



**Some thoughts on designing primers:**

1. Primers should be 17-28 bases in length
2. Base composition should be 50.60% (G+C)
3. Primers should end (3') in a G or C, or CG or GC: this prevents "breathing" of ends and increases efficiency of priming
4. Tms between 55.80°C are preferred

## Primer design steps:

The screenshot shows the Ensembl genome browser search interface. At the top, the Ensembl logo is followed by navigation links: BLAST/BLAT, BioMart, Tools, Downloads, Help & Documentation, Blog, and Mirrors. Below this is a search bar with a dropdown menu set to 'Human' and a search input field containing 'SDHA'. A 'Go' button is to the right of the search field. Below the search bar, there is a suggestion: 'e.g. BRCA2 or rat 5:62797383-63627669 or coronary heart disease'. The main content area is divided into several sections: 'Browse a Genome' with a description of the Ensembl project; 'Popular genomes' with icons and names for Human (GRCh38), Mouse (GRCm38), and Zebrafish (Zv9); and a grid of tool tiles including 'ENCODE data in Ensembl', 'Variant Effect Predictor', 'Gene expression in different tissues', 'Find SNPs and other variants for my gene', 'Retrieve gene sequence', and 'Compare genes across species'. A 'Log in to customize this list' link is also present.

**Figure 1.** Go to ensemble data base and select (human) and type the gene name *SDHA* in this example.

The screenshot shows the Ensembl genome browser interface for the SDHA gene. The top navigation bar includes the Ensembl logo, links for BLAST/BLAT, BioMart, Tools, Downloads, Help & Documentation, Blog, and Mirrors, along with a search bar and a 'Login/Register' link. The main header indicates the species is Human (GRCh38), the location is 5,218,241-256,700, and the gene is SDHA.

The left-hand menu, titled 'Gene-based displays', lists various options. The 'Sequence' option is highlighted with a red box. Below it, there are buttons for 'Configure this page' (also highlighted with a red box), 'Add your data', and 'Export data'.


The main content area displays the following information for the SDHA gene (ENSG00000073578):

- Description:** succinate dehydrogenase complex, subunit A, flavoprotein (Fp) [Source:HGNC Symbol;Acc:HGNC:10680]
- Synonyms:** FP, SDH2, SDHF
- Location:** [Chromosome 5: 218,241,256,700](#) forward strand.
- INSDC coordinates:** chromosome:GRCh38:CM000667.2:218241:256700:1
- Transcripts:** This gene has 20 transcripts (splice variants) [Show transcript table](#)

The 'Marked-up sequence' section includes a 'Download sequence' button and a 'BLAST this sequence' button. Below this is a 'Key' section with a 'Features' box containing the text 'All exons in this region'. The sequence itself is displayed with exons highlighted in yellow.

```
>chromosome:GRCh38:5:217641:257300:1
GTGGGAAAGCTGCAAAACAGTGTAGAGCTTCTGAGGCTCCTCCTGTTGCTGATGATG
CAGAGATTACAGATTTCTCTGCTCGGGCTCAGCATGCACCCCTCCATAGCTAGGCTGGGG
GCTAACGTTAGTAAGAAATTTACTGATTTACTGATACCTACAGAAATTTTCACATTTTCA
CAAAGTAGAGTCACTTTAAGTTTGTATACATGTGACCTCCACAAACACACAAATFACAA
AATAGTTTACTTAATCAGGGATGTTAATAAACAAGAGTTCGACAGGCCCCGGGCGGTC
CCGCCACCTCCGGGACGAGGCCCCCTCCGACCCCAHTCCCTCGGGTCTGGGGCTTC
AGGCCCTTCGGTCTGGGGGATCCGGAGAACACCCACGGGGCTTTAAAATGTTGGTGC
CCACCCTCCCGGAAACAGGGCCCGCTCTACCTCGGTCCGGGAGGCGGGACCTCAGGG
TTCCCTTAACGCCACCGTCCGGGGTCCGGTTTCCGACGGGCGGGGCCCTCAGTA
CCCGCTCCGGGCTGGCATGGTGGCAGGGGGATGTCGCCACTGCAGCCCCGCTCGAC
TCCGGGTGGTGGCAGGGGGGTATCCCCCTCCCGCCAGCTGCACCCGGTGTGGT
GGCAGGGCCAGCTCCGCCAGGACTGGCGGACTGGCGGGCCACACAGACATGTC
```

**Figure 2.** Selects (sequence) from the list on the left then select (configure this page) to show variations.



**Figure 3.** The entire sequence (Exons and introns) of the whole gene will appear, copy the gene sequence and paste in new word document to make the designing process much easier.

```

>chromosome:GRCh38:5:217641:257300:1
217641 GTGGGGAAGGTGCAAACAGTGCTAGAAAGCTTCTGAGGCCCTCCTCCTGGTGGTGATGAATG 217700
217701 CAGAGATTACRAGAYTCTCTGCTTCGGGCTCAGCATGCACCCTCCATAGCTAGGCTGGGG 217760
217761 GCTAACGTTAGTAAGAATTTACTGATTTACTGAGTACCTACACGARGTTTCACATTTTCA 217820
217821 CAAAGTAGAGTCACTTTAAGTTTGTTTACATGTGAMCTCCACAAACACACACRATTACAA 217880
217881 AATAGTTTAGTTAATCAAGGGATGTATAATAAACAAGAGGYCGGAGGCCCCGGGCGCTC 217940
217941 CCGCCACCTCCGCGGACGAGCGCCCTCCGACCCATTCCCTGCGGTCTGGGCGTTC 218000
218001 AGGCCCTTCGGTCTGGGCKATCCCGGAGAACCACCACGGGGCTTTAAAAATGTTGGTGC 218060
218061 CCACCACCTCCCGGAACAGGGCCCGCTCTACCTCGGTCGGGGAGCGCGGRACCTCAGCG 218120
218121 TTCCTTAACGCCACCCTCCGCGGGTCCGCTTTGCGCAGGCGCGGCCCCCACTTCAGTA 218180 EX1
218181 CCCGCTCCGGGCGTGGCATGGTGCAGGCGCGATGTCCCCACTGCAGCCCCGCTCGAC 218240
218241 YCCGGCGTGGTGCAGGCGCGGTATCCCCCTCCCGCCAGCTCGACCCCGGTGTGGT 218300
218301 GCGCAGGCGCAGTCTGCGCAGGGACTGGCGGGACTGCGCGGCGGCRACAGCRGACVTGTC 218360
218361 GGGGGTCCGGGRCTGTGCGGCTGCYRAGCGCTCGGCGCTGGCGCTGGCCAAGGCGT 218420
218421 GAGTCCGTGCCGCGGACCGGGGCGGGGCAAGCGGGGCCGAGGCGGCGGTAGGAGCGGGA 218480
218481 CGRTCCCCAGCGGGTCCGAGCGGAGCGGGCGCCGGGTCCCTGCGCCCTCTGTCCCGGAT 218540
218541 CGRGAAGGGGCTGAGAGAGCCCTGGGCCGTGCGAGGGGAAGCCCGGGGCGGACTCGGG 218600
218601 GACCCGGGAGCTCGGTCCCTTAGTAGATAGTCCGCGTCCGGGTGAAGGTCACAACCCCGC 218660
218661 GGGCTTGCTGGGCGTCCCCTCCGCGCYTTGGTCTYGGGCCTGGGGTCTGGGACCCCGCG 218720
218721 GGCTGAGGTAGCCCCTCGCCTCAGTGCTTGGCAGGTGGACTCGGGGAGGAGTCGTGTCTG 218780
218781 CCCAAGGTCACCCGGGCGGATAGCGGCGGTGGCCGCGCTGGCTGGGCTGGGCCTCTGCC 218840
218841 GCCCTCTGTGCGGGTTGTCCTGAGGAGCAGCCCGCAGCCCGTGGGTGGGGCCGGCGGGG 218900

```

**Figure 4.** Select the sequence of the exon and a wide range of the sequence both before and after the exon or select the entire sequence (**Exon 1** will be presented as an example) ,place square brackets [ ] both before and after the beginning and end of the first fragment by 60 -100bp (whatever is between the square brackets is untouchable). Select the fragment and about 300-400bp before and after the brackets (the more area you select before and after the fragment, the more options you provide for the primer selection). Sometimes we have to deal with a relatively big exon, In this case the exon will be divided in to two fragments (1a &1b).

```

>chromosome:GRCh38:5:217641:257300:1
217641 GTGGGGAAGGTGCAAAACAGTGCTAGAAGCTTCTGAGGCCCTCCTCRTGGTGGTGATGAATG 217700
217701 CAGAGATTACRGAGYTCTCTGCTTCGGGCTCAGCATGCACCTCCATAGCTAGGCTGGG 217760
217761 GCTAACGTTAGTAAGAATTTACTGATTTACTGAGTACCTACACGARRGTTTCACATTTCA 217820
217821 CAAAGTAGAGTCACTTTAAGTTTGTTCATGTGAMCTCCACAAACACACACRATTACAA 217880
217881 AATAGTTTAGTTAATCAAGGGATGTATAATAAACAAGAGGYCGGAGGCCCGGGCCGCTC 217940
217941 CGGCCACCTCCGCGGACGAGGCCCCMCTCCGACCCATTCCCTGCGGTCTGGGGCTTC 218000
218001 AGGCCCTTCGGTCTGGGKATCCCGGAGAACCACCACGGGGCTTTAAAAATGTTGGTGC 218060
218061 CCACCACCTCCCGGAACAGGGCCCGCTCTACCTCGGTGGGGAGCGCGRACCTCAGCG 218120
218121 TTCCCTTAACGCCACCGTCCGCGGGTCCGCTTTGCGCAGGCGCGATGTCCCCCACTGCAGCCCGCTCGAC 218180 EX1
218181 CCGCTCCGGCGTGGCATGGTGCAGGCGCGATGTCCCCCACTGCAGCCCGCTCGAC 218240
218241 YCGGGCTGGTGCGCAGGCGCGGTATCCCCCCTCCCCGCCAGCTCGACCCCGGTGTGGT 218300
218301 GCGCAGGCGCAGTCTGCGCAGGGACTGGCGGGACTGCGCGCGGRACAGRGACYVTGTC 218360
218361 GGGGGTCCGGRCCTGTGCGGGCTGCYRAGCGCTCGGCGCTGGCGCTGGCCAAGGCGGT 218420
218421 GAGTCCGTGCCCGGACCGGGGCGGGGCAGGCGGGGCCGAGGCGGCGGTAGGAGCGGGA 218480
218481 CGRTCCCCAGCGGGTCCGAGCGGAGCGGGCGCGGGTCCCTGCGCCCTCTGTCCCGGGAT 218540
218541 CGKGAAGGGGCTGAGAGAGCCCTGGGCGGTGCGAGGGGAAGCCCGGGGCGGACTCGGG 218600
218601 GACCCGGGGAGCTCGGTCCTTAGTAGATAGTCCGCGTCCGGGTGAAGGTCACAACCCCGC 218660
218661 GGGCTTGCTGGGCGTCCCTCCGCGCTTTGGTCGGGCCTGGGGTCTGGGACCCCGCG 218720
218721 GGCTGAGGTAGCCCTCGCCTCAGTGCCTGGCAGGTGGACTCGGGGAGGAGTCTGTGTG 218780
218781 CCAAGGTACCCGGGCGGATAGCGGCCGTGGCGCCCTGGCTGGGCTGGGCCTCTGCC 218840

```

Primer3web

ersion 4.0.0 - Pick primers from a DNA sequence.

[disclaimer](#)   [code](#)  
[cautions](#)

Select the [Task](#) for primer selection generic

Paste source sequence below (5'->3', string of ACGTNacgtn -- other letters treated as N -- numbers and blanks ignored). FASTA format ok. Please N-out undesirable sequence (vector, ALUs, LINEs, etc.) or use a [Mispriming Library \(repeat library\)](#) NONE

```

218241 YCCGGCGTGGTGCAGGCGCGGTATCCCCCTCCCCGCCAGCTCGACCCCGGTGGT 218300
218301 GCGCAGGCGCAGTCTGCGCAGGGACTGGCGGGACTGCGCGCGGCRACAGCRGACTGTC 218360
218361 GGGGGTCCGGRCCTGTGCGGGCTGCRYAGCGCTCGGCGCTGGCGCTGGCCAGGGCGGT 218420
218421 GAGTCCGTGCCCGGACCGGGGCGGGGCAAGCGGGGCCGAGGCGGCGGTAGGAGCGGGA 218480
218481 CGRTCCCCAGCGGGTCCGAGCGGAGCGGGGCGGGTCCCTGGGCGCTCTGTCCCGGGAT 218540
218541 CGKGAAGGGGCTGAGAGAGCCCTGGGCGGTGCGAGGGGAAGCCCGGGGCGGACTCGGG 218600

```

<input checked="" type="checkbox"/> Pick left primer, or use left primer below	<input type="checkbox"/> Pick hybridization probe (internal oligo), or use oligo below	<input checked="" type="checkbox"/> Pick right primer, or use right primer below (5' to 3' on opposite strand)

Pick Primers
Download Settings
Reset Form

**Figure 5.** Go to primer 3 data base, and paste in the blank empty window and select pick primers. (At this point, in computer language, you are ordering primer 3 to select primers for the selected fragment anywhere out the square brackets).



```

>chromosome:GRCh38:5:217641:257300:1
217641 GTGGGGAAGGTGCAAACAGTGCTAGAAAGCTTCTGAGGCCCTCCTC RTGGTGGTGATGAATG 217700
217701 CAGAGATTACR GAGY TCTCTGCTTCGGGCTCAGCATGCACCCTCCATAGCTAGGCTGGGG 217760
217761 GCTAACGTTAGTAAGAATTTACTGATTTACTGAGTACCTACACGAR GTTTCACATTTTCA 217820
217821 CAAAGTAGAGTCACTTTAAGTTTGTTTACATGTGAM CTCCACAAACACACACR ATTACAA 217880
217881 AATAGTTTAGTTAATCAAGGGATGTATAATAAAACAAGAGGY CGGAGGCCCGGGCCGCTC 217940
217941 CCGCCACCTCCGCGGACGAGCGCCCTCTCCGACCCCATTCCTGCGGTCTGGGCGTTC 218000
218001 AGGCCCTTCGGTCTGGGCK ATCCCGGAGAACCACCCACGGGGCTTTAAAAATGTTGGTGC 218060
218061 CCACCACCTCCCCG GAACAGGGCCCGCTCTAC CTGGTCGGGGAGCGCGGR ACCTCAGCG 218120
218121 TTCCCTTAACGCCACCGTCCGCGGGTCCGCTTTGCGCAGGCGCGGCCCCCACT [ CAGTA 218180 EX1
218181 CCCGCTCCGGGCGTGGCATGGTGCAGGCGCGATGTCCCCACTGCAGCCCCGCTCGAC 218240
218241 YCCGGCGTGGT GCGCAGGCGCGGTAT CCCCC CTCCCCGCCAGCTCGACCCCGGTGTGGT 218300
218301 GCGCAGGCGCAGTCTGCGCAGGACTGGCGGGACTGCGCGGCGGCR ACAGCRGAC YTGTC 218360
218361 GGGGGTCCGGGR CCTGTGCGGGCTGC YRAGCGCTCGGCGCCTGGCGCTGGCCAAGGCGGT 218420
218421 GAGTCCGTGCCGCGGACCGGGCGGGGCA GG CGGGGGCCGAGGCGGCGGTAGGAGCGGGA 218480
218481 CGR TCCCCAGCGGGTCCGAGCGGAGCGGGCGCCGGGTCCCTGCGCCCTCTGTCCC] GGGAT 218540
218541 CGKGAAGGGCTGAGAGAGCCCTGGGCCGTGCGAGGGGAAGCCCGGGGCGGACTCGGG 218600
218601 GACCCG EGGAGCTCGGTCCTTAGTAG ATAGTCCGCGTCCGGGTGAAGGTCAACAACCCCGC 218660
218661 GGGCTTGCTGGGCGTCCCCTCCGCCGY TTGGTC YGGGCTGGGGTCTGGGA CCCC GCG 218720
218721 GGCTGAGGTAGCCCCTCGCCTCAGTGCTGGCAGGTGGACTCGGGGAGGAGTCGTGTCTG 218780
218781 CCCAAGGTCACCCGGGCGGATAGCGGCCGGTGGCCGCCCTGGCTGGGCTGGGCCTCTGCC 218840
218841 GCCCTCTGTGCGGGTTGTCCTGAGGAGCAGCCCGCAGCCCGTGGGTGGGGCCGGCGGGGC 218900
218901 GGGTGAR ACCGCCCGGGTGGGTGCGAGGAGTGGCCGGGCTCGGCCGGTGGGCGTCC RGT 218960
218961 GGAAGCGTGGCGCGCCGAGCTTAGGCTGTCAGTTGCGCTTCCAGAAAGCGCAAATCT 219020

```

**Figure 7.** Finally, highlight both forward and reverse primers to make sure that they cover the whole gene.

## H.W

Order	Mapped genes	OMIM ID	OMIM Title	Gene Map Disorder	Expression in the nervous system and muscles
1	SDHA( <i>5p15</i> )	#252011	Succinate Dehydrogenase Complex, Subunit A	-Mitochondrial respiratory chain complex II -Leigh syndrome -Cardiomyopathy, dilated -Parangliomas 5	↑
2	PDCD6( <i>5p15.33</i> )	*601057	Programmed Cell Death 6	-Neurodegenerative disorders -Apoptosis -Cancer -Autoimmune disease -Congenital malformations -Immunodeficiency	↑
3	AHRR( <i>5p15.3</i> )	*606517	Aryl-Hydrocarbon Receptor Repressor	-Ahh inducibility -Cancer lung -Liver diseases -Tumors -Cancer	↑
4	EXOC3( <i>5p15.33</i> )	*608186	Exocyst Complex Component	NA	↑
5	SLC9A3( <i>5p15.3</i> )	*182307	Solute Carrier Family 9, Subfamily A (NHE3, Cation Proton Antiporter 3), Member 3	-Metabolic acidosis -Acidosis -Hypokalemia -Adenoma -Inflammatory bowel diseases -Diabetic nephropathies -Renal failure -Proteinuria -Colon cancer -Cystic fibrosis	↑
6	TPPP( <i>5p15.3</i> )	*608773	Tubulin Polymerization Promoting Protein	-Lewy body dementia -Multiple system atrophy	↑
7	TRIP13( <i>5p15.33</i> )	*604507	Thyroid Hormone Receptor Interactor 13	-Lung cancer	↑
8	NKD2( <i>5p15.3</i> )	*607852	Naked Cuticle Homolog 2 ( <i>Drosophila</i> )	-Hyperopia	↑
9	SLC12A7( <i>5p15</i> )	*604879	Solute Carrier Family 12 (Potassium/Chloride Transporter), Member 7	-Renal tubular acidosis	↑
10	SLC6A19( <i>5p15.33</i> )	*608893	Solute Carrier Family 6 (Neutral Amino Acid Transporter), Member 19	-Hartnup disorder (HND) -Hyperglycinuria (HG) -Iminoglycinuria (IG) , digenic -Aminoaciduria	↑
11	SLC6A18( <i>5p15.33</i> )	*610300	Solute Carrier Family 6 (Neutral Amino Acid Transporter), Member 18	-Iminoglycinuria -Hyperglycinuria -Iminoglycinuria, digenic	↑

Order	Mapped genes	OMIM ID	OMIM Title	Gene Map Disorder	Expression in the nervous system and muscles
12	TERT( <i>5p15.33</i> )	+187270	Telomerase Reverse Transcriptase	-Aplastic anemia (AA) -Dyskeratosis congenita, autosomal dominant, 2 (DKCA2) -Pulmonary fibrosis, and/or bone marrow failure, telomere-related, 1 (PFBMFT1) -Dyskeratosis congenita, autosomal recessive, 4 (DKCB4) -Pulmonary fibrosis, idiopathic (IPF) -Melanoma, cutaneous malignant 9 (CMM9) -Laryngeal squamous cell carcinoma -Benign meningioma -Cancer, Tumors	↑
13	CLPTM1L( <i>5p15.33</i> )	*612585	CLPTM1-Like Protein	-Cleft lip -Lung cancer	↑
14	SLC6A3( <i>5p15.3</i> )	*126455	Solute Carrier Family 6 (Neurotransmitter Transporter), Member 3	Parkinsonism-dystonia infantile (PKDYS)	↑
15	LPCAT1( <i>5p15.33</i> )	*610472	Lysophosphatidylcholine Acyltransferase 1	NA	↑
16	MRPL36( <i>5p15.3</i> )	*611842	Mitochondrial Ribosomal Protein L36	NA	↑
17	NDUFS6( <i>5p15.33</i> )	*603848	NADH Dehydrogenase (Ubiquinone) Fe-S Protein 6, 13kDa (NADH-Coenzyme Q Reductase)	-Mitochondrial complex I deficiency -Leigh disease	↑
18	IRX4( <i>5p15.3</i> )	*606199	Iroquois Homeobox 4	NA	↑
19	IRX2( <i>5p15.33</i> )	*606198	Iroquois Homeobox 2	-Soft tissue sarcoma -Squamous cell carcinoma of the head and neck	↑
20	C5orf38( <i>5p15.33</i> )	*610522	Chromosome 5 Open Reading Frame 38 (Coordinated Expression To IRXA2 Protein)	-Hypertension	↑
21	IRX1( <i>5p15.3</i> )	*606197	Iroquois Homeobox 1	NA	↑