

Lab sheet #5**-Qualitative Analysis of Renal Calculi-****Objectives:**

- Identification and qualitative analysis of kidney stones.

Method:**(1) Test for Uric acid:**

1. Put a small amount of the sample.
2. Add 5-7 drops of concentrated nitric acid (*Carefully*).
3. Heating in a water bath.

➔ (The positive result is yellow to orange color on the inner surface of the test tube).

(2) Test for carbonate:

1. Add 0.5 ml con. hydrochloric acid (2M HCL) to a small portion of the sample.

➔ (Gas bubbles will indicate the presence of carbonate).

(3) Test for oxalate:

1. Heat a part of the sample with 2 ml dilutes sulphuric acid (2M H₂SO₄) for 1 min.
2. Add 2 drops (one by one) of, potassium permanganate (KMnO₄) solution and Mix.

➔ (The decolorization and evolution of bubbles will confirm the presence of oxalate).

(4) Test for phosphates:

1. Dissolve a little of the sample in about 1.5 ml of concentrated nitric acid (HNO₃).
2. Add an equal volume (1.5 ml) of ammonium molybdate solution.
3. Heat to boiling.

➔ (If phosphates are present, a yellow precipitate of ammonium phosphomolybdate is obtained).

(5) Test for calcium:

1. Dissolve small amount of the sample by heating with 2 ml dilute hydrochloric acid (2M HCL).
2. Add 1 ml ammonium oxalate.

→ (A white precipitate of calcium oxalate shows the presence of calcium).

(6) Test for magnesium:

1. On a few amounts of magnesium, add 1 ml of titan followed by 1 ml potassium hydroxide—*to be strongly alkaline*—.

→ (An orange to red color indicates the presence of magnesium).

Results:

	Observation	Type of stone/s
Uric acid		
Carbonate		
Oxalate		
Phosphates		
Calcium		
Magnesium		