## Lab (2)

# Quantitative amino acids estimation by ninhydrin method

## **Materials:**

#### Chemical

Standard amino acid stock solution (100  $\mu$ g/ml), 8% w/v of ninhydrin reagent, 50% v/v ethanol, distilled water.

#### Equipment and Glassware

Test tubes, rack, pipette, pipette pump, aluminum foil, plastic cuvette, spectrophotometer, water bath.

#### **Protocol:**

#### 1. Set up 7 test tubes as following:

Tube	Standard amino acid solution (100 µg/ml) (ml)	Distilled water (ml)	Unknown sample	Ninhydrin reagent (ml)
Blank	-	4		
Α	1.2	2.8		
В	1.6	2.4		
С	2	2		1
D	2.4	1.6		1
Ε	2.8	1.2		
Unknown sample	-	-	4	

- 1. Mix the contents of the tubes by vertexing/shaking the tubes.
- 2. Cover the tubes with aluminium foil.
- 3. Place all the test tubes in 80°C for 5 minutes.
- 4. Cool the test tubes in cold water, the add 1 ml of ethanol to each test tube and mix well.
- 5. Record the absorbance of all tubes against the blank at 570 nm using a colorimeter (spectrophotometer).
- *↓ PAUSE AND THINK → What the blank should contain? Why?*
- 6. Calculate the amino acid concentration for each standard amino acid solution using  $C1 \times V1 = C2 \times V2$  formula.
- Plot standard curve for absorbance against amino acids concentration (μg/ml) using results for solutions (A-E).
- 8. From the standard curve, estimate the concentration of the amino acids present in your unknown sample.

### **Results:**

Test tube	Amino acid concentration [µg/ml]	Absorbance at 570 nm
Blank		
Α		
В		
С		
D		
E		
Unknown		
sample		

## **Supporting materials:**

- A video explains the mechanism of reaction between ninhydrin and amino acid: <u>https://www.youtube.com/watch?v=P-iK4QUU9t0</u>
- A video shows the practical steps of the experiment step by step: -Part 1: <u>https://www.youtube.com/watch?v=wvVYF8T7uiE</u>
  -Part 2: <u>https://www.youtube.com/watch?v=B4EgYeFqV5Q</u>