




# **Characterization of The DNA by Spectrophotometric Assay and Melting Temperature (T<sub>M</sub>)**

---

# After DNA Extraction... What NEXT ?



# 1. Quantifying DNA concentration

- Is determined by measuring absorbance at **260 nm**. Why?
- At 260 nm double-stranded DNA has specific absorption coefficient of  $0.02 (\mu\text{g/ml})^{-1}\text{cm}^{-1}$ .
- So:  
→ **Concentration of DNA** =  $(A_{260} / \epsilon L) \times \text{Dilution Factor (DF)}$ .

Beer-Lambert Law:

$$A = \epsilon cl$$



## 2. DNA purity:

### 1. To detect nucleic acid purity from proteins contamination:

→ Calculate  $A_{260}/A_{280}$

- Highly purified DNA samples have a  $A_{260}/A_{280}$  nm ratio of (1.8-1.9).


→ What if the ration is below 1.8? What that means?

### 2. To detect nucleic acid purity from carbohydrates, peptides, ethanol or any organic compounds:

→ Calculate  $A_{260}/A_{230}$

- Purified DNA samples have a  $A_{260}/A_{280}$  nm ratio of (2-2.2).

# DNA and protein absorption spectrum:






## 3. GC content:

- **Denaturation:** is when the double-stranded DNA (dsDNA) unwinds {dissociated "melted"} and separates into single-stranded (ssDNA) by heat or altered pH, which breaks the hydrogen bonds between complementary bases (A = T and G≡C).
- Hyperchromic and hypochromic effect.
- The **melting temperature (T<sub>m</sub>)** is the temperature at which **50% of the DNA is unpaired** (denatured).
- GC content can be calculated by generating T<sub>m</sub> profile (DNA melting curve).


$$\%(G+C) = 2.44 (T_m - 69.3)$$

# DNA melting curve:



**FIGURE 4.4 DNA melting curve.** A melting curve of DNA showing  $T_m$  (the melting temperature) and possible molecular conformations for various degrees of melting.

# Relationship between $T_m$ and GC%:



**FIGURE 4.5** Effect of G-C content on DNA melting temperature.  $T_m$  increases with increasing percent of G + C.

What do you notice in the relation between GC content and  $T_m$ ?





# Practical Part

---



## Aim:

- Determination the concentration and purity of extracted DNA using UV spectrophotometer.
- Determination of DNA melting temperature and GC content percentage.

## Principle:

- DNA and proteins have a maximum absorbance at 260 and 280 respectively.
- dsDNA will be separated to ssDNA by heat (denaturation).
- O.D at 260 nm will increase during denaturation... Why?
- Temperature for midpoint of denaturation gives  $T_m$ . Why it is important to know  $T_m$  of DNA?
- The DNA of each species has a specific denaturation curve.. Why?



# Results:

- As in the lab sheet