

BCH 447

Estimation of Glutathione in Plasma

- Objectives:

✓ To draw the standard curve of glutathione by given known amount of glutathione assay procedure using spectrophotometric technique.

✓ To estimate the amount of glutathione in plasma sample.

Glutathione


- -It occurs naturally as a tripeptid (gamma-glutamylcysteinylglycine)
- -it has two form : Reduced form(GSH) and Oxidized form (GSSG)
- It has nucleophilic and reducing properties, those properties due to the presence of

Sulfahydryl group or thiol group (-SH).

- **<u>Nucleophelic</u>, means GSH will loss the H ion then become nucleophilic (contain negative charge) can react with the compound which carrying positive charge.
- **reducing , means it will give the other protein the H ion to maintain them in reducing form.

Functions of Glutathione

➤ It acts as anti-oxidant , blood booster , and cell detoxifier in the body.




- Physiological values of intracellular GSH generally range from 1 to 10 mM.
- Glutathione deficiencies have been linked to many forms of <u>cancer</u>.

Assay Principle

•The principle of the assay is based on the oxidation of the reduced form of glutathione (GSH) by the aromatic disulphide compound [5,5-dithiobis-2-nitrobenzoic acid (DTNB)] to form GS-TNB and the aromatic thiol, 5thio-2nitrobenzoic acid (TNB).

•Yellow color is formed. It can be measured at 412nm and is proportional to the amount of glutathione present in the sample.

Note: This method is also useful to determine the GSSG.First, GSSG will be converted to GSH by glutathione reductase(GR) and NADPH assay then it will react with DTNB reagent



- Method:

1- Prepare serial GSH concentration, as the following table:

Tube No.	GSH stock Solution (ml)	Phosphate Solution (ml)	Total Volume (ml)	GSH Conc (mg/dl)
Blank	0	3	3	0
1	0.6	2.4	3	2
2	1.2	1.8	3	4
3	1.8	1.2	3	6
4	2.4	0.6	3	8
5	3	0	3	10

2- For standard curve (Tubes from blank to 5)

<u>In a separate test tube</u>, take 0.5 ml of solution+ 2ml of Phosphate Solution + 0.25ml DTNB+ 0.25ml H₂O.

3- For Sample: Take 0.3 ml of sample +2ml of Phosphate Solution + 0.25ml DTNB+ 0.45ml H₂O.

**Incubate all tubes for 10min at 37°C→ Read absorbance at 412 nm **

- Results:

Tube No.	GSH concentration (μg/ml)	Absorbance At 412 nm
1	20	
2	40	
3	60	
4	80	
5	100	
Blank	0	
Sample	From Sandard Curve	

- How to convert concentration unit from mg/dl to μg/ml?
- 2 mg/dl to μ g/ml 2 x 1000/100= 20 μ g/ml

- Calculations:

- Plot the standard curve and determine glutathione concentration from the graph expressed as μg/ml.
- Calculate the glutathione concentration in Total volume of extracted plasma (= 3 ml).
 - Example: 0.3 ml of sample \rightarrow 0.2 µg
 - 3 ml of plasma \rightarrow ?
- Calculate glutathione concentration in mg by converting μg to mg.

Normal range of glutathione concentration = 3.344-4.84 mg