BCH 302 practical :

[1]Amino acid

-Qualitative tests of a.a

[2]Protein

-Qulitative chemical reactions of a.a protein functional groups

[3]Carbohydrates

- -Physical properties
- -Chemical properties
- -Qualitative analysis of CHO

[4]Lipid

-Qualitative tests of lipids

[5]Spectrophotometric DNA quantification.



Amino Acid

Mona Sh. AL-Harbi BCH302

0

*Amino acids play central roles both as building blocks of proteins and as intermediates in metabolism. There are 20 natural amino acid all of them involved in protein synthesis . All of them are L-alfa amino acids.



*All amino acids found in proteins have this <u>basic structure</u>, differing only in the structure of the R-group or the side chain. *The simplest, and smallest, amino acid found in proteins is <u>glycine</u> for which the R-group is hydrogen (H). *All these amino acids are found in solutions in their ionized form (Zwitter ion),

*The charge of a.a depend on PH of medium where they are located .





*Polar amino acids are more soluble in water than non-polar, due to presence of amino and carboxyl group which enables amino acids to accept and donate protons to aqueous solution, and therefore, to act as acids and bases. A molecule that functions as such is known as <u>an amphoteric</u>.

*<u>isoelectric point (pI)</u>: The pH value at which concentration of anionic and cationic groups are equal (i.e the net charge of this molecule equals zero)

*Amino acids are able to rotate polarized light either to the left (livo) L. or to the right (dextro) D, since they have an asymmetric C atom (a carbon atom linked to 4 different groups), glycine which lacks asymmetric C atom (has 2 H+ on a alfaC) is an exception.

Qualitative tests of amino acids:

[1]Ninhydrin test:

A test often used to detect -L- alfa amino acids.which containing free amino group.

Principle:







Contining free amino group



Dose not Contining free amino group

[2]Xanthoproteic test

This test is used to differentiate between aromatic amino acids which give positive results and other amino acids. **Principle:**



Positive result — yellow

[3] Millon's test

This test is specific for tyrosine, the only amino acid containing a phenol group, a hydroxyl group attached to benzene ring.

Principle:

In Milon's test, the phenol group of tyrosine is first nitrated by nitric acid in the test solution. Then the nitrated tyrosine complexes mercury ions in the solution to form a <u>brick-red solution</u> or <u>precipitate</u> of nitrated tyrosine, in all cases, appearance of red color is positive test. [4]Detection of amino acids containing sulfhydral group (- SH)- Lead Sulfite Test

This test is specific for –SH containing amino acid (Cystein).



PbS(black precipitate)



Report:

- -Title of the experiment
- -Name of the test
- -Objective
- -Principle
- -Material and method: as in the lab. Sheet
- -Result
- -Calculation
- -Discussion
- -Questions