




Laboratory Report (109 chem)

Experiment (8): Carbohydrates

Student Names: Section No:

<i>Compounds Reactions</i>	Observation						Conclusion
	Glucose	Fructose	Lactose	Maltose	Sucrose	Starch	
<p><u>Molisch Test:</u> 1 ml sugar solution + 1ml reagent (α-naphthol in 1ml ethanol) + 1ml conc. H_2SO_4 by drop carefully inside the wall of tube</p>							
<p><u>Reduction of Fehling Solution (A+B):</u> Fehling Solution (0.5 ml of Fehling A+ 0.5 ml of Fehling B) add to sugar solution</p>	+ve Red ppt of copper oxide appears	+ve Red ppt of copper oxide appears	+ve Red ppt of copper oxide appears	+ve Red ppt of copper oxide appears	-Ve	-Ve	
<p><u>Barfoed Test:</u> Small amount of sugar + 1ml Barfoed reagent+ heat 5min</p>	+ve Red ppt of copper oxide appears	+ve Red ppt of copper oxide appears	-Ve Red ppt. only after prolonged heating (more than 15 min) as it is hydrolyzed to glucose and galactose.	-Ve Red ppt. only after prolonged heating (more than 15 min) as it is hydrolyzed to glucose	-Ve	-Ve	
<p><u>Furfural test</u> Amount of sugar solution +3ml HCl + heat and hold paper over the mouth of tube</p>	-Ve violet color appears after some time.	+ve Violet color	-Ve violet color appears after some time.	-Ve violet color appears after some time.	+ve Violetcolor	-Ve Faint violet color after a long time	

<p><u>Iodine Test</u> <u>(I₂ Solution)</u> 1 mL sugar solution + Add 1 mL of I₂ solution</p>	-Ve	-Ve	-Ve	-Ve	-Ve		
<p><u>Tollen's reagent</u> <u>Ammoniacal silver Nitrate</u> 0.5 ml AgNO₃+ 2drop 10% NaOH (black ppt)+5drops NH₄OH add to sugar solution + heat(1-2min)</p>							
<p><u>Osazone formation:</u> Take solid sugar, phenylhydrazine hydrochlorid e, and sodium acetate in the ratio of 1:2:3. - Add 5 drops of water. - Place in a hot water bath for 10-15 minutes.</p>	<p>Yellow ppt. yellow needles aggregated in the form of sheaves.  Glucosazone or Fructosazone</p>			-Ve	-Ve		
<p><u>Water Solubility Test:</u> Amount of sugar + 1ml water</p>							

Question:

1) What is the chemical formula for carbohydrates?