
Qualitative tests of Lipids -I-


BCH302 [Practical]

Lipids

- It can be defined as **nonpolar** organic compound insoluble in polar solvent, but **soluble in organic solvents** such as benzene ,ether, chloroform.
- Lipids are esters of long chain **fatty acids** and alcohols.
- Fatty acids (F.A) are lipids' building blocks.
- There are two types of fatty acids:
- 1. Saturated fatty acids.
- 2. unsaturated fatty acids.
- Biological role of lipids:
- 1. It presents in cell membranes.

2. An essential source of energy in the body. It give **more energy** than carbohydrate and proteins.






Types of fatty acids



Saturated lipids only




Mixed saturated and unsaturated






Saturated




Monounsaturated

saturated fatty acid



unsaturated fatty acid




Classification of Lipids:

• Lipids can be divided according to their chemical composition to:

I) Simple lipids.II) Compound (conjugated) lipids.III) Derived lipids .


I) Simple lipids :

• These compounds are: esters of fatty acids with alcohol.




A. Neutral lipids:

Esters of fatty acids with alcohol (glycerol), e.g.: Triacylglyceride.



B. Waxes:

Esters are formed from fatty acids and long chain alcohol which have only one hydroxyl group (mono hydroxyl alcohol) and are larger than Glycerol (high molecular weight) e.g.: Beeswax



II) Compound (conjugated) lipids :

• Lipids are linking with <u>other compounds:</u>


a)Phospholipids:


Most phospholipids contain diglyceride, and phosphate group.

b)Glycolipids: Lipids with a carbohydrate attached.

c)Proteolipids :


Any of a group of proteins to which a lipid molecule is attached.





III) Derived lipids :


• They are substances that are soluble in lipid or derived from the lipids by hydrolysis; for examples, cholesterol and fat soluble vitamins (A, K, E and D).





Practical part

Qualitative tests of lipids



Experiment 1 : Solubility test

Objective:

• To test the solubility of oils in different solvent.

Principle:

- Fats are not dissolved in water due to their nature, which is non-polar (hydrophobic), but it is soluble in organic solvents such as chloroform, benzene, and boiling alcohol.
- Different lipids have ability to dissolve in different organic solvent.
- This property enable us to separate a mixture of fat from each other for example:
- Phosphatidelipid can not dissolve in acetone.
- Cerebroside and sphingomyline can not dissolve in the ether

Experiment 1 : Solubility test

Method:

- 1. Place 0.5ml of oil in 6 test tubes clean, dry containing 4ml of different solvents (acetone, chloroform and ether and ethanol, cold ethanol and hot water).
- 2. Shake the tubes thoroughly, then leave the solution for about one minute.
- 3. Note if it separated into two layers, the oil are not dissolve; but if one layer, homogeneous transparent formed, oil be dissolved in the solvent.

Results:

Tube	Solvent	Degree of solubility
A	chloroform	
В	ether	
С	ethanol	
D	diluted acid	



dissolved not dissolved

Experiment 2 : Saponification test


• TAG can be hydrolyzed into their component fatty acids and alcohols. This reaction can also be carried out in the laboratory by a process called saponification, where the hydrolysis is carried out in the presence of a strong base (such as NaOH or KOH).

Objective:

• To form the soap.

Principle:


- Saponification is a process of hydrolysis of oils or fat with alkaline and result in glycerol and salts of fatty acids (soap).
- Soap can be defined as **mineral salts of fatty acids**



• Soap works on <u>emulsification of oils and fats in the water</u> as it works to reduce the attraction surface of the solution.



• When soap molecules are dissolved in water, the water-repelling hydrophobic tails cluster together while hydrophilic heads surround them arranging themselves in a spherical form toward water molecules.




Experiment 2 : Saponification test

Method:

- 1. Place 2 ml of oil in a large flask.
- 2. Add 4 ml of alcoholic potassium hydroxide.
- 3. Boil the solution for 3 minutes. After this period, make sure it is perfectly saponification process, by taking a drop of the solution and mix with the water if oil separated indicates that the non-completion of the saponification. In this case, continued to boil until all the alcohol evaporates.
- 4. Take the remaining solid material (soap) and add about 30 ml of water and keep it for the following tests.
- 5. Shake the solution after it cools and noted to be thick foam.

Results:



Experiment 3 : Testing the separation of soap from the solution by salting out

Objective:

• To investigate the effect of NaCl on soap solubility.

Principle:

- To get the soap out of solution by salting out, when added solid sodium chloride to the solution until saturation; separated soap in the form of insoluble and floats above the surface.
- The NaCl solution provides Na+ and Cl- ions that bind to the polar water molecules, and help separate the water from the soap.
- This process is called salting out the soap.

Experiment 3 : Testing the separation of soap from the solution by salting out

Method:

- 1. Place about 10 ml of soap in the beaker.
- 2. Then add small amounts of sodium chloride in batches, stirring until saturated solution.

Results:

Tube	Observation
Soap + NaCl	



Experiment 4: Test formation insoluble fatty acids salt (insoluble soaps)

Objective:

• To investigate the effect of different cations on soap solubility.

Principle:

- Working calcium, magnesium, lead or iron ions leads to the **deposition of soap** making it <u>insoluble in water</u>. These ions **replace the sodium or potassium present in soap**.
- Hard water contains significant quantities of Ca2+, Mg2+ and some Fe3+ which react with the charged ends of the soaps to form insoluble salts of fatty acid.
- The insoluble salts of fatty acid that Ca2+ and Mg2+ form with soap anions cause <u>white precipitate</u> to from.

K⁺ soap + Ca⁺²sulfate=> Ca⁺² soap + K⁺ sulfate.

Potassium soap + calcium sulfate \rightarrow calcium soap + potassium sulfate. (a white precipitate from calcium stearte or oleate is formed).

17

Experiment 4: Test formation insoluble fatty acids salt (insoluble soaps)

Method:

- 1. Add about 4 ml of distilled water to 2 ml of soap in two test tubes.
- 2. Add to the first tube a few drops of calcium chloride, to second tube MgCl.

Results:

Tube	Observation
CaCl ₂	
MgCl ₂	



Howe Work:

- Why lipids are not soluble in acid ?
- What is the chemical composition of your soap?
- Why potassium hydroxide is used in saponification test?