Spect. Deter. Of Total Iron

Introduction:

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Faculty Member [BSc, MSc, DSc, Major in "Instrumental Analysis"-Expert & Advanced Major in:
FAAS FAES GC GC-MS HPLC IEC ICP-OES ICP-MS ICP-R/E FTIR & NMR]

It is well known that IRON is found in surface and drinking water and plays a vital role in biological systems, and also it is the most abundant metal ion in the human body. Many proteins use iron for oxygen transport, electron transport, and as a catalyst in oxidoreductase reactions.

In this experiment, we present a UV-visible spectroscopy assay for detection of total Fe.

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
about 1,535°C;; b.p. about 2,750°C;; sp. gr. 7.87 at 20°C;; valence +2, +3, +4, or +6. Iron is

biologically significant. Because iron is a component of hemoglobin, a red oxygen-carrying

Faculty Member [BSc. MSc. DSc. Major in "Instrumental Analysis"-Expert & Advanced Major in
FAAS FAES GC GC-MS HPLC IEC ICP-GES ICP-MS ICP-R.E FTIR & NMR
أ.د. طارق العصيمي
0114670404

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Iron is a lustrous, ductile, malleable, silver-gray metal found in Group 8 of the periodic table. It is known to exist in four distinct crystalline forms (see allotropy). The most common is the α -form, which is stable below about 770°C;, and has a body-centered cubic crystalline structure; it is often called ferrite. Iron is attracted by a magnet and is itself easily magnetized (see magnetism). It is a good conductor of heat and electricity. It displaces hydrogen from hydrochloric or dilute sulfuric acid, but becomes passive (loses its normal chemical activity) when treated with cold nitric acid.

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
اتصل بي
Turkou@outlook.com 0114670404

Iron Common Compounds:


Iron forms such compounds as oxides, hydroxides, halides, acetates, carbonates, sulfides, nitrates, sulfates, and a number of complex ions.

Compound Name	Formula	Compound Name	Formula
Iron(III) Carbonate	Fe ₂ (CO ₃) ₃	Iron(III) Perchlorate	Fe(ClO ₄) ₃
Iron(II) Sulfate Heptahydrate	FeSO ₄ .7H ₂ O	Iron(III) Sulfite	Fe ₂ (SO ₃) ₃
Iron(III) Oxide	Fe ₂ O ₃	Iron(II) Chromate	FeCrO ₄
Iron(III) Acetate	Fe(C ₂ H ₃ O ₂) ₃	Iron(II) Carbonate	FeCO ₃
Iron(II) Nitrate	Fe(NO ₃) ₂	Iron(II) Sulfite	FeSO ₃
Iron(II) Phosphate	Fe ₃ (PO ₄) ₂	Iron(II) Cyanide	Fe(CN) ₂
Iron(II) Acetate Tetrahydrate	Fe(CH ₃ COO) ₂ .4H ₂ O	Iron(II) Chloride	FeCl ₂
Iron(II) Nitrite	Fe(NO ₂) ₂	Iron(III) Nitrite	Fe(NO ₂) ₃
Iron(III) Sulfate	Fe ₂ (SO ₄) ₃	Iron(II) Oxide	FeO
Iron(III) Chlorate	Fe(ClO ₃) ₃	Iron(III) Nitrate	Fe(NO ₃) ₃
Iron(III) Phosphate	FePO ₄	Iron(III) Permanganate	Fe(MnO ₄) ₃
Iron(III) Chloride Hexahydrate	FeCl ₃ .6H ₂ O	Iron(III) Bromate	Fe(BrO ₃) ₃
Iron(III) Chloride	FeCl ₃	Iron(III) Sulfide	Fe ₂ S ₃
Iron(III) Hydroxide	Fe(OH) ₃	Iron(III) Cyanide	Fe(CN) ₃
Iron(II) Sulfate	FeSO ₄	Iron(II) Iodate	Fe(IO ₃) ₂
		Iron(II) Perchlorate	Fe(ClO ₄) ₂

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
Conversion of Iron from Fe(III) to Fe(II).
Faculty Member [BSc. VSC Post. MSc] in Instrumental Analysis - Expert & Assisted Major in:
FAAS FAES GC GC-MS HPLC IEC ICP-GES ICP-MS ICP-R/E FTIR & NMR

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
Importance of 1,10-Phenanthroline(1,10-Phen.):

Complex Compound:

Complex or coordination compound is formed when metal ion (central ion) reacts with molecules or ions which have free electron pairs (ligand). Metal ion and ligand are bonded with a polar covalent bond in which both electrons are given by the ligand [Please refer to next experiment for further information about Complex].

Experimental:

1- Prepare 50ml Of [20ppm] (Fe^{2+}) from $[(\text{NH}_4)_2\text{Fe}(\text{SO}_4)_2 \cdot 6\text{H}_2\text{O}]_{(s)}$ using distilled water.

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7- WAIT up to (5min.) to form the Complex(colored) -not in blank-.

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8- Fill all with distilled water.

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9- Take (10ml) of the unknown sample and follow steps (3-8).

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10- Move to the next Laboratory, adjust the λ_{max} to (508nm) and:

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Find the concentration of $[\text{Fe}]$ in the Unknown Sample.

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Turku@outlook.com


0114670404

Calculate the percentage compatibility/incompatibility.

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Calibration Graphs:

Comparison and Calculation of Percentage Compatibility:

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