



قسم الكيمياء الحيوية
Biochemistry Department
College of Science - King Saud University


BCH 471

Experiment (11)

**Separation of Serum Proteins by Cellulose
Acetate Membrane Electrophoresis**

Electrophoresis

- Its principle is that the charged particles of a sample migrate in an applied electrical field.
- It is applied for the separation and characterization of proteins, nucleic acids and subcellular-sized particles like viruses and small organelles.



Advantages of cellulose acetate membrane

- The virtual elimination of trailing because of the very small amount of adsorption.
- Well-defined bands are obtained on an almost colorless background making accurate quantitation possible.
- The volume of serum required is very small and with the small scale technique only a half to two hours is required for the separation.
- Subsequent staining, washing and drying is rapid so that the entire procedure usually takes no more 1-2 hours.

-Rate of migration in an electric field depends:

1- Net Charge

↑ net charge = ↑ migration

2- Size

↑ Size = ↓ migration

3-Electric field strength


↑ electric field strength = ↑ migration

4- Temperature of operation

↑ Temp = ↑ migration

Effect of pH and buffer on protein charge

- Proteins are amphoteric compounds and are therefore either positively or negatively charged
- Isoelectric Point - pH where there is no net charge in molecule.




Barbitone buffer (pH 8.60, 0.06mol/L)

Albumin	pI=4.9
α1-globulin	pI=5.28
α2-globulin	pI=5.82
β-globulin	pI=6.8
γ-globulin	pI=7.3

(Lowest mwt)

(Highest mwt)



1=albumin;

2, 3, 4, 5 is $\alpha 1-$, $\alpha 2-$, $\beta-$ and γ -globulin;

6 is origine